WECK METAL LIGATING CLIPS – ENVIRONMENTAL CONDITIONS

MR Conditional

The Weck Hemoclip, Hemoclip Plus, and Horizon Ligating Clips were determined to be MR-conditional.

Non-clinical testing was conducted using the Stainless Steel Hemoclip Plus Large Ligating Clip, which was selected in consideration of the fact that it represents the largest mass of material and the material has the highest magnetic susceptibility, such that the MRI findings apply to those similar implants that have lower masses and/or made from materials that have lower magnetic susceptibilities. This non-clinical testing demonstrated that the Hemoclip, Hemoclip Plus, and Horizon Ligating Clips are MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the Hemoclip Plus Ligating Clips produced the following temperature rise during MRI performed for 15-min of scanning (i.e., per pulse sequence) in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system:

Highest temperature change +1.6°C

Therefore, the MRI-related heating experiments for Hemoclip Plus Ligating Clips at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.6°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Hemoclip Plus Ligating Clips. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of the Hemoclip Plus™ Ligating Clips.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	571-mm ²	364-mm ²	1,109-mm ²	877-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

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