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HMC Premedical S.p.A. declares the above mentioned medical devices are conforming to the essential requirements of the Medical Device Directive consolidated with the requirements of the 2007/47/EC reported in the Annex1 and conforming to the applicable rules of product and process. The Technical File and supporting documentations are retained under the premises of the Manufacturer and at disposition of the Competent Authorities and Notified Body. The Manufacturer HMC Premedical S.p.A. has developed procedures for Post-Market Surveillance and for Vigilance according to provision of Regulation (EU) 2017/745, and it is the sole responsible for this Declaration of Conformity.

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(Harmonized European applicable rules)

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UVA PIT System

Operator's Manual CS (Closed System)

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

Contents

1	Warnings and Cautions	5
2	System Overview	7
2.1	Intended use.....	7
2.2	Note.....	7
2.3	Indications for use	7
2.4	Contra-indications	8
2.5	Side effects	8
3	Description of the UVA PIT System	9
3.1	Components of the UVA PIT System	9
3.2	Location and function of the System components	10
3.2.1	Power switch	10
3.2.2	Operator Interface.....	10
3.2.3	Irradiation chamber.....	10
3.2.4	Peristaltic Pump.....	10
3.2.5	UVA sensor	11
3.2.6	Smart card	11
3.3	Lifting and Carrying of the System	11
4	Description of the PIT KIT	12
4.1	Components	12
4.2	Pouch Label	13
4.3	Storage & use of the PIT KIT	14
4.4	Disposal of the PIT KIT	14
5	Photoactivator 8-MOP.....	14
6	Preparation of the System	15
7	Preparation of the PIT KIT	16
8	The Treatment.....	17
8.1	Priming	17
8.2	Buffy Coat Cell Suspension Transfer.....	19
8.3	Data Entry.....	19
8.4	Addition of 8-MOP.....	20
8.5	Irradiation.....	22
8.6	Rinsing Lines	23
8.7	End of Procedure.....	23
9	Disconnection.....	24
10	Return of Buffy Coat Cell Suspension.....	24
11	Manual termination of the Procedure	25

12	Alarms and Safety features	26
12.1	Overview.....	26
12.2	Mechanical control.....	27
12.3	Temperature control	28
12.4	UVA control	28
12.5	Pressure control	29
12.6	Other features for safety and protection.....	30
12.7	Pause function.....	31
13	Cleaning.....	32
14	Maintenance and Service.....	32
14.1	Lamp replacement.....	33
14.2	Preventive Maintenance/Safety check	33
14.3	Provision of Documentation.....	33
14.4	Lifetime of the device.....	34
14.5	Disposal of the device	34
14.6	Fuse	34
15	Signs and symbols	35
15.1	General	35
15.2	PIT System	35
15.3	PIT KIT	35
15.4	Packaging.....	36
16	Technical specification	37
16.1	UVA PIT System	37
16.2	EMC UVA PIT SYSTEM.....	38
16.3	PIT KIT	39
17	Manufacturer and Certifications.....	40
17.1	Manufacturer	40
17.2	Authorized Representative.....	40
18	Registrations / Certifications.....	40
18.1	UVA PIT System	40
18.2	PIT KIT	40
19	Limited warranty	41
20	Documentation of the Procedure	42

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

This manual must not be copied and must not be given to non-users.

Copyright PIT Medical Systems GmbH.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

1 Warnings and Cautions

	<p>WARNINGS</p> <p>Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.</p>
	<p>CAUTIONS</p> <p>Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.</p>
	<p>IMPORTANT</p> <p>Instructions including important information. To avoid damage, injury, or improper use of the product.</p>

Consult this manual prior to use.

PIT Medical is only responsible for the basic safety, reliability, and performance of the UVA PIT System if it used in accordance with this operator's manual. Operator's Manual must be kept with the device.

Please read and observe these warnings and safety instructions before use.

- Please note that these operator's manual is a general guide for the use of the product. Medical matters must be addressed by a physician.
- Consult the indications for use, contraindications, warnings and cautions before using - see chapter [System Overview](#)
- Do not modify this equipment without authorization from the manufacturer.
- Only use properly grounded receptacles and circuits designed to accommodate this device.
- Position the system securely using all four rubber feet.
- To avoid any risk of electrical shock, this equipment must be connected to an electrical power supply with a protective earth (PE).
- The equipment should be used only with the power cord that is supplied by the manufacturer.
- Do not use multiple plug extension.
- Prevent liquids from entering the system. The saline bag may have only a maximum of 250 ml. Please use one for priming and one for rinsing.
- Do not operate the system in the presence of high oxygen concentration, anesthetic gases, or other easy flammable gases.
- For any unexpected device condition, shut off the device using the main switch on the back.
- Locate the UVA PIT System in a room with a maximum temperature of 28°C and humidity less than 80%.
- The device is stationary used in medical environment.
- The device is not in physical contact with the patient. PIT Medical Systems recommends the use of the UVA PIT System should occur at least 1.5 m away from the patient.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

- The system has two fans for cooling on each side. Do not cover the fans and allow a minimum space of 40cm around all sides of the system.
- The UVA PIT system must remain in an upright (display in front) position during use.
- To avoid any risk of crushed fingers, do not reach into rotating parts. Visually inspect the motor cover before use.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the UVA PIT System, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation .
- The UVA PIT System is an “eye-safe” device when operated properly, however, as with any light source, never look directly at the light output.
- For the operator: There is no need to use UVA protective glasses during normal use.
- For the patient: UVA protective glasses should be used mandatory if the patient takes oral 8- MOP. To reduce the risk of any increased residual amount of 8-MOP after reinfusion of the buffy coat cell suspension, the patient should wear UV-blocking sunglasses and protect from the sun 24 to 48 hours after ECP treatment. Follow instructions in the package/ patient information leaflet from the pharmaceutical manufacturer of the 8-MOP.**The working instructions for the procedure, including the information of the patient information leaflet of the manufacturer of the 8-MOP, and recent guidelines must be followed strictly.**
- During the whole procedure, sterile working is mandatory.
- PIT Medical can only guarantee the safe functioning of the system when used in combination with the original PIT KIT accessory. **Only the disposable KIT (exclusively manufactured by HMC REF PIT-KIT/CS is approved to be used with this system.** Other disposables may not lead to the expected result and lead to the expiry of the regulatory approval of the product.
- **Please consult the instructions for use on the pouch label of the PIT KIT and read the warnings and safety instructions before use.**
- Do not use the PIT KIT if:
 - The sterile package is damaged.
 - The package seal is open.
 - The PIT KIT is wet.
 - The expiration date has passed.
- The Procedural KIT is intended for a single patient’s entire therapy. Reuse on another patient might cause cross contamination.
- Always check the patient ID (barcode sticker) against the ID (barcode sticker) on the Recirculation bag (RC-bag) before beginning the reinfusion.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

2 System Overview

The UVA PIT System has been developed to enhance the UVA irradiation.

The UVA PIT System incorporates technical advances in the fluid irradiation functions:

- Continuous recirculation
- Offline recirculation
- UVA control
- Simple setup and input of parameters via touch screen
- Automatic energy calculation
- Temperature control and UV monitor function

A pre-connected, single use disposable has been specifically designed for the use with the UVA PIT System. Only the **PIT KIT** should be used with the UVA PIT System.

2.1 Intended use

The UVA PIT System provides Photo Immune Therapy on autologous blood.

The UVA PIT System is CE certified as a medical system.

The UVA PIT System is indicated for the use in the UVA irradiation (in the presence of the photoactive substance 8-MOP) of extracorporeal circulating leukocyte enriched blood.

Collection of the leukocyte-enriched blood to be treated with the UVA PIT System can be performed with any of the existing standard equipment available on the market.

PIT Medical Systems recommends the use of the UVA PIT System should occur at least 1.5 m away from the patient.

2.2 Note

Extracorporeal photopheresis (ECP) should be performed only by persons with special training (nurses and doctors) and in institutions with access to the necessary equipment and materials for this treatment.

Therapy using psoralen/8-MOP and UV irradiation should take place under constant supervision by a physician with the appropriate training.

2.3 Indications for use

- Cutaneous T-cell lymphoma (CTCL)
- Chronic Graft-versus-Host disease (cGvHD)

Please consider actual clinical and scientific data, publications, and guidelines for other possible applications of ECP.

The Operator's Manual concerning this matter is updated regularly.

The decision to treat a patient with ECP should be made by a specialist/physician.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

2.4 Contra-indications

As of the date of this publication, the ECP is contraindicated for patients with the following conditions based on single adverse events reported in connection with each condition.

- Hypersensitivity to Psoralen/8-MOP
- Comorbidities, which can lead to increased photosensitivity
- Aphakia
- Coagulation disorders such as heparin-induced thrombocytopenia (HIT)
- Cardiovascular diseases
- Low hematocrit values
- Incompatibility of an extracorporeal blood volume loss, e.g. in case of severe cardiac, hepatic or renal insufficiency or a low body weight (for more information concerning body weight, see: Messina et al., 2003 and Kanold et al., 2007). This issue depends on the used apheresis device. In general, the patient's total blood volume determines the safe extracorporeal blood volume, which should not exceed 15% of blood volume (Howell et al., 2015).
- Pregnancy and lactation (at present, insufficient information is available)
- In addition, patients with Cold agglutinin disease (CAD) or with an abnormal 8-MOP-metabolism or -transport might be contra-indicated for ECP. Both reported conditions were individual single cases (Mask-Bull et al., 2014; Vagace et al., 2007).
- The instructions within the patient information leaflet of the 8-MOP of the pharmaceutical manufacturer (warning and precautions, possible side effects, and interaction with other medicinal products) must be strictly followed.

Messina et al., Br J Haematol. 2003 Jul;122(1):118-27; Kanold et al., Transfusion. 2007 Dec;47(12):2276-89; Howell et al., Transfus Med. 2015 Apr;25(2):57-78; Mask-Bull et al., J Am Acad Dermatol. 2014 Mar;70(3): e61-2; Vagace et al., Arch Dermatol. 2007 May;143(5):622-5

2.5 Side effects

Adverse reactions / side effects of ECP are generally very low.

- Transient hypotension
- Low to moderate temperature increases
- Increased erythema
- Fatigue and lethargy that usually last for 24 to 48 hours after each treatment
- Anemia, especially during long-term ECP applications by inadequate and incomplete reinfusions and rarely due to hemolysis
- Reversible hypertension
- Infections
- Thrombocytopenia

As of the date of this publication, there are no reports of serious side effects of grade III or IV on the WHO classification scale.

3 Description of the UVA PIT System

3.1 Components of the UVA PIT System


Figure 1: Components of the UVA PIT System

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

3.2 Location and function of the System components

3.2.1 Power switch

The device must be plugged into electrical outlet for the ON/OFF power switch to function.

The power switch is located on the rear side of the instrument above the power cord connection.

IMPORTANT

The unit must not be installed in such a way that the ON/OFF switch is difficult to operate, or the power supply cable cannot be removed from the unit.

3.2.2 Operator Interface

The operator interface consists of a touch screen.

The touch screen displays the following information:

- Procedure status
- Remaining energy to be delivered (in Joules)
- Temperature (in °C)
- Status of the lamp energy

When the touch screen has been inactive for a specified amount of time (not during a procedure), it goes blank. Touch the screen to reactivate the display.

Audible and visible alarms are used to alert you to special operating conditions.

3.2.3 Irradiation chamber

The UVA PIT System controls the real UVA output energy, which delivers the optimum UVA energy for effective irradiation. UVA output energy during an ECP procedure is targeted to be 2 J/cm². The UVA PIT System light source consists of six lamps (three at the top of the irradiation chamber and three at the bottom).

Approximately 500 hours of life can be expected from the lamps. The intensity of the lamp energy is electronically monitored, and the values are tracked by the UVA PIT System and constantly provided in the irradiation algorithm.

The UVA PIT System automatically compensates for the aging of the lamps. The lamp intensity is displayed on the touch screen, which indicates when the lamps should be replaced. See chapter [Lamp Replacement](#).

Inside the photo activation chamber is an acrylic glass which is consisting of a special UVA translucent material. Its function is to hold the irradiation bag. Do not attempt to replace this glass. Replacement of this glass may interfere with the UVA transluence for the treatment of the buffy coat cell suspension.

3.2.4 Peristaltic Pump

The peristaltic pump continuously recirculates the buffy coat cell suspension through the irradiation bag for a uniform exposure to the UVA light.

The peristaltic pump is located at the front of the system and moves the fluid through the disposable set lines. The pump itself consists of two parts: the drive containing the electronics, which is inside the unit, and the rotor, protected by an outer cover. The rotor is designed to gently move the blood and does not damage the structure of the cells.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

3.2.5 UVA sensor

The UVA sensor, which is located in the photoactivation chamber, detects the amount of energy delivered by the UVA lamps. The sensor determines the actual energy that is being delivered by the UVA lamps and provides data back to the system to adjust the energy, as necessary.

3.2.6 Smart card

Each box of PIT KITs (25 per box) is equipped with a smart card that contains 25 procedures. When you press **Start** on the touch screen (after entering data and adding 8-MOP), one procedure is deducted from the smart card. When you press **OK** at the end of the procedure, the touch screen displays the remaining number of procedures on the smart card.

3.3 Lifting and Carrying of the System

The UVA PIT System weighs about 21kg. Lift it with one hand placed on each side between the device feet.

IMPORTANT

Before lifting/carrying the device:

- Turn off the device at the main switch.
- Disconnect the power cable from the device.
- Remove the pole by lifting it vertically.
- Close the motor cover.
- Close the irradiation tray.
- Remove the smart card.

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	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

4.3 Storage & use of the PIT KIT

Store the PIT KIT away from moisture and within the temperature range of -5°C to $+35^{\circ}\text{C}$.

Always follow the instructions on the carton box and on the pouch label- see also the chapters:

[Warnings and Cautions](#), [Signs and symbols](#).

IMPORTANT

Open the sterile packaging just before use and handle the PIT KIT using aseptic techniques.

4.4 Disposal of the PIT KIT

The used disposable should be considered as biological hazard. Follow the policies and procedures of your institution's infection control plan for the handling and disposal of blood contaminated materials.

5 Photoactivator 8-MOP

The UVA PIT System uses 8-MOP (8-Methoxypsoralene) as a photoactivator.

8-MOP is a light sensitive molecule and attaches to the DNA of cells. The 8-MOP molecule is activated by UVA light and cross-links the DNA to inhibit replication. Apoptosis is induced within the treated buffy coat cell suspension.

IMPORTANT


Follow the instructions in the package/patient information leaflet from the manufacturer of the 8-MOP.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

6 Preparation of the System


- To start the system, ensure that the power cord is connected to the wall socket, then turn on the device using the main switch power switch located on the back.

A control signal (beep) sounds. The touch screen displays “**Insert Card**” and the installed firmware version.

- With the chip facing up, insert the smart card into the smart card reader on the right side of the device.

If the card is empty or invalid, the screen is displayed:


If the card is valid, the touch screen displays the number of available treatments (procedures) (25 for a new card):


7 Preparation of the PIT KIT

- Inspect the PIT KIT. If there are loose parts, holes or any other irregularities, **do not use the PIT KIT**.
- Attach one of the two barcode stickers to the appropriate area of the RC bag (recirculation bag) and the other one to the related document or to the patient file.
- **Close all 5 pinch clamps (red, blue, white, yellow, green), 2 slide clamps (white) and the roller clamp (blue).**
- Position the RC bag on the lower hook of the pole (see above Fig. 1, No. 3).
- Hang saline bag on the upper right hook and spike with the saline line (Fig 2, No. 16).
- Prepare the collected cell bag.
- Remove the cover from the spike on the buffy coat line (Fig. 2, No. 1; the spike connector is red). Spike the cell product bag with the spike or connect it with a sterile connection device and place it on the left upper hook (Fig 1, No. 2).
- Slide out the irradiation tray.
- Place the irradiation bag onto the irradiation plate with the pump segment on the right front side. Slip the bag under the pressure sensor, as shown below. Secure the bag holes over the four pins shown below, so it stays in place.

- Open the pump cover
- Insert the pump segment into the pump (Fig. 1 No. 6) with the red dot marking to the left side and then turn the handle of the pump manually to insert the tube into the pump housing.

- Before closing the motor cover, **check the correct position of the pump segment/tube inside the pump: The tube should not be twisted and should be flush in the guideway/guiding groove.**
- Close the pump cover carefully and make sure that the **cover is closed completely and aligned to the pump housing.**
- Slide the tray back into place.


	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

8 The Treatment


8.1 Priming

- After the PIT KIT is positioned on the irradiation tray, the pump is loaded, and the pump cover and tray are closed, press **START**.



The device runs for approximately 12 seconds to adjust the pump segment of the PIT KIT. Meanwhile this screen is displayed:



Next, the system prompts you to check the lines/tubes and the pump segment to be sure they are properly seated.

- If the lines, especially the pump segment, are properly seated, press **OK** to continue to the priming screen.
- If the lines or the pump segment are not correctly inserted/adjusted (e.g., if they are twisted), press **NOK** (= not OK).

If you press **NOK**, the device prompts you to remove the tube from the pump and reload the pump segment.


- Reload the pump segment, then press **OK**.

The pump check/adjustment is performed again. (Refer to the steps above).


If you press **STOP**, the procedure is aborted. [Note: No procedure is deducted from the smart card.]

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

- Select whether you will prime the set (Yes or No).


- If you have selected **No**, continue with the next section: [Buffy Coat Cell Suspension Transfer](#).
- If you have pressed **Yes**, open the **blue** and the **white** pinch clamps.
- Select the priming volume: 40ml or 80ml.

The pump runs counterclockwise and stops automatically after priming the PIT KIT with the volume of saline you selected.


- During the countdown, you can press **Skip** to skip the remainder of the priming or **Pause** to pause the priming. Information about the pause function can be found in the [Pause function](#) section.

- After the pump has stopped, or if you have chosen to skip the priming, close the **white** pinch clamp.

The tubing contains 40 or 80 ml of saline solution, depending on which volume you chose.

This screen is displayed:

Remark: If you have selected 80 ml priming volume, you have the option of priming the reinfusion line.

- Open the **yellow** pinch clamp and the roller clamp (Fig. 2, No. 7, 9).
- Remove the cap from the male luer lock (Fig. 2, No. 10).
- Let some saline solution pass through the reinfusion line.
- **Close** the **yellow** pinch clamp and the roller clamp.
- **Attach** the cap firmly to close the luer lock.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

IMPORTANT

Once the system has been primed, **do not** turn it off.

- Before continuing to the next step, cell transfer, inspect the PIT KIT again for overall integrity and for leakage, especially at the welded components (tubing connections) and the pump segment.
- Open the tray and inspect the IR bag. If you chose to prime the PIT KIT, the IR bag should be partly or completely filled with saline solution, depending on whether you chose 40ml or 80ml of priming volume.
- If an IR bag does not open during the priming, repeat the priming. If you detect any other irregularities, discard the KIT and use a new one.

8.2 Buffy Coat Cell Suspension Transfer

Close the **blue** pinch clamp and open the **red** pinch clamp.


- Transfer all cells by gravity to the RC bag.
- Close the **red** pinch clamp.

The RC bag has two latex free sampling ports (Fig. 2, No. 11). Use these ports to take a sample, if needed. Continue with data entry.

8.3 Data Entry


The touch screen displays a default volume of 150 ml.

- Use - or + button to adjust the **volume to match the volume of the collected buffy coat cell suspension**. It can range from 50 to 800 ml. If it is less than 50 ml, enter 50 ml.

- Press **OK**.


The touch screen displays a default hematocrit % of 1.0. Use the + and – buttons to adjust the **hematocrit % to match the hematocrit % of the collected buffy coat cell suspension**. It can range from **0.5 % to 6 %**. If it is less than **0.5 %**, enter **0.5 %**.

- Press **OK**.


The touch screen displays the data you have entered.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

- If the data is correct, press **Yes** to confirm.
- If the data is not correct, press **No** and follow the steps to re-enter the volume and hematocrit %.

If you press **Yes**, this screen is displayed:

8.4 Addition of 8-MOP

The touch screen displays the minimum of 8-MOP that should be added based on the volume and hematocrit % you have entered. However, you should strictly follow the instructions in the package/patient information leaflet from the manufacturer of the 8-MOP.

As a calculation example, taken from a PIL of 8-MOP with an initial concentration of 20 µg/ml Methoxypsoralen. The dosage is calculated according to the treatment volume, using the following formula:

$$\text{Apheresis collection volume in ml} \times 0.017 = \text{volume of recommended 8-MOP in ml.}$$


Figure 3: Schematic diagram, detail of the UVA PIT KIT

- Attach one of the two stickers **8-MOP inserted** to the appropriate area of the RC bag (recirculation bag) and the other one to the related document or to the patient file.
- Please shake the bag gently after adding 8-MOP.

The **Sample taking ports** are provided for taking a sample, if necessary (see below and Fig. 2, No. 12).


Figure 4: Schematic diagram, detail of the UVA PIT KIT

IMPORTANT

Before continuing with the next step:

- Check all lines.
- Make sure that the **white** slide clamps are closed.
- Make sure the **red** and the **white** pinch clamps are closed.
- Open the **green** and **blue** pinch clamp on the PIT KIT.


	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

8.5 Irradiation

- To start the procedure, press **Start**.

The system calculates the required energy based on the specific volume and hematocrit % you entered for the collected buffy coat cell suspension and starts the irradiation phase.

The touch screen displays the total required energy in Joules for this buffy coat cell suspension:

The UVA lamps start to irradiate the buffy coat cell suspension and the touch screen displays the remaining energy to be delivered (see the large numbers at the top of the display).

The pump runs counterclockwise.


The buffy coat cell suspension is circulated from the RC bag via the pump through the IR bag and back to the RC bag.

Touch screen:

- On the left side of the touch screen, the chamber temperature is displayed in °C. See additional information on: [Alarms and security features: Temperature control](#).
- On the right side, the irradiation UV power, expressed as energy per area (E/A), is displayed as a color-coded bar. The E/A-bar may be green, yellow or red, depending on the actual lamp level. For this, see [Alarms and security features: UVA control](#).
- You can press **Pause** to pause the irradiation. The pause function is described in the [Alarms and security features: Pause function](#) section.

When irradiation is complete (the calculated energy has been delivered), the device stops irradiation and continues to the next step.


The screen is displayed:

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

8.6 Rinsing Lines


- Close the **green** pinch clamp.
- Open the **white** pinch clamp.

- Press **Rinsing**.

The pump runs counterclockwise and purges all cells to the RC bag.

An additional volume of about 50ml saline solution will be added to the whole buffy coat cell suspension (for complete rinsing).


- To pause the rinsing, press **Pause**.

See the description in the [Alarms and security features: Pause function](#) section.

- To skip rinsing, press **Skip**.
- To skip rinsing during the countdown, press **Skip** for at least 3 seconds.

8.7 End of Procedure

After you either skip rinsing or allow it to complete, this screen is displayed:

- Press **OK**. The device returns to the **Start** screen and displays the number of treatments (procedures) remaining on the smart card.
- Close the **white** and the **blue** pinch clamp.

9 Disconnection

IMPORTANT

Ensure that all 5 pinch clamps and the 2 slide clamps are closed.

- Separate the RC bag from the rest of the PIT KIT CS by using tube sealers.
- The RC bag is now ready for further handling.
- Remove the IR bag from the IR chamber and the pump segment from the pump.
- To dispose of the PIT KIT, follow your facility's SOP for handling hazardous blood contaminated materials.

10 Return of Buffy Coat Cell Suspension

The reinfusion line (with blood filter, roller clamp a luer connection; see Fig. 5, No. 8, 9 and 10) is designated to return the treated cells in a safe and sterile manner.

- Prepare the patient for reinfusion according to your facility's SOP.
- Hang the RC bag on a pole. Remove the paper tape and uncoil the line. Using the luer connector, connect the reinfusion bag to the patient access according to your facility's SOP.
- Open the **yellow** pinch clamp (Fig. No.7) and allow the cells return by gravity. Adjust the speed by using the roller clamp. Follow your facility's SOP for reinfusion.


Figure 5: Schematic diagram, detail of the UVA PIT KIT CS

IMPORTANT

When the reinfusion is complete, see chapter [Disposal of the PIT KIT](#).

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

11 Manual termination of the Procedure

IMPORTANT



If the power fails, the system stops immediately and returns to the **Start** screen. The operator or physician should determine whether to return the cells manually to the patient or perform the procedure again.

Follow these general guidelines to end the procedure manually:


- Close the **blue** pinch clamp.
- Make sure the **red** and the **white** pinch clamps are closed too.
- Remove the IR bag from the IR chamber.
- Remove the pump segment out of the pump.
- Make sure the **green** pinch clamp is open.
- Lift the IR bag above the liquid level of the RC bag and let the cells flow to the RC bag. You may apply some slight pressure to the IR bag.
- Close the **green** pinch clamp , if the IR bag is empty.
- The RC bag could be separated from the whole PIT KIT by using tube sealers. As during the whole treatment, pay attention to sterility.
- Proceed with chapter [Return of Buffy Coat Cell Suspension](#).

12 Alarms and Safety features

12.1 Overview

Message	Description	Solving
	<p>This message appears if the drawer is open or not completely closed</p>	<p>To rectify the error message, please close the drawer completely.</p>
	<p>This message appears if the pump cover is not closed. A red light is illuminated on the pump.</p>	<p>To rectify the error message, please close the motor cover completely. A green light is illuminated on the pump.</p>
	<p>At the time the IR chamber temperature exceeds 43°C, the device terminates the process.</p>	<p>By pressing "OK", the device will go into the start mode and will show the remaining number of treatments, but the message "Temperature Error" is still displayed. To delete the error message, power cycle the device.</p>
	<p>If the UVA intensity drops below a certain value, the E/A-bar changes from green to yellow indicating that the lamp level is low, and the lamps have to be replaced soon. If the lamp level drops further below a threshold, the E/A-bar changes to red.</p>	<p>Please call your local representative in case this happens.</p>

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

	<p>The UV intensity is below the minimum value for longer than 30s. No more treatment can be performed.</p>	<p>Please call your local representative in case this happens.</p>
	<p>The UVA PIT KIT is under too high pressure and there is no circulation.</p>	<p>At this point, the operator should check the status of the manual valves / clamps and inspect the lines for a kink or an obstruction.</p> <p>Then press OK.</p>

12.2 Mechanical control

If the irradiation tray is not in place, the UVA lamps switch off and the pump stops.


The screen is displayed:

- After closing the irradiation tray, press **START** to continue.

The irradiation tray's position is controlled by a magnetic switch. This magnetic switch enables the correct operation of the device and is not an interlock to prevent radiation emission.


If the pump cover is open, this screen is displayed:

- Close the motor cover and make sure it is closed completely and aligned to the pump housing. Press **START** to continue.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN


12.3 Temperature control

If the temperature reaches 38 °C, the numbers in the Celsius field turn red (see above). As soon as the temperature exceeds 40 °C, the device starts to beep continuously, and the UVA irradiation stops. The fans and the pump continue to operate in order to cool down the IR chamber.

After the temperature drops below 38 °C again, the beeping stops, and the device automatically continues the procedure.

If the IR chamber temperature exceeds 43 °C, the device terminates the procedure, and this screen is displayed:


In this case, the product should be tested for hemolysis. The deactivation of the 8-MOP might have not been complete. A physician or specialist must determine whether the buffy coat cell suspension can be reinfused.

- Press **OK** to return to **Start** screen. The touch screen displays the remaining number of procedures, but the message “Temperature Error” is still displayed. To clear the error message, power cycle the device.

12.4 UVA control

An internal UVA sensor continuously controls the emitted UVA irradiation. The processor deducts the emitted energy from the total calculated energy. This control loop ensures to provide a consistent process.

If the UVA intensity drops below a certain value, the E/A-bar changes from green to yellow indicating that the lamp level is low, and the lamps must be replaced soon. Please call your local representative in case this happens.

The device will allow you to continue the treatment, as well as a limited number of further treatments are still possible.

In order to remind the operator, the indication “Lamp Level Low” will be shown during pause and in the next start screen (after power cycle, the message is deleted).


	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

If the lamp level drops further below a threshold, the E/A-bar changes to red.

After 30 seconds below this threshold the procedure will be aborted, and this screen is displayed:


In this case, the lamps must be replaced **immediately** in order to continue working.

If the procedure is completed, paused, or aborted within the 30 seconds after the E/A-bar has turned red, the message "Replace Lamps" is also displayed.

See Chapter [Lamp replacement](#).

12.5 Pressure control

The UVA PIT system is equipped with a pressure sensor that monitors the pressure in the PIT KIT tubing.

If you forget to open the clamp, or if a line is kinked or obstructed, the device stops the pump to avoid excessive pressure in the recirculation system.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

In this situation, the device starts to beep continuously, and this screen is displayed:

IMPORTANT

Check the status of the manual valves/clamps and inspect the lines/tubing set for a kink or an obstruction.

- Once the error is resolved, press **OK** to continue. The pressure sensor is for correct operation and not an interlock.

12.6 Other features for safety and protection


- The IR bag is fixed at four pins on the irradiation tray to ensure it is properly seated.
- The front slit is covered to avoid UVA exposure of the operator.
- A Smart card enforces correct use of the device.
- Six long life UVA lamps have a homogeneous irradiation profile; they have a life span of approximately 500 hours.
- A special peristaltic pump protects cells from damage.
- Two high performance cooling fans maintain temperature.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

12.7 Pause function


Any time **Pause** is displayed on the touch screen, you have the option to pause the operation of the device. Depending on the phase of the procedure, different screens are displayed after you press **Pause**.

As shown below, you may have the option to press **START** to continue the procedure during the irradiation phase, **Priming** to continue priming, or **Rinsing** to continue rinsing. You can also press **STOP** or **Skip**.

- During the irradiation phase of the procedure, **STOP** aborts the procedure. Press **STOP** for more than three seconds, then release.

This screen is displayed:

- During priming, **STOP** aborts the procedure. Press **STOP** for more than three seconds, then release. The device returns to the **Start** screen. Note: No procedure is deducted from the smart card if the procedure is aborted during priming.
- During rinsing, **Skip** ends the procedure. Press **Skip** for more than three seconds, then release. The touch screen displays "Treatment Done".

Opening the irradiation tray or pump/motor cover has a similar effect as using the pause function. The screens that are displayed after you open the motor cover are shown in the [Alarms and Safety features](#) chapter.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

13 Cleaning

WARNING

Always disconnect the device from power before cleaning.

Exterior of the device: As required; with soapy water

IMPORTANT

- Do not use chemicals
- Do not use excess water

Irradiation chamber: As required; with soapy water

IMPORTANT

- Do not clean with chemicals or grinding detergents.
- Do not use excess water
- Do not pour fluid into the photo activation chamber.
- Only use lint-free soft-tissues.

If disinfection should be required, only alcohol-free and aldehyde-free disinfectants suitable for acrylic glass should be used. **Do not use excess of disinfectant and do not directly spray onto touch screen.**

14 Maintenance and Service

If the device shows an error or does not perform properly, please call your local representative or PIT Medical Systems. The device must not be opened by non-trained personnel with no training certification from PIT Medical Systems. If this occurs, the distributor and manufacturer will not be responsible for any malfunction, errors or any other incidences caused by the device. The device must undergo complete, rigorous testing before further use.

WARNINGS

- Only persons authorized by PIT Medical Systems may perform service inside the protective covers of the UVA PIT System. Dangerous voltages are present inside the system. Authorized persons should perform maintenance after the system is shut down and disconnected from power.
- Turn off and unplug the UVA PIT System before servicing. Never perform maintenance on the UVA PIT System during a procedure.

Maintenance procedures are to be performed only by fully trained personnel.

This Service technician must have the following qualifications:

- Principles of electrical engineering
- Knowledge in dealing with UVA radiation
- Completed product and service training at PIT Medical Systems GmbH

Contact PIT Medical Systems or your local representative once a year or every 500 procedures (whatever comes first).

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

For a detailed technical, service and maintenance description please refer to the Service & Maintenance Manual. If you do not have the document, ask the manufacturer for it.

14.1 Lamp replacement

The UVA lamps should be replaced, at the latest, when indication on the touch screen. See also [UVA control](#).

The emitted lamp energy is displayed during irradiation on the right side of the touch screen as E/A (energy/area). Check the lamp energy only when the device is warmed up (after at least seven minutes of irradiation).

When the E/A-bar changes from green to yellow (indicating that the lamp level is low), contact your local representative to replace the lamps.

The device allows you to continue the procedure, as well as a limited number of further procedures.

If the energy drops further below a threshold, the E/A-bar changes to red. After 30 seconds, the touch screen turns red and displays the alarm “**Error – UV Level Low**”, the lamps turn off and the procedure is aborted.

In order to continue working, the lamps must be replaced.

Only authorized personnel should replace the lamps.

Use only lamps provided by PIT Medical Systems.

14.2 Preventive Maintenance/Safety check

This device must be checked annually.

Related to the Medical Device Directive, the operator is responsible for initiating the preventive maintenance and safety check and calling the service engineer.

PIT Medical Systems will not take any responsibility and will decline every liability if preventive maintenance and safety checks are not performed in time.

The service procedure shall perform an inspection shall ensure that the PE (Protective Earth) is correctly connected.

14.3 Provision of Documentation

If additional documents such as circuit diagrams, component lists, descriptions, calibration instructions or support are required, please contact us. We will provide you with these documents.

Contact details can be found under chapter [Manufacturer](#).

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

14.4 Lifetime of the device

The intended life span of the device is seven years.

14.5 Disposal of the device

The lamps should be removed and disposed of according to your facility's SOP gas containing lamps.

The device must be disposed of as electronic waste.

14.6 Fuse

The device is equipped with a replaceable fuse:

Fuse:	T2A
- Rated Voltage	250VAC
- Rated Current	2.0A
- Release Speed	Time-Lag slow-blow
- Breaking Capacity	35A – 200A

The fuse is located at the main switch.




	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

15 Signs and symbols

15.1 General

	This symbol indicates the name and the address of the manufacturer.
	This symbol indicates the date of manufacture (four digits for the year).
	This symbol indicates that the device is in conformance with the Medical Device Directive 93/42/ EEC.

15.2 PIT System

	This symbol indicates to not dispose the device together with household refuse (for EU only).
	This symbol indicates a CAUTION or WARNING associated with the device.
	This symbol indicates manufacturer's serial number.
	This symbol indicates to refer to instruction manual/booklet/operator's manual.
	This symbol indicates to not reach into rotating parts.
	This symbol indicates the minimum and maximum ambient temperature during storage / transport
	This symbol indicates the minimum and maximum ambient air pressure during storage / transport / operation.
	This symbol indicates the maximum humidity during operation.
	This symbol indicates that this connection may only be used for service purposes.

15.3 PIT KIT

	The symbol indicates that the product must not be used after the end of the month.
	This symbol must be accompanied by the manufacturer's code for the batch. Synonyms: "lot number" - "batch number"
	To indicate that the device is provided sterile and has been sterilized using irradiation.
	This symbol identifies the product as non-reusable. Synonyms of "non-reusable" are: "single use" and "use only once"
	This symbol must be accompanied by the code identifying the product.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

	This symbol identifies the product as not re-sterilizable for re-use
	Latex free
	DEHP free
	Medical Device
	To indicate that caution is necessary when operating the device.

15.4 Packaging

	This symbol indicates to keep the goods away from sunlight.
	To indicate that the transport package shall be kept away from rain and in dry conditions. Transport packaging design and materials handling
	To indicate correct upright position of the transport package.
	Indicates a medical device that can be broken or damaged is not handled carefully.
	To indicate that hooks shall not be used for handling the transport package.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

16 Technical specification

16.1 UVA PIT System

Physical	
Dimensions (Height x Width x Depth)	64cm x 30 cm x 23 cm
Weight:	21 kg
Recommended operating space:	Minimum of 40 cm clearance on all sides
Power Cord Length:	1, 8 m
Electrical	Power supply: 230 VAC, 50/60 Hz;
	Power consumption: 300 VA
Recirculation pump	Flow rate: 80 ml/min
Environmental conditions: Operation	Temperature: 18°C - 28°C Humidity: < 80% Air pressure: 700hPa to 1060 hPa
Environmental conditions: Transport	Temperature: -20°C - 50°C Humidity: < 80% Air pressure: 300hPa to 1060 hPa
Environmental conditions: Storage	Temperature: -20°C - 50°C Humidity: < 80% Air pressure: 300hPa to 1060 hPa
Switch-on cycle:	Max. ON-D 90min / Min. OFF-D 20min
Protection class	Class I
IP-Class:	IP X0
Regulations & norms	IEC 60601-1-2
	IEC 61010-1
	IEC 60601-1

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

16.2 EMC UVA PIT SYSTEM

Environment		Professional health care facility except in the vicinity of active devices of the HF surgical devices and MRT devices
Emissions		
Class		CISPR 11, Class B, Group 1
Harmonic Distortion (IEC 61000-3-2)		Comply with IEC 61000-3-2 class a
Voltage fluctuation and Flicker (IEC 61000-3-3)		Comply with IEC 61000-3-3
Immunity		
Electrostatic discharge (IEC 61000-4-2)		CD +/- 8 kV, AD +/- 15 kV
Radiated RF EM Fields (IEC 61000-4-3)		80 – 2700 MHz: 3 V/m
Proximity Wireless fields (IEC 61000-4-3)		385 – 5785 MHz: 9 – 28 V/m
Electrical Fast Transients and bursts (IEC 61000-4-4)		+/- 2 kV: Power supply lines
Surges (IEC 61000-4-5)		+/- 1 kV Differential mode +/- 2 kV Common mode
Conducted disturbances, induced by RF fields (IEC 61000-4-6)		3 V rms and 6 V rms in ISM bands
Rated Power-frequency Magnetic Field (IEC 61000-4-8)		30 A/m
Voltage dips (IEC 61000-4-11)		0% U_T ; ½ period At 0, 45, 90, 135, 180, 225, 270 and 315 degree 0% U_T ; 1 period and 70% U_T ; 25/30 periods Single phase: at 0 degree
Voltage interruptions (IEC 61000-4-11)		0% U_T ; 250/300 periods

Essential performance: By applying risk analysis, PIT Medical has determined that the UVA PIT System does not have any performance characteristics that could lead to destruction of blood cells caused by UV underexposure or overexposure. As such, the UVA PIT System is not considered to have essential performance as defined in the IEC 60601-1-2 standard. However, the UVA PIT System was tested per IEC 60601-1-2 ensuring that the treatment could be maintained during testing.

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

16.3 PIT KIT

Packaging Unit	Box with 25 KITs <ul style="list-style-type: none"> - Single packed - sterile - 1 Smart Card (for 25 procedures)
Weight of box	6 kg
Size of box	58 cm x 40 cm x 26 cm

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

17 Manufacturer and Certifications

17.1 Manufacturer

PIT Medical Systems GmbH
Am Farrnbach 1
90556 Cadolzburg
Germany

17.2 Authorized Representative

Mr. Waldfried Weber
Am Farrnbach 1
90556 Cadolzburg
Germany

18 Registrations / Certifications

18.1 UVA PIT System

DIMDI	UVA PIT System	DE/CA59/5728/2016-R/Hd
UMDN	Irradiators, Blood	17-437
CND	HEMATOLOGY AND HEMOTRANSFUSION DEVICES. DEVICES AND KITS FOR EXTRACORPOREAL PHOTOCHEMOTHERAPY OR ...	B030204
GMDN	In vitro ultraviolet blood irradiator	58075
EAN	UVA PIT System	42 60197191 004
Notified Body MDC	UVA PIT System	0483
Quality Management System according EN ISO 13485		
Council Directive/Medical Device Directive (MDD) 93/42/EEC and subsequent amendments verified by conformity valuation procedure. Annex II - excluding Section 4.		

18.2 PIT KIT

Ministry of Health, Italy	PIT KIT	54D1B01000; 54C1G05000
CND	HEMATOLOGY AND HEMOTRANSFUSION DEVICES. DEVICES AND KITS FOR EXTRACORPOREAL PHOTOCHEMOTHERAPY OR ...	B030204
GMDN	Off-line blood photoimmunotherapy bag	62222
EAN	PIT KIT	42 60197190 014
Notified Body IMQ	PIT KIT	0051
Quality Management System according EN ISO 13485		
Council Directive/Medical Device Directive (MDD) 93/42/EEC and subsequent amendments verified by conformity valuation procedure. Annex II - excluding Section 4.		

	Document name	Operator's Manual UVA PIT System CS		
	Document No.	Operators_manual_UVA_PIT_CS	Language	EN

19 Limited warranty

PIT Medical Systems warrants its device to conform to the specification contained herein and against defects in workmanship or material under normal use and service for a period of one year from the date of startup.

PIT Medical Systems' sole obligation is to repair or replace defective components sold under this warranty at its own discretion.

The equipment installation and setup are performed by PIT Medical Systems or an authorized representative.

This warranty does not apply to electrical fuses, lamps, filters, or other consumable items. PIT Medical Systems will not warrant parts found to be defective due to shipping damage, misuse, use of unauthorized parts or lack of maintenance as specified within the UVA PIT System Operator's Manual.

This warranty does not apply if the device is operated by an unauthorized or non-trained operator. Operators need to be authorized by PIT Medical Systems.

This warranty extends to the original purchaser of device from PIT Medical Systems and is not transferred upon resale.

Buyer's sole and exclusive remedy for breach of PIT Medical Systems warranty obligation shall be limited to the cost of equipment and in no event include consequential damages or liability for personal injury or property damages.

The foregoing warranty is in lieu of all other warranties, expressed or implemented and no other warranties exist including, without limitation, any warranty of merchantability or fitness for a particular purpose.

If a defect develops with your UVA PIT System, please contact PIT Medical Systems or your local representative.

Am Farrnbach 1
 90556 Cadolzburg, Deutschland
 Tel: +49-(0)9103-68269-80
 E-Mail: office@pitmed.com

20 Documentation of the Procedure

As an example, PIT Medical Systems GmbH provides this form for an internal documentation of an ECP.

User Protocol

Please complete the following document.

UVA PIT	System S/N:	<input type="text"/>	Software version:	<input type="text"/>
	PIT KIT Lot-No.:	<input type="text"/>		
8-MOP	Manufacturer:	<input type="text"/>	Batch No.:	<input type="text"/>
	Form and dosage:	<input type="text"/>		
Apheresis unit	Name of company and device:	<input type="text"/>		
Patient	ID / Name / Date of Birth:	<input type="text"/>		
Date / Time	Date:	<input type="text"/>	Start of procedure (time):	<input type="text"/>
	ECP procedure	Priming:	<input type="text"/>	Calculated energy:
	BC Volume:	<input type="text"/>	Lamp energy:	<input type="text"/>
	HCT:	<input type="text"/>	Chamber temperature:	<input type="text"/>
	Treatment done:	<input type="text"/>	Aborted:	<input type="text"/>
			End of procedure (time):	<input type="text"/>

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List of devices covered by this notification

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PIT Medical Systems GmbH | Am Farrnbach 1 | 90556 Cadolzburg
TO WHOM IT MAY CONCERN

24.07.2023

LETTER OF AUTHORIZATION

We, **PIT Medical Systems GmbH, Am Farrnbach 1, 90556 Cadolzburg, Germany**, do hereby authorize Medservice Sztandar Sp. k. (Limited Partnership)
ul. Bojkowska 59 C, 44 - 100 Gliwice, Polska
as our exclusive distributor for the products:

UVA PIT System©
DATA LOG
PIT KIT© (disposable Kits) and its variant PIT KIT/CS©

and technical service in the territory of Poland.

PIT Medical Systems GmbH has assumed all responsibilities for the UVA PIT system product made by Med Tech Solutions GmbH, Am Farrnbach 1, 90556 Cadolzburg.
Since 2018, the PIT KIT has been exclusively manufactured for PIT Medical Systems GmbH by HMC Pre Medical SPA.
Since 2018, the UVA PIT System is manufactured by PIT Medical Systems GmbH.

If you have any questions regarding this matter, please reach us at office@pitmed.com.

Yours Faithfully

Petra Weber
CEO