
# Datasheet AquaWTU

## Compact reverse osmosis system

### ISO standards

The AquaWTU is the ideal system for fulfilling the requirements of the ISO dialysis water quality standards.

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### Medical Device Class IIb according Medical Device Directive 93/42/EEC

#### Key features & functionality

- Dead-space free design
- Detailed error reporting
- Integrated water pretreatment
- Fine filter (10 µm)
- Alternating softener
- UV lamp
- Activated carbon filter: 5 µm
- AUTO START/STOP programming
- Soft start function
- Semi-automatic chemical disinfection
- Semi-automatic decalcification
- Initial self-testing of all safety-relevant actuators and sensors
- High-speed membrane flushing interval during rinse programme
- Internal leakage sensor

#### Options

- Dialysis water tank
- Remote control with status indicator LEDs
- 3/4" raw water pressure reducer
- Permeate ring main pressure reducer upstream of haemodialysis system
- Integrated 20 µm prefilter
- External volt-free contact for error indicator (staff call)

## Specifications

Haemodialysis devices	1–5 devices
Permeate capacity	<b>AquaWTU 125/250:</b> 125–250 L/h
Efficiency/yield	Up to 75 %
Permeate pressure	Approx. 2–3 bar
Dimensions in mm (h × w × d)	925 × 510 × 580 (incl. optional prefilter 720 mm)
Weight	<b>AquaWTU 125/250:</b> 110 kg
Concentrate pressure	max. 17 bar
Noise level	Noise level in SUPPLY mode: 49–50 dB (distance of 1 m) (depending on system capacity and features)

## Electrical supply

Different voltage types possible	– 110 V/60Hz – 120 V/60 Hz – 220 V/230 V; 50 Hz – 220 V/230 V; 60 Hz – 240 V/50 Hz
Electrical supply/operating current	Type <b>AquaWTU 125/250:</b> – 3.9 A nominal current/voltage type (220–240 V); – 7.5–8.2 A nominal current/voltage type (110–120 V)
Power consumption	~0.95 kVA
Radiated heat/loss	<b>AquaWTU 125/0.50 kW</b> <b>AquaWTU 250/0.65 kW</b>
Overcurrent protection	Overcurrent protection: overload circuit breaker B10/C10, depending on specified current RCCD 30 mA recommended
Socket	Standard (Schuko plug) US-NEMA, NEMA 5–15, 3 phases US-NEMA, 5–15, 3 phases, type SJT, hospital grade UK, British Standard BS 1363 Swiss, SEV 1011 Australia, AS 3112
Type of protection against electric shock	Protection class I
Applied parts classification	Type B
Degree of protection against ingress of liquids	Drip-proof
Leakage currents	According to EN 60601-1
Overvoltage category	II
Pollution severity	II
Material group	III b
Operating mode	Continuous operation (standby)

## Product water quality

Bacteria (CFU) and endotoxins (EU)	> 99 %
Total dissolved solids	> 96 %

### Water supply

Feed pressure	Dynamic 3.5–6 bar
Minimum inlet flow	<b>AquaWTU 125:</b> min. 600 L/h <b>AquaWTU 250:</b> min. 900 L/h
Permeate connection	– DN 10 permeate feed – DN 10 permeate return
Inlet connection	¾ inch external thread
Drain water connection	– DN 10 concentrate drain – DN 10 permeate drain – DN 10 softener drain
Maximum length of the ring main DN 10	The length of the ring main is restricted: max. 30 m

### Operating conditions

Water hardness	< 30 °dH
Iron	< 0.1 mg/L
Manganese	< 0.05 mg/L
Chloride	< 100 mg/L
Silicate	< 25 mg/L
Total chlorine	< 0.1 mg/L
Conductivity	< 2000 µS/cm
Total salt content	< 1500 mg/L
pH	6–8
Silt density index	< 5
Feed water temperature	min. 5 °C/max. 35 °C
Atmospheric pressure	Ambient pressure: 700–1150 hPa
Ambient temperature range	+5 °C to +40 °C
Relative humidity	Up to 80 % at 20 °C (non-condensing)

### External connection options

Tank control	Tank control option
Remote control	Optional remote control and optional volt-free contact alarm

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---

Type	Material	
Polymers	PP	Polypropylene
	PSU	Polysulfone
	PVC	Polyvinyl chloride
Rubber	EPDM	Ethylene propylene diene monomer
		Silicone
Metals	1.4571	Stainless steel
	1.4404	Stainless steel
	Ti	Titanium

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

Fresenius Medical Care

# Datasheet AquaUNO

## Single station reverse osmosis system

### ISO standards

The AquaUNO is the ideal system for fulfilling the requirements of the ISO dialysis water quality standards.

**ISO 23500-1** – Guidance for the preparation and quality management of fluids for haemodialysis and related therapies

**ISO 23500-2** – Water treatment equipment for haemodialysis applications and related therapies

**ISO 23500-3** – Water for haemodialysis and related therapies

Medical Device Class IIb according Medical Device Directive 93/42/EEC

### Key features & functionality

- Dead-space free design
- Detailed error reporting
- Semi-automatic chemical disinfection and decalcification
- Initial self-testing of all safety-relevant actuators and sensors

### Options

- AquaCART (alternate portable water pretreatment system)
- AquaSTOP (leakage detector)
- 3/4" raw water pressure reducer
- Pressure switch option for automatic START/STOP function
- Tubing holder
- Different ring main lengths: 2 x 10 m
- Portable softener (to reduce the raw water hardness below 1° dH)
- START/STOP of reverse osmosis with haemodialysis device (any brand), controlled by hydraulic sensors (if pressure switch option is installed)
- U-shaped permeate connection (different hydraulic connectors available)

## Specifications

Haemodialysis devices	1 device
Permeate capacity with < 25°dH	Approx. 50 L/h at 10°C
Permeate output with < 1°dH	– Approx. 50 L/h at 5°C – Approx. 60 L/h at 10°C – Approx. 70 L/h at 15°C
Permeate pressure	Approx. 2–3 bar
Efficiency	Up to 50 % by 100% consumption
Dimensions in mm (h x w x d)	455 x 515 x 200
Weight (empty/filled)	(Approx. 30 kg/approx. 32 kg)
Noise level	Noise level in SUPPLY mode: 44–46 dB (A); (distance of 1 m); (depending on system capacity and features)

## Electrical supply

	Voltage	Frequency	Max. power consumption	Max. operating current
<b>Different voltage types possible</b>				
	220 V	50/60 Hz	400 VA	10 A
	230 V	50/60 Hz	400 VA	10 A
	240 V	50/60 Hz	400 VA	10 A
	100 V	50/60 Hz	475 VA	10 A
	110 V	50/60 Hz	475 VA	10 A
	115 V	50/60 Hz	475 VA	10 A
	120 V	50/60 Hz	475 VA	10 A
<b>Overcurrent protection for current</b>				
Protection/220 V to 240 V/current 6 A specification 9–10 A	Overload circuit breaker B6/C6 (depending on specified current)			
Protection/100 V to 120 V/current 10 A	Overload circuit breaker B10/C10 (depending on power supply system)			
<b>Residual current circuit breaker</b>	RCCD 30 mA recommended			
<b>Socket</b>	Standard (Schuko plug) US-NEMA, NEMA 5–15, 3 phases US-NEMA, 5–15, 3 phases, Typ SJT, hospital grade UK, British Standard BS 1363 Swiss, SEV 1011 Australia, AS 3112			
<b>Type of protection against electric shock</b>	Protection class I			
<b>Applied parts classification</b>	Type B			
<b>Degree of protection against ingress of liquids</b>	Drip-proof			
<b>Leakage currents</b>	According to EN 60601-1			
<b>Overvoltage category</b>	II			
<b>Pollution severity</b>	II			
<b>Material group</b>	III b			
<b>Operating mode</b>	Continuous operation (standby)			

### Product water quality

Bacteria (CFU) and endotoxins (EU)	> 99 %
Total dissolved solids	> 96 %

### Water supply

Feed pressure	Max. 6 bar
Minimum inlet flow	– 1.5 bar dynamic at 150 L/h feed and max. 50 L/h permeate consumption – 2.0 bar dynamic at 150 L/h feed and more than 50 L/h permeate consumption
Permeate connection	Nipple (standard) or coupling with DN6 tubing connection
Inlet connection	G ¼" union nut
Drain water connection	DN6 tubing connection
Permeate ring main	Standard: 2 × 2 m Option: 2 × 10 m

### Operating conditions

Water hardness	< 25 °dH
Iron	< 0.1 mg/L
Manganese	< 0.05 mg/L
Chloride	< 100 mg/L
Silicate	< 25 mg/L
Total chlorine	< 0.1 mg/L
Conductivity	< 2000 µS/cm
Total salt content	< 1500 mg/L
pH	6–8
Silt density index	< 3
Inlet volume	Min. 120 L/h
Feed water temperature	Min. 5 °C/max. 35 °C
Atmospheric pressure	Ambient pressure: 700–1150 hPa
Ambient temperature range	+5 °C to +35 °C
Relative humidity	Up to 80 % at 20 °C (non-condensing)

### External connection options

Interface:	Interface for connection to a 4008S/5008/5008S/6008 haemodialysis device
------------	--

### Transport and storage conditions

Storage temperature range	+5 °C to +40 °C (protect from frost)
Storage time	Storage time of preserved system: maximum 12 months
Atmospheric pressure	Ambient pressure: 500–1150 hPa
Relative humidity	Up to 80 % at 20 °C (non-condensing)

## Materials in contact with dialysis water

---

Type	Material	
Polymers	PP	Polypropylene
	PA	Polyamide
	PSU	Polysulfone
	PPO	Polyphenylene oxide
	FEP	Fluorinated ethylene propylene
	PVDF	Polyvinylidene fluoride
	PBT	Polybutylene terephthalate
Rubber	EPDM	Ethylene propylene diene monomer
		Silicone
Metals	1.4571	Stainless steel
	1.4404	Stainless steel
	Ti	Titanium
Ceramics	Al <sub>2</sub> O <sub>3</sub>	Ceramic

**FRESENIUS  
MEDICAL CARE**

# Datasheet AquaC CART Basic

Transporting vehicle with prefiltration

## ISO standards

The **AquaC CART Basic** is the ideal system for fulfilling the requirements of the ISO dialysis water quality standards.

**ISO 23500-2** – Water treatment equipment for haemodialysis applications and related therapies

Accessory according to Medical device Directive 93/42/EEC, covered by Declaration of Conformity of AquaC UNO H and AquaUNO

## AquaC CART Basic

The **AquaC CART Basic** water pretreatment system is a portable unit for the pretreatment of drinking water for the operation of an AquaC UNO H or AquaUNO type reverse osmosis system.

## Key features & functionality

- Easy to move
- Wheel-stop function
- Service position position for reverse osmosis
- Convenient position for use
- Space-saving design
- Easy-to-clean materials
- Inlet with raw water shut-off valve and manometer
- Particle filter 20" / 50  $\mu\text{m}$
- Fine filter 20" / 20  $\mu\text{m}$
- Activated carbon filter 20" / 5  $\mu\text{m}$
- Fine filter 20" / 1  $\mu\text{m}$
- Soft water sampling with manometer

## Options

- Individual filter equipment available

## Specifications

Dimensions in mm (h × w × d)	850 × 510 × 760
<b>Weight</b>	
Empty	~ 29 kg
Filled	~ 34 kg
<b>20" particle filter:</b>	
Filtration	50 µm
<b>20" fine filter</b>	
Filtration	20 µm
<b>20" activated carbon filter</b>	
Filtration	5 µm
<b>20" fine filter</b>	
Filtration	1 µm

## Water supply

<b>Water inlet pressure (dynamic)</b>	
min.	1 bar
max.	6 bar
<b>Raw water connection</b>	Rp ¾"
<b>Internal connections</b>	Rp ¾"
<b>Reverse osmosis connections</b>	Tubing connection

## Transport and storage conditions

<b>Storage temperature range</b>	+5 °C to +40 °C (protect from frost)
<b>Transport</b>	Keep upright, handle with care!
<b>Storage</b>	Store upright!
<b>Relative humidity</b>	Up to 80 % at 20 °C (non-condensing)

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## 1.10 AquaC CART Basic

POS	Description
A	Particle filter (50 $\mu\text{m}$ )
B	20" particle filter (20 $\mu\text{m}$ )
C	20" activated carbon filter (5 $\mu\text{m}$ )
D	20" particle filter (1 $\mu\text{m}$ )

### 1.10.1 Part number – AquaC CART Basic

Part number	Description
6343491	AquaC CART Basic Cart with water pretreatment system (supplied without reverse osmosis)

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