stryker®

Operations Manual


stryker[®]

Operations Manual

Symbols

English
EN

| | |
|---|---|
|  | Refer to instruction manual/booklet |
|  | Operating instructions / Consult instructions for use |
|  | General warning |
|  | Caution |
|  | Warning; crushing of hands |
|  | Warning; crushing of feet |
|  | Do not insert lift pole |
|  | Headboard and footboard orientation |
|  | Do not store items under the bed. |
|  | Gatch positioning |
|  | Catalogue number |
|  | Serial number |
|  | CE mark |
|  | Manufacturer |
|  | Mass of equipment with safe working load |
|  | Safe working load |
|  | Maximum patient weight |
|  | Adult patient |

Symbols

English
EN

| | |
|---|---|
|  | Direct current |
|  | Alternating current |
|  | Dangerous voltage |
|  | Unit provides terminal for connection of a potential equalization conductor. The potential equalization conductor provides direct connection between the unit and potential equalization busbar of the electrical installation. |
|  | Protective Earth terminal |
| IPX4 | Protection from liquid splash |
|  | Type B applied part |
|  | In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately. Contact your local distributor for disposal information. See the maintenance manual for servicing instructions and recyclable components. |

Table of Contents

English
EN

| | |
|--|------|
| Warning/Caution/Note Definition | 1-3 |
| Summary of safety precautions | 1-4 |
| Pinch points | 1-7 |
| Introduction | 1-8 |
| Product description | 1-8 |
| Indications for use | 1-8 |
| Expected service life | 1-9 |
| Contraindications | 1-9 |
| Warranty | 1-9 |
| Specifications | 1-9 |
| Product illustration | 1-12 |
| Applied parts | 1-13 |
| Contact information | 1-13 |
| Serial number location | 1-14 |
| Setup | 1-15 |
| Operation | 1-16 |
| Plugging or unplugging the battery cable | 1-16 |
| Plugging or unplugging the product | 1-17 |
| Charging the battery | 1-17 |
| Storing the battery long term | 1-18 |
| Storing the power cord | 1-18 |
| Transporting the product | 1-19 |
| Applying or releasing the brakes | 1-19 |
| Applying or releasing steer lock | 1-21 |
| Applying or releasing the fifth wheel (option) | 1-22 |
| Activating and resetting the CPR release | 1-23 |
| Removing or replacing the headboard | 1-23 |
| Removing or replacing the footboard | 1-24 |
| Raising or lowering the lower leg section | 1-25 |
| Raising or lowering the siderails | 1-26 |
| Nurse control panel (outside siderail) (option) | 1-27 |
| Patient control panel (inside siderail) (option) | 1-28 |
| Patient control pendant (option) | 1-30 |
| Nurse control pendant | 1-31 |
| Extending the bed extender (option) | 1-32 |
| Installing the bed extender bolster mattress | 1-33 |
| Extending or retracting the linen tray (option) | 1-34 |
| Inserting or removing a cassette from the X-ray cassette holder (option) | 1-34 |
| Accessories | 1-36 |
| Installing the IV pole | 1-36 |
| Adjusting the IV pole | 1-37 |
| Installing the lifting pole | 1-37 |
| Installing the lifting pole handle | 1-38 |

Table of Contents

English
EN

Installing the oxygen bottle holder 1-39

Installing the Foley bag basket 1-39

Cleaning..... 1-41

 Preparing the product for cleaning 1-41

 Cleaning..... 1-41

 Cleaning the siderails 1-42

Disinfecting 1-43

Preventive maintenance 1-44

EMC information 1-46

Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

WARNING

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

CAUTION

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

Note: Provides special information to make maintenance easier or important instructions clearer.

Summary of safety precautions

English
EN

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

WARNING

- Only use the input voltage and frequency as rated on the product.
 - Always allow the product to reach room temperature before you begin setup or test functional operations, to prevent permanent product damage.
 - Do not use this product if it has identifiable failures, defects, malfunctions, or damage.
 - Do not use this product under any conditions if use would cause injury to the operator or patient.
 - Only operate the product when all operators are clear of the mechanisms.
 - Do not attach the power cord to any parts of the product.
 - Always unplug the power cord and call maintenance if unanticipated motion occurs.
 - Do not store items under the bed.
 - Do not use the bed without the litter covers.
 - To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
 - Always store the power cord to avoid the risk of entanglement, damage to the power cord, or potential shock hazards. If the power cord is damaged, remove the product from service and contact the appropriate maintenance personnel.
 - Always allow enough clearance between the head end of the product and the adjacent wall, so you can unplug the power cord from the wall outlet during an emergency.
 - Always store the power cord before transporting the product.
 - Always disconnect the power cord from the wall outlet if you detect overheating of the battery, control cables, or pendants. Do not use the product again until it has been inspected, serviced, and confirmed to work as intended by authorized maintenance personnel.
 - Always replace the battery after it surpasses its expected service life.
 - Do not open a dead battery.
 - Do not throw the battery into a fire.
 - Do not spill liquid onto the battery or submerge the battery in liquid.
 - Always unplug the battery cable from the control box before you store the product for a long period of time.
 - Always lock the siderails in the full up position with the sleep surface horizontal in the lowest position when transporting a patient.
 - Always keep limbs, hands, fingers, and other body parts clear of mechanisms and gaps.
 - Always make sure that there are no obstacles near the product. Injury to the patient, operator, bystanders or damage to the frame or surrounding equipment could occur if you collide with an obstacle.
 - Do not attempt to move the product laterally. This may cause the product to tip.
 - Do not move the product after you apply the brakes.
 - Always apply the brakes when a patient is getting into or out of the product to avoid instability.
 - Always apply the brakes when the patient is unattended.
 - Always make sure that all persons and equipment are away from the area below and around the backrest before you activate the CPR release. The CPR release is for emergency use only.
 - Always properly orient the headboard when replacing the headboard to avoid entrapment.
 - Always properly orient the footboard when replacing the footboard to avoid entrapment.
 - Always make sure that all persons and equipment are away from the area below and around the leg rest before you lower the lower leg section.
 - Always lock the siderails unless a patient's condition requires extra safety measures.
 - Always lock the siderails in the full up position when the patient is unattended.
-

Summary of safety precautions

English
EN

WARNING (CONTINUED)

- Do not use siderails as restraint devices to keep the patient from exiting the product. The operator must determine the degree of restraint necessary to make sure that the patient is safe.
 - Do not sit on the siderails.
 - Always lock bed motion controls when the patient is unattended.
 - Never store the nurse control panel within the patient's reach.
 - Do not sit on the bed extender. This may cause the product to tip.
 - Always lock the bed extender before you place weight on the bed extender.
 - Always retract the linen tray (option) before you put the product in motion.
 - Always retract the linen tray (option) when not in use.
 - Do not use the product for X-ray procedures if it does not have the radiolucent backrest (option).
 - Do not use accessories to support patient limbs or other body parts.
 - Do not clean, service, or perform maintenance while the product is in use.
 - Always power off and unplug the power cord before cleaning, servicing, or performing maintenance.
 - Always power off the product and unplug the power cord from the wall outlet when large spills occur near the circuit boards, cables, and motors. Remove the patient from the product, clean up the fluid, and have service personnel inspect the product. Fluids can cause unpredictable operation and decreased functionality of any electrical product. Do not return the product to service until it is completely dry and has been thoroughly tested for safe operation.
 - Do not spray cleaners directly onto the battery, control boxes, actuators, cables, or other electric equipment.
 - Do not use abrasive powder, steel wool, or similar materials that may damage the product surface.
 - Do not use **Virex**® TB for product disinfecting.
 - Do not use acid-based chemicals or flammable chemicals, such as gasoline, diesel, or acetone for cleaning purposes.
 - The cleaners and disinfectants must not be highly alkaline or acidic (pH value 6-8).
 - Do not directly spray or saturate the siderail control panel, patient control pendant, or nurse control pendant with cleaners.
 - Do not use sharp objects to clean the siderail control panel.
-

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
 - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
 - To minimize the risk of any electromagnetic interference, the product design follows the standard IEC 60601-1-2. To avoid problems, use the bed in accordance with the EMC/EMI requirements in the EMC section of this operations manual.
 - Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in.) to any part of the SV2, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
 - Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
-

Summary of safety precautions

English
EN

CAUTION (CONTINUED)

- Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Do not place objects into any gaps of the product.
- Always plug the product into a wall outlet (regulated AC-power source) when not in use to maintain a sufficient battery charge and to maximize product performance while operating on battery power.
- Always replace batteries that have corrosion at the terminals, display cracking, have expanded or bulging sides, or no longer can maintain a full charge.
- Only use authorized batteries when replacing the batteries. Use of unauthorized batteries may lead to unpredictable system performance.
- Do not place or store heavy objects on the product.
- Do not squeeze or pinch the power cord in the bed frame.
- Do not use the siderails as a push or pull device. Always move the product using the integrated handles in the headboard and footboard.
- Always remove the patient lift pole before transporting the product.
- Do not use the IV pole as a push or pull device.
- Always make sure that the IV pole is at a low height during transport.
- Always apply the brake to prevent unintended movement.
- Do not apply the brake pedal to stop a moving product.
- Always hang the nurse control pendant onto a foot end siderail or the linen tray (option) before removing the footboard.
- Do not raise the lower leg section while the bed extender is in use. This is to avoid the situation where the product does not support the lower legs of a taller patient.
- Always place the patient control pendant safely on the mattress while the pendant is in use.
- Always hang the patient control pendant onto the siderail when the pendant is not in use.
- Do not squeeze or pinch the pendant cord in the bed frame.
- Always place the nurse control pendant onto the footboard.
- Do not remove the footboard after you extend the bed extender.
- The safe working load of the linen tray is 15 kg.
- Only use authorized accessories for this product. The use of unauthorized accessories may result in product damage or injury to the operator or patient. Stryker is not responsible for any damage or injury that may result from the misuse of the product or the use of unauthorized accessories.
- Always make sure that accessories are locked into position.
- Do not allow accessories to interfere with mechanical or electric mechanisms of the product.
- Do not use the lifting pole as a push or pull device.
- Do not use the oxygen bottle holder as a push or pull device.
- Always turn the oxygen bottle holder in toward the bed before transporting a patient.
- Do not strike the oxygen bottle holder while transporting a patient.
- The safe working load of each Foley hook is 2 kg.
- Do not steam clean, pressure wash, ultrasonically clean, or immerse any part of the product in water. Exposure to water may damage the internal electric parts. These methods of cleaning are not recommended and may void this product's warranty.

Summary of safety precautions

English
EN

⚠ CAUTION (CONTINUED)

- Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning. Some cleaning products are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, you may leave a corrosive residue on the surface of the product that could cause premature corrosion of critical components. Failure to follow these cleaning instructions may void your warranty.

Pinch points


Figure 1-1: SV2 pinch points

Introduction

English
EN

This manual assists you with the operation or maintenance of your Stryker product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Notes

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- Stryker continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact Stryker Customer Service or Technical Support at + 90 (352) 321 43 00 (pbx).

Product description

SV2 is an AC-powered bed with a battery backup system. **SV2** is for patients to receive treatment for extended periods of time at hospitals and care centers. **SV2** has four electrical actuators that allow it to adjust to numerous positions, including CPR, Trendelenburg, Reverse Trendelenburg, and chair positions. **SV2** is equipped with retractable siderails, removable headboard and footboard, and options and accessories that assist with the care of the patient.

SV2 is an electromechanical MedSurg and ICU bed with DC-powered actuators and controls to adjust the patient sleep surface. The patient sleeping surface consists of four sections: the backrest, seat, upper leg section, and lower leg sections. Siderails are split, with two siderails on the head end, and two siderails on the foot end. The siderails secure in the full up position. When unlatched, siderails open outside and move to the lowest position.

You can actuate electromechanical functions with the siderail control panel, patient control pendant, and nurse control pendant. The control box consists of logic controls and a power supply that power and control signals to all four actuators via a distribution box. The siderail control panels, patient control pendant, and nurse control pendant controls are also controlled by the control box via distribution box.

The bed is equipped with two pairs of actuators (four actuators total). The first pair below the litter surface control the backrest down and up functions, and upper leg down and up functions. The second pair of actuators below the undercarriage control the litter down and up functions, Trendelenburg, and Reverse Trendelenburg.

Additional bed mechanisms allow for manual CPR, knee gatch motion, and bed length extension. The bed is also equipped with brake and steer control for the casters. Casters help in emergency or non-emergency intra-hospital transport of a patient on the bed.

Indications for use

SV2 is for use by human adult patients in a MedSurg and ICU setting requiring the support of a hospital bed. Use this product with a patient sleep surface.

Operators for the bed include healthcare professionals (such as nurses, nurse aides, and medical doctors), service or maintenance personnel, patients, and bystanders who can use bed motion functions.

SV2 is for use in medical, surgical, and critical care healthcare environments, including hospitals, institutions and clinics.

The **SV2** bed frame, litter mounted accessories, and mattresses can come in contact with human skin.

Indications for use (Continued)

The **SV2** bed frame is not intended to be used with an oxygen tent, in the presence of flammable anesthetics, or to support more than one individual at a time.

Expected service life

SV2 has a ten year expected service life under normal use conditions and with appropriate periodic maintenance.

The battery has a one year expected service life under normal use conditions.

The casters have a two year expected service life under normal use conditions.

The optional fifth wheel has a two year expected service life under normal use conditions.

Contraindications

None known.

Warranty

Warranty may vary by country. Contact your local Stryker Medical representative for additional information.

Specifications

| | | |
|---|---|-----------------------|
|  | Safe working load | 250 kg |
| | Note: Safe working load indicates the sum of the patient, mattress, and accessory weight | |
|  | Maximum patient weight | 215 kg |
| Product weight | | 160 kg |
| Overall product size | Length | 2200 mm (±10 mm) |
| | Length (with bed extender - option) | 2510 mm (±10 mm) |
| | Width | 990 mm (±10 mm) |
| Product height (without mattress) | Low | 375 mm (+15 / -25 mm) |
| | High | 755 mm (±10 mm) |
| Under product clearance | | 150 mm |
| Caster size (single and optional dual-casters) | | 150 mm |
| Product angle indicator | | 0° - 15° |
| Backrest angle indicator | | 0° - 90° |
| Backrest angle | | 0° - 60° |

Introduction







English
EN

Specifications (Continued)

| | |
|-------------------------------------|--|
| Trendelenburg/Reverse Trendelenburg | 0° - 12° |
| Gatch angle | 0° - 30° |
| Electrical requirements | |
| Battery | 24 VDC, 10 amps, Model BA1812 |
| Control box | 100-240 VAC, 50 Hz - 60 Hz nominal, 5 amps |
| Electrical classification | Class 1 when product is plugged into mains power Internally powered when the product is unplugged |
| Duty cycle | 2 mins of actuation and 18 mins idle |

Class I Equipment: Equipment that protects against electrical shock and does not solely rely on basic insulation, but which includes an additional safety precaution that is provided for the connection of the equipment to the protective earth conductor in the fixed wiring of the installation that accessible metal parts cannot become live in the event of a failure of basic insulation.

| | |
|------------------------------|---------------------------|
| Compatible mattresses | |
| 7002-2-012 | 2000 mm x 860 mm x 120 mm |
| 7002-2-014 | 2000 mm x 860 mm x 140 mm |
| 7002-2-512 | 2000 mm x 860 mm x 120 mm |
| 7002-2-514 | 2000 mm x 860 mm x 140 mm |
| 7002-5-012 | 2000 mm x 860 mm x 120 mm |
| 7002-5-014 | 2000 mm x 860 mm x 140 mm |
| 7002-5-512 | 2000 mm x 860 mm x 120 mm |
| 7002-5-514 | 2000 mm x 860 mm x 140 mm |
| 7002-5-712 | 2000 mm x 860 mm x 120 mm |

| Environmental conditions | Operation | Storage and transportation |
|--------------------------|---|---|
| Temperature |  (5 °C) to (38 °C) |  (-10 °C) to (50 °C) |
| Relative humidity |  20% to 90% |  20% to 90% |
| Atmospheric pressure |  800 hPa to 1060 hPa |  800 hPa to 1060 hPa |

Specifications (Continued)

Specifications listed are approximate and may vary slightly from product to product or by power supply fluctuations.

Stryker reserves the right to change specifications without notice.

| Standards applied | |
|--|--|
| IEC 60601-1:2012 | Medical electrical equipment - Part 1: General requirements for basic safety and essential performance |
| IEC 60601-1-2:2014 | Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests |
| IEC 60601-2-52:2009 + A1:2015 | Medical electrical equipment - Part 2-52: Particular requirements for the basic safety and essential performance of medical beds |
| IEC 60601-2-54:2009 + A1:2015* *Only applicable when the product is equipped with the radiolucent backrest option | Medical electrical equipment - Part 2-54: Particular requirements for the basic safety and essential performance of X-ray equipment for radiography and radioscopy |

WARNING

Only use the input voltage and frequency as rated on the product.


CAUTION

- To minimize the risk of any electromagnetic interference, the product design follows the standard IEC 60601-1-2. To avoid problems, use the bed in accordance with the EMC/EMI requirements in the EMC section of this operations manual.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in.) to any part of the SV2, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Introduction

English
EN

Product illustration

| | | | |
|---|---------------------------------------|---|--|
| A | Accessory sleeve | L | Linen tray (option) |
| B | Backrest | M | Lower leg section |
| C | Backrest indicator | N | Roller bumper |
| D | Bed angle indicator | O | Seat section |
| E | Bed extender (option) | P | Siderails |
| F | Brake/steer pedal | Q | Siderail control panel (Inside siderail) (option) |
| G | Casters (Dual-wheel casters optional) | R | Siderail control panel (Outside siderail) (option) |
| H | CPR release | S | Siderail latch |
| I | Foley hooks | T | Upper leg section |
| J | Footboard | U | Nurse control pendant (option) |
| K | Headboard | V | Patient control pendant (option) |

Applied parts

Figure 1-2: Type B applied parts

Contact information


Contact your Stryker Customer Service:

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Have the serial number (A) of your product available. Include the serial number in all written communication.

Serial number location

WARNING

- Always allow the product to reach room temperature before you begin setup or test functional operations, to prevent permanent product damage.
 - Do not use this product if it has identifiable failures, defects, malfunctions, or damage.
 - Do not use this product under any conditions if use would cause injury to the operator or patient.
 - Only operate the product when all operators are clear of the mechanisms.
 - To avoid the risk of electric shock, this equipment must be connected to a supply mains with protective earth.
 - Always store the power cord to avoid the risk of entanglement, damage to the power cord, or potential shock hazards. If the power cord is damaged, remove the product from service and contact the appropriate maintenance personnel.
 - Always allow enough clearance between the head end of the product and the adjacent wall, so you can unplug the power cord from the wall outlet during an emergency.
 - Always store the power cord before transporting the product.
 - Do not attach the power cord to any parts of the product.
 - Always unplug the power cord and call maintenance if unanticipated motion occurs.
 - Do not store items under the bed.
 - Do not use the bed without the litter covers.
-

CAUTION

Do not place objects into any gaps of the product.

Note: Product has suitable mains to isolate its circuits electrically from the supply mains on all poles at the same time.

Before placing the product into service, make sure that these components work:

1. Visually inspect the product for any signs of shipping damage.
2. Make sure that the product and all components and accessories have arrived.
3. Depress the brake pedal and make sure that the brake, steer, and neutral positions operate.
4. Raise and lower the siderails to make sure that they move, stow, and lock securely in the full up position.
5. Plug the battery cable into the control box ([Plugging or unplugging the battery cable on page 1-16](#)).
6. Plug the power cord into a protective earthed outlet ([Plugging or unplugging the product on page 1-17](#)).
7. Press each button on the siderail control panel, nurse control pendant, and patient control pendant (option) to make sure that each function operates ([Nurse control pendant on page 1-31](#)).
8. Make sure that the battery is fully charged.
9. Make sure that the CPR release handle operates.
10. Make sure that optional accessories are installed and operate as described in this manual.

Plugging or unplugging the battery cable

To plug the battery cable into the control box:

1. Remove the upper leg rest cover (A) ([Figure 1-3 on page 1-16](#)).
2. Locate the control box ([Figure 1-4 on page 1-16](#)).
3. Connect the battery cable to the control box.
4. Press the battery cable lock to lock the battery cable into the control box (A) ([Figure 1-4 on page 1-16](#)).


Figure 1-3: Removing the upper leg rest cover

To unplug the battery cable from the control box:

1. Remove the upper leg rest cover (A) ([Figure 1-3 on page 1-16](#)).
2. Locate the control box ([Figure 1-4 on page 1-16](#)).
3. Unlock the battery control cable to the control box (A) ([Figure 1-4 on page 1-16](#)).
4. Disconnect the battery from the control box.
5. Using tape, secure the battery cable lock to the litter frame ([Figure 1-5 on page 1-16](#)).


Figure 1-4: Locking or unlocking the battery cable

Figure 1-5: Unplugging the battery cable from the control box

Plugging or unplugging the product

WARNING

- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- Always store the power cord to avoid the risk of entanglement, damage to the power cord, or potential shock hazards. If the power cord is damaged, remove the product from service and contact the appropriate maintenance personnel.
- Always allow enough clearance between the head end of the product and the adjacent wall, so you can unplug the power cord from the wall outlet during an emergency.
- Always store the power cord before transporting the product.

To plug in the product, plug the power cord into a protective earthed outlet.

To unplug the product, grasp the mold near the outlet and pull the cord in a direction parallel to the floor (not at an angle).

Charging the battery

WARNING

- Always disconnect the power cord from the wall outlet if you detect overheating of the battery, control cables, or pendants. Do not use the product again until it has been inspected, serviced, and confirmed to work as intended by authorized maintenance personnel.
- Always replace the battery after it surpasses its expected service life.
- Do not open a dead battery.
- Do not throw the battery into a fire.
- Do not spill liquid onto the battery or submerge the battery in liquid.

CAUTION

- Always plug the product into a wall outlet (regulated AC-power source) when not in use to maintain a sufficient battery charge and to maximize product performance while operating on battery power.
- Always replace batteries that have corrosion at the terminals, display cracking, have expanded or bulging sides, or no longer can maintain a full charge.
- Only use authorized batteries when replacing the batteries. Use of unauthorized batteries may lead to unpredictable system performance.

SV2 is equipped with a battery backup system that charges when the product is plugged into a wall outlet. The battery backup system allows the operator to use the product when the product is unplugged, during a power failure, or while transporting a patient. Battery backup functionality activates when you unplug the product.

Always check battery backup function according to the preventive maintenance checklist (see [Preventive maintenance on page 1-44](#)). Always replace the battery if it does not perform as intended during preventive maintenance.

Charging the battery (Continued)

To charge the battery, connect the product to a wall outlet. The battery has a full charge within 10 to 12 hours.

Storing the battery long term

WARNING

Always unplug the battery cable from the control box before you store the product for a long period of time.

CAUTION

Do not place or store heavy objects on the product.

Store the battery per the environmental conditions listed in the specifications section (see [Specifications on page 1-9](#)).

To store the battery:

1. See [Plugging or unplugging the product on page 1-17](#).
2. See [Plugging or unplugging the battery cable on page 1-16](#).

Storing the power cord

WARNING

- Always store the power cord before transporting the product.
- Always store the power cord to avoid the risk of entanglement, damage to the power cord, or potential shock hazards. If the power cord is damaged, remove the product from service and contact the appropriate maintenance personnel.

CAUTION

Do not squeeze or pinch the power cord in the bed frame.

To store the power cord, wrap the power cord around the cord wrap (A) underneath the head end of the product ([Figure 1-6 on page 1-18](#)).

Figure 1-6: Storing the power cord

Transporting the product

WARNING

- Always store the power cord before transporting the product.
- Always lock the siderails in the full up position with the sleep surface horizontal in the lowest position when transporting a patient.
- Always keep limbs, hands, fingers, and other body parts clear of mechanisms and gaps.
- Always make sure that there are no obstacles near the product. Injury to the patient, operator, bystanders or damage to the frame or surrounding equipment could occur if you collide with an obstacle.
- Do not attempt to move the product laterally. This may cause the product to tip.

CAUTION

- Do not use the siderails as a push or pull device. Always move the product using the integrated handles in the headboard and footboard.
- Always remove the patient lift pole before transporting the product.
- Do not use the IV pole as a push or pull device.
- Always make sure that the IV pole is at a low height during transport.

To transport the product:

1. Lock the siderail control panel and patient control pendant functions (see [Nurse control pendant on page 1-31](#)).
2. Unplug the power cord from the wall outlet.
3. See [Storing the power cord on page 1-18](#).
4. Store the control pendants.
5. Retract the linen tray (see [Extending or retracting the linen tray \(option\) on page 1-34](#)).
6. Lower the IV pole.
7. Turn in the oxygen bottle holder toward the product.
8. Raise and lock the siderails in the full up position (see [Raising or lowering the siderails on page 1-26](#)).
9. Release the brakes (see [Applying or releasing the brakes on page 1-19](#)).
10. Push the product from the headboard or footboard.

Applying or releasing the brakes

WARNING

- Do not move the product after you apply the brakes.
- Always apply the brakes when a patient is getting into or out of the product to avoid instability.
- Always apply the brakes when the patient is unattended.

CAUTION

- Always apply the brake to prevent unintended movement.
- Do not apply the brake pedal to stop a moving product.

Applying or releasing the brakes (Continued)

The brake pedals are located on each caster.

To apply the brakes, depress the red pedal. The brake pedal locks all four casters. This holds the product in place ([Figure 1-7 on page 1-20](#)).

Figure 1-7: Applying the brakes

To release the brakes, depress the green pedal until the pedal is in the neutral position ([Figure 1-8 on page 1-20](#)). This releases all four casters and allows you to freely move the product.

Figure 1-8: Releasing the brakes

Applying or releasing steer lock

The steer pedals are located on each caster.

To apply the steer caster, depress the green pedal ([Figure 1-9 on page 1-21](#)). This allows you to maneuver the product forward and backward in a straight path.

Figure 1-9: Applying the steer pedal

To release the steer caster, depress the red pedal until the pedal is in the neutral position ([Figure 1-10 on page 1-21](#)). This releases the right side caster on the foot end and allows you to freely move the product.

Figure 1-10: Releasing the steer pedal

Applying or releasing the fifth wheel (option)

The steer pedals are located on each caster.

To apply the fifth wheel, depress the green pedal ([Figure 1-11 on page 1-22](#)). This lowers the fifth wheel and allows you to maneuver the product forward and backward in a straight path.

Figure 1-11: Applying the fifth wheel

To release the fifth wheel, depress the red pedal until the pedal is in the neutral position ([Figure 1-12 on page 1-22](#)). This retracts the fifth wheel and allows you to freely move the product forward, backward, and from side to side.

Figure 1-12: Releasing the fifth wheel

Activating and resetting the CPR release

WARNING

Always make sure that all persons and equipment are away from the area below and around the backrest before you activate the CPR release. The CPR release is for emergency use only.

If the backrest is raised and you need quick access to the patient, you can position the product to 0° by activating the CPR release.

You can find the instant CPR release lever at the head end section on both the left and right sides of the backrest.

To activate the CPR release:

1. Grasp either lever (A) and pull outward (Figure 1-13 on page 1-23).
2. Guide the backrest down to the flat position.

To reset the backrest motor after activating the CPR release, press the backrest down button on the control panels or the CPR button on the nurse control panel.

Figure 1-13: CPR release

Removing or replacing the headboard

WARNING

Always properly orient the headboard when replacing the headboard to avoid entrapment.

You can remove the headboard for patient accessibility and cleaning.

To remove the headboard:

1. Hang the patient control pendant onto the head end siderail.
2. Grasp the handles and lift the headboard straight up and off the product (Figure 1-15 on page 1-24).

To replace the headboard:

1. Align the curved headboard corners with the foot end of the bed (Figure 1-14 on page 1-24).
2. Align the headboard pegs (A) with the plastic sleeves (B) at the head end of the product (Figure 1-15 on page 1-24).

Removing or replacing the headboard (Continued)

3. Lower the headboard until it seats into the plastic sleeves (B) ([Figure 1-15 on page 1-24](#)).

Figure 1-14: Headboard orientation


Figure 1-15: Removing or replacing the headboard

Removing or replacing the footboard

WARNING

Always properly orient the footboard when replacing the footboard to avoid entrapment.

CAUTION

Always hang the nurse control pendant onto a foot end siderail or the linen tray (option) before removing the footboard.

You can remove the footboard for patient accessibility and cleaning.

To remove the footboard:

1. Hang the nurse control pendant onto a foot end siderail or the linen tray (optional).
2. Grasp the handles and lift the footboard straight up and off the product ([Figure 1-17 on page 1-25](#)).

To replace the footboard:

1. Align the curved footboard corners with the head end of the bed ([Figure 1-16 on page 1-25](#)).
2. Align the footboard pegs with the plastic sleeves at the foot end of the product ([Figure 1-17 on page 1-25](#)).
3. Lower the footboard until it seats into the plastic sleeves ([Figure 1-17 on page 1-25](#)).

Removing or replacing the footboard (Continued)

Figure 1-16: Footboard orientation

Figure 1-17: Removing or replacing the footboard

Raising or lowering the lower leg section

WARNING

Always make sure that all persons and equipment are away from the area below and around the leg rest before you lower the lower leg section.

CAUTION

Do not raise the lower leg section while the bed extender is in use. This is to avoid the situation where the product does not support the lower legs of a taller patient.

You can raise or lower the lower leg rest manually.

To raise the lower leg section:

1. Grasp the lower leg section with both hands.
2. Raise the lower leg section to the desired height.
3. Release the lower leg section to lock the section into place.

To lower the lower leg section:

1. Grasp the lower leg section with both hands.
2. Raise the lower leg section all the way up to unlock the lower leg section.
3. Guide the lower leg section back down onto the litter.

Raising or lowering the siderails

WARNING

- Always lock the siderails in the full up position with the sleep surface horizontal in the lowest position when transporting a patient.
- Always keep limbs, hands, fingers, and other body parts clear of mechanisms and gaps.
- Always lock the siderails unless a patient's condition requires extra safety measures.
- Always lock the siderails in the full up position when the patient is unattended.
- Do not use siderails as restraint devices to keep the patient from exiting the product. The operator must determine the degree of restraint necessary to make sure that the patient is safe.
- Do not sit on the siderails.

CAUTION

Do not use the siderails to move the product. Always move the product using the integrated handles in the headboard and footboard.

You must raise and lower the siderails with both hands. Siderails only lock in the full up position.

When raising the siderails, listen for the “click” that indicates that the siderail is locked in the raised position. Pull on the siderail to make sure that it is locked into position.

To raise the siderails, grasp and lift the siderail.

To lower the siderails, lift the yellow release latch (A) ([Figure 1-18 on page 1-26](#)) and guide the siderail down.



Figure 1-18: Raising or lowering the siderails

Nurse control panel (outside siderail) (option)

WARNING

Always lock bed motion controls when the patient is unattended.

| | Name | Function |
|---|-----------------------|--|
| A | Emergency CPR | Overrides the control panel lockout to achieve the flat position at low height. Also available if the control panels are turned off. |
| B | Litter lockout LED | Illuminates when you lock the litter section |
| C | Litter up | Raises the litter |
| D | Litter down | Lowers the litter |
| E | Backrest lockout LED | Illuminates when you lock the backrest section |
| F | Backrest up | Raises the backrest |
| G | Backrest down | Lowers the backrest |
| H | Upper leg lockout LED | Illuminates when you lock the upper leg section |
| I | Upper leg up | Raises the upper leg section |
| J | Upper leg down | Lowers the upper leg section |
| K | Trendelenburg | Places the product into the Trendelenburg position (head down with foot up) |

Operation

English
EN

Nurse control panel (outside siderail) (option) (Continued)

| | Name | Function |
|----------|------------------------------|--|
| L | Reverse Trendelenburg | Places the product into the Reverse Trendelenburg position (head up with foot down) |
| M | Plug indicator | Illuminates when the product is plugged in |
| N | Battery charge indicator | Illuminates amber when you connect the product to a wall outlet and the batteries are recharging. The battery has a full charge within 10 to 12 hours. When the battery is charged, the LED no longer illuminates. |
| O | Low height indicator | Illuminates green when the product is within 2 cm from its lowest height position |
| P | Examination position | Flattens the litter and raises the litter to the highest height |
| Q | One-button Vascular position | Overrides the control panel lockout to achieve 12° Trendelenburg |


Patient control panel (inside siderail) (option)

WARNING

Always lock bed motion controls when the patient is unattended.

Patient control panel (inside siderail) (option) (Continued)

| | Name | Function |
|---|-----------------------|--|
| A | Autocontour up | Raises the backrest and the upper leg section at the same time |
| B | Autocontour down | Lowers the backrest and the upper leg section at the same time |
| C | Litter lockout LED | Illuminates when you lock the litter section |
| D | Litter up | Raises the litter |
| E | Litter down | Lowers the litter |
| F | Backrest lockout LED | Illuminates when you lock the backrest section |
| G | Backrest up | Raises the backrest |
| H | Backrest down | Lowers the backrest |
| I | Upper leg lockout LED | Illuminates when you lock the upper leg section |
| J | Upper leg up | Raises the upper leg section |
| K | Upper leg down | Lowers the upper leg section |
| L | Egress | Lowers the litter, lowers the upper leg section, and raises the backrest so the patient can enter and exit the product |
| M | Chair position | Place the product into the chair position |

Operation

English
EN


Patient control pendant (option)

WARNING

Always lock bed motion controls when the patient is unattended.

CAUTION

- Always place the patient control pendant safely on the mattress while the pendant is in use.
- Always hang the patient control pendant onto the siderail when the pendant is not in use.
- Do not squeeze or pinch the pendant cord in the bed frame.


| | Name | Function |
|---|------------------|--|
| A | Backrest up | Raises the backrest |
| B | Backrest down | Lowers the backrest |
| C | Upper leg up | Raises the upper leg section |
| D | Upper leg down | Lowers the upper leg section |
| E | Autocontour up | Raises the backrest and the upper leg section at the same time |
| F | Autocontour down | Lowers the backrest and the upper leg section at the same time |
| G | Litter up | Raises the litter |
| H | Litter down | Lowers the litter |

WARNING

- Always lock bed motion controls when the patient is unattended.
- Never store the nurse control panel within the patient's reach.

CAUTION

- Always place the nurse control pendant onto the footboard.
- Always hang the nurse control pendant onto a foot end siderail or store in the linen tray (option) before removing the footboard.
- Do not squeeze or pinch the pendant cable in the bed frame.

| | Name | Function |
|----------|-----------------------------------|--|
| A | Emergency CPR | Overrides the control panel lockout to achieve the flat position at low height. Also available if the control panels are turned off. |
| B | Litter lockout/Litter lockout LED | Enables or disables locks for litter motion. Illuminates when you lock the litter section. |
| C | Litter up | Raises the litter |

Nurse control pendant (Continued)

| | Name | Function |
|---|--|--|
| D | Litter down | Lowers the litter |
| E | Backrest up lockout/Backrest lockout LED | Enables or disables locks for the backrest. Illuminates when you lock the backrest. |
| F | Backrest up | Raises the backrest |
| G | Backrest down | Lowers the backrest |
| H | Upper leg lockout/Upper leg lockout LED | Enables or disables locks for the upper leg section. Illuminates when you lock the upper leg section. |
| I | Upper leg up | Raises the upper leg section |
| J | Upper leg down | Lowers the upper leg section |
| K | Battery charge indicator | Illuminates amber when you connect the product to a wall outlet and the batteries are recharging. The battery has a full charge within 10 to 12 hours. When the battery is charged, the LED no longer illuminates. |
| L | Examination position | Flattens the litter and raises the litter to the highest height |
| M | Plug indicator | Illuminates when the product is plugged in |
| N | Chair position | Place the product into the chair position |
| O | One-button Vascular position | Overrides the control panel lockout to achieve 12° Trendelenburg |
| P | Autocontour up | Raises the backrest and the upper leg section at the same time |
| Q | Autocontour down | Lowers the backrest and the upper leg section at the same time |
| R | Trendelenburg | Places the product into the Trendelenburg position (head down with foot up) |
| S | Reverse Trendelenburg | Places the product into the Reverse Trendelenburg position (head up with foot down) |
| T | Egress | Lowers the litter, lowers the upper leg section, and raises the backrest so the patient can enter and exit the product |

Extending the bed extender (option)

WARNING

- Do not sit on the bed extender. This may cause the product to tip.
- Always lock the bed extender before you place weight on the bed extender.

Extending the bed extender (option) (Continued)

CAUTION

- Do not remove the footboard after you extend the bed extender.
- Do not raise the lower leg section while the bed extender is in use. This is to avoid the situation where the product does not support the lower legs of a taller patient.

The bed extender allows you to extend the length of the product by 31 cm.

To extend the bed extender:

1. Pull and turn each knob 90° to unlock the bed extender ([Figure 1-19 on page 1-33](#)).
2. Grasp the footboard handles.
3. Pull the footboard to extend the bed extender ([Figure 1-20 on page 1-33](#)).
4. Pull and turn each knob 90° to lock the bed extender into place.

Note: When locking the bed extender, listen for the “click” that indicates that the bed extender is locked. Push and pull on the footboard to make sure that the bed extender is locked.


Figure 1-19: Unlocking the bed extender

Figure 1-20: Extending the bed extender

Installing the bed extender bolster mattress

For mattress specifications, see the MA series mattress manual.

The recommended bed extender bolster mattresses are:

| Compatible bolster mattresses | Dimension |
|-------------------------------|--------------------------|
| 7002-4-018 | 330 mm x 710 mm x 180 mm |
| 7002-4-020 | 330 mm x 710 mm x 200 mm |
| 7002-4-518 | 330 mm x 710 mm x 180 mm |
| 7002-4-520 | 330 mm x 710 mm x 200 mm |

To install the bed extender bolster mattress:

1. See [Extending the bed extender \(option\) on page 1-32](#).
2. Place the bed extender bolster mattress between the mattress and the footboard.
3. Press down on the bed extender bolster mattress to secure it into place.

Extending or retracting the linen tray (option)

The linen tray is an optional built-in storage unit that can store a patient's clothes, laundry, or the nurse control pendant. You can find the linen tray at the foot end of the product.

WARNING

- Always retract the linen tray (option) before you put the product in motion.
- Always retract the linen tray (option) when not in use.

CAUTION

The safe working load of the linen tray is 15 kg.

To extend the linen tray, grasp the plastic linen tray and pull the linen tray out toward you.

To retract the linen tray, grasp the plastic linen tray and push the linen tray into the frame.

Figure 1-21: Storing the nurse control pendant

Inserting or removing a cassette from the X-ray cassette holder (option)

WARNING

Do not use the product for X-ray procedures if it does not have the radiolucent backrest (option).

Inserting or removing a cassette from the X-ray cassette holder (option) (Continued)

The SV2 may include an optional radiolucent backrest to allow X-ray images to be taken while the patient is on the bed.

You can take X-ray images by inserting an X-ray cassette into the housing that is located behind the backrest. You do not need to move the patient to insert an X-ray cassette or take an X-ray.

X-ray guide dimensions: 390 mm x 660 mm x 16 mm

The recommended X-ray cassette dimensions are:

- 385 mm x 385 mm x 15 mm
- 460 mm x 383 mm x 15 mm

To insert an X-ray cassette:

1. See [Removing or replacing the headboard on page 1-23](#).
2. Slide the X-ray cassette into the X-ray cassette holder.
3. Adjust the patient to the desired position.

To remove an X-ray cassette:

1. Slide the X-ray cassette out of the X-ray cassette holder.
2. See [Removing or replacing the headboard on page 1-23](#).

Accessories

English
EN

These accessories may be available for use with your product. Confirm availability for your configuration or region.

CAUTION

Only use authorized accessories for this product. The use of unauthorized accessories may result in product damage or injury to the operator or patient. Stryker is not responsible for any damage or injury that may result from the misuse of the product or the use of unauthorized accessories.

| Name | Number | Safe working load |
|--|--------|--------------------|
| IV pole | MM017 | Each IV hook: 2 kg |
| Patient lift pole | MM003 | 75 kg |
| Upright oxygen bottle holder (120 mm diameter) | MM006 | 7.5 kg |
| Foley bag basket | MM029 | 4 kg |

Installing the IV pole

WARNING

Do not use accessories to support patient limbs or other body parts.

CAUTION

- Always make sure that accessories are locked into position.
- Always make sure that the IV pole is at a low height during transport.
- Do not use the IV pole as a push or pull device.
- Do not allow accessories to interfere with mechanical or electric mechanisms of the product.

You can install the IV pole into any of the four accessory sleeves on the corners of the bed. The IV pole has a telescopic pole that extends to provide a second height position.

To install the IV pole:

1. Insert the IV pole into one of the four accessory sleeves ([Figure 1-22 on page 1-36](#)).
2. Rotate and lock the IV pole into the accessory sleeve ([Figure 1-23 on page 1-36](#)).


Figure 1-22: Installing the IV pole

Figure 1-23: Locking the IV pole

Adjusting the IV pole

To adjust the IV pole:

1. Turn the telescoping knob counter clockwise to unlock the IV pole ([Figure 1-24 on page 1-37](#)).
2. Grasp the IV pole.
3. Raise the IV pole to the desired height.
4. Turn to the telescoping knob clockwise to lock the IV pole ([Figure 1-24 on page 1-37](#)).


Figure 1-24: Adjusting the IV pole

Installing the lifting pole

The lifting pole assists the patient with changing position in bed.

WARNING

Do not use accessories to support patient limbs or other body parts.

CAUTION

- Always make sure that accessories are locked into position.
- Always remove the lifting pole before transporting the product.
- Do not use the lifting pole as a push or pull device.
- Do not allow accessories to interfere with mechanical or electric mechanisms of the product.

You can install the lifting pole into either of the two accessory sleeves on the head end of the bed.

To install the lifting pole:

Installing the lifting pole (Continued)

1. Insert the lifting pole into one of the two accessory sleeves ([Figure 1-25 on page 1-38](#)).


Figure 1-25: Installing the lifting pole

2. Rotate and lock the lifting pole into the accessory sleeve.

Installing the lifting pole handle

To install the lifting pole handle, place the black grip of the lifting pole between the two stoppers on the lifting pole ([Figure 1-26 on page 1-38](#)).


Figure 1-26: Installing the lifting pole handle

Installing the oxygen bottle holder

WARNING

Do not use accessories to support patient limbs or other body parts.

CAUTION

- Always make sure that accessories are locked into position.
- Do not use the oxygen bottle holder as a push or pull device.
- Always turn the oxygen bottle holder in toward the bed before transporting a patient.
- Do not strike the oxygen bottle holder while transporting a patient.
- Do not allow accessories to interfere with mechanical or electric mechanisms of the product.

To install the oxygen bottle holder:

1. Insert the oxygen bottle holder into one of the two accessory sleeves near the head end ([Figure 1-27 on page 1-39](#)).


Figure 1-27: Installing the oxygen bottle holder

2. Rotate and lock the oxygen bottle holder into the accessory sleeve ([Figure 1-28 on page 1-39](#)).

Figure 1-28: Locking the oxygen bottle holder

Installing the Foley bag basket

WARNING

Do not use accessories to support patient limbs or other body parts.

Installing the Foley bag basket (Continued)

CAUTION

- The safe working load of each Foley hook is 2 kg.
- Do not allow accessories to interfere with mechanical or electric mechanisms of the product.

To install the Foley bag basket, hook the basket onto the Foley hooks ([Figure 1-29 on page 1-40](#)).


Figure 1-29: Installing the Foley bag basket

Preparing the product for cleaning

Cleaning and disinfecting are two separate processes. Clean before disinfecting to make sure that the cleaning agent is effective.

To prepare the product for cleaning:

1. Raise the litter to its highest height.
2. Lock the siderail control panel and patient control pendant functions (see [Nurse control pendant on page 1-31](#)).
3. Unplug the power cord from the wall outlet.
4. See [Storing the power cord on page 1-18](#).
5. See [Applying or releasing the brakes on page 1-19](#).
6. Remove the mattress.

Cleaning

WARNING

- Do not clean, service, or perform maintenance while the product is in use.
- Always power off and unplug the power cord before cleaning, servicing, or performing maintenance.
- Always power off the product and unplug the power cord from the wall outlet when large spills occur near the circuit boards, cables, and motors. Remove the patient from the product, clean up the fluid, and have service personnel inspect the product. Fluids can cause unpredictable operation and decreased functionality of any electrical product. Do not return the product to service until it is completely dry and has been thoroughly tested for safe operation.
- Do not spray cleaners directly onto the battery, control boxes, actuators, cables, or other electric equipment.
- Do not use abrasive powder, steel wool, or similar materials that may damage the product surface.
- Do not use **Virex**® TB for product disinfecting.
- Do not use acid-based chemicals or flammable chemicals, such as gasoline, diesel, or acetone for cleaning purposes.
- Do not directly spray or saturate the siderail control panel, patient control pendant, or nurse control pendant with cleaners.
- The cleaners and disinfectants must not be highly alkaline or acidic (pH value 6-8).

CAUTION

- Do not steam clean, pressure wash, ultrasonically clean, or immerse any part of the product in water. Exposure to water may damage the internal electric parts. These methods of cleaning are not recommended and may void this product's warranty.
- Always make sure that you wipe each product with clean water and thoroughly dry each product after cleaning. Some cleaning products are corrosive in nature and may cause damage to the product if you use them improperly. If you do not properly rinse and dry the product, you may leave a corrosive residue on the surface of the product that could cause premature corrosion of critical components. Failure to follow these cleaning instructions may void your warranty.

To clean product surfaces:

1. Using a clean, soft, damp cloth, wipe product surfaces with a mild soap and water solution to remove foreign material.
2. Wipe product surfaces with a clean, dry cloth to remove any excess liquid or cleaning agent.
3. Dry thoroughly.

Cleaning the siderails

WARNING

- Do not directly spray or saturate the siderail control panel, patient control pendant, or nurse control pendant with cleaners.
 - Do not use sharp objects to clean the siderail control panel.
 - Do not use abrasive powder, steel wool, or similar materials that may damage the product surface.
 - Do not use **Virex® TB** for product cleaning.
 - Do not use acid-based chemicals or flammable chemicals, such as gasoline, diesel, or acetone for cleaning purposes.
-

To clean the siderails:

1. Raise the siderail.
2. Latch the siderail.
3. Use a clean, soft, damp cloth to wipe down the siderail and the siderail control panel.
4. Allow the siderail control panel to dry thoroughly.

Suggested disinfectants:

- Quaternary cleaners without glycol ethers (active ingredient - ammonium chloride)
- Chlorinated bleach solution (5.25% - less than 1 part bleach to 100 parts water)
- 70% Isopropyl alcohol

Always follow the disinfectant's instructions for appropriate contact time and rinsing requirements.

Avoid oversaturation and make sure that the product does not stay wet longer than the chemical manufacturer's guidelines for proper disinfecting.

To disinfect the product:

1. Thoroughly clean and dry the product before you apply disinfectants.
2. Apply recommended disinfectant solution by spray or pre-soaked wipes
Note: Make sure that you follow the disinfectant's instructions for appropriate contact time and rinsing requirements.
3. To disinfect mechanisms, lift the backrest and leg rest up to the highest height.
4. Wipe product surfaces and mechanisms with a clean, dry cloth to remove any excess liquid or cleaning agent.
5. Allow the product to dry completely before returning to service.

Preventive maintenance

English
EN

At a minimum, check all items listed during annual preventive maintenance for all Stryker Medical products. You may need to perform preventive maintenance checks more frequently based on your level of product usage.

Remove product from service before performing preventive maintenance. Preventive maintenance should only be performed by trained or certified personnel. See the maintenance manual for servicing and maintenance instructions.

Inspect the following items:

- _____ All welds and all fasteners are secure
- _____ Tubing or sheet metal for bends or breaks
- _____ Casters are free of debris
- _____ Casters are secure and swivel
- _____ Casters lock securely by depressing the brake pedal
- _____ Locking steer caster applies and releases
- _____ Steer pedal latches
- _____ Backrest operates
- _____ Litter up and down operates
- _____ Trendelenburg and Reverse Trendelenburg operates
- _____ IV pole is intact and operating (optional)
- _____ Accessory sleeves are not damaged or cracked
- _____ Bed extender extends and locks (option)
- _____ Headboard, footboard, and siderail panels for cracks or splits
- _____ All covers are not damaged and do not have sharp edges
- _____ Radiolucent backrest is clean and not cracked (option)
- _____ Cassette holder is clean and not cracked (option)
- _____ Underbed light operates
- _____ CPR release operates
- _____ Siderails move, latch, and stow
- _____ All functionality on all control panels
- _____ Batteries for replacement
- _____ Batteries for corrosion at the terminals, cracking, expanded or bulging at the sides, or can no longer maintain a full charge
- _____ Lower leg section moves, latches, and stows
- _____ Pendants for any physical damage
- _____ Power cord not worn or frayed
- _____ Cables not worn or pinched
- _____ All electrical connections tight
- _____ All grounds secure to the frame
- _____ Ground Impedance Check ($\leq 0.2 \text{ Ohm}$)
- _____ Leakage current: Normal Polarity, No Ground, L2 Active ($\leq 300 \text{ }\mu\text{A}$)
- _____ Leakage current: Normal Polarity, No Ground, No L2 ($\leq 600 \text{ }\mu\text{A}$)
- _____ Leakage current: Reverse Polarity, No Ground, L2 Active ($\leq 300 \text{ }\mu\text{A}$)
- _____ Leakage current: Reverse Polarity, no Ground, No L2 ($\leq 600 \text{ }\mu\text{A}$)
- _____ Enclosure is free from wear, tear, stresses and mechanical damage
- _____ High potential test 1500 VAC (trip current not more than 10 mA)
- _____ No rust or corrosion of parts
- _____ Control boxes are not damaged or cracked
- _____ Actuator functionality
- _____ Labels for legibility, proper adherence, and integrity

Preventive maintenance

| |
|------------------------|
| Product serial number: |
| Completed by: |
| Date: |

English
EN

EMC information

English
EN

Guidance and Manufacturer's Declaration - Electromagnetic Emissions

SV2 is intended for use in an electromagnetic environment specified below. The customer or the user of **SV2** should assure that it is used in such an environment.

| Emissions Test | Compliance | Electromagnetic Environment |
|--|------------|--|
| RF Emissions CISPR 11 | Group 1 | SV2 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF Emissions CISPR 11 | Class A | SV2 is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic Emissions IEC 61000-3-2 | Class A | |
| Voltage Fluctuations Flicker Emissions IEC 61000-3-3 | Complies | |

Note: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

SV2 is suitable for use in the electromagnetic environment specified below. The customer or the user of **SV2** should assure that it is used in such an environment.

| Immunity Test | IEC 60601 Test Level | Compliance Level | Electromagnetic Environment Guidance |
|---|--|--|--|
| Electrostatic Discharge (ESD) IEC 61000-4-2 | ± 8 kV contact ± 15 kV air | ± 8 kV contact ± 15 kV air | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrostatic fast Transient/burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input/output lines | ± 2 kV for power supply lines ± 1 kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | Mains power quality is that of a typical commercial and/or hospital environment. |

(Continued)

| Guidance and Manufacturer's Declaration - Electromagnetic Immunity | | | |
|--|--|--|---|
| Voltage dips, voltage variations and short interruptions on power supply input lines IEC 61000-4-11 | 0%U _T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0%U _T for 1 cycle 70%U _T (30% dip in U _T) for 25 cycles 0% U _T for 250 cycles | 0%U _T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0%U _T for 1 cycle 70%U _T (30% dip in U _T) for 25 cycles 0% U _T for 250 cycles | Mains power quality should be that of a typical commercial and/or hospital environment. If the user of SV2 requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery. |
| Power frequency magnetic fields IEC 61000-4-8 | 30 A/m | 30 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment. |
| Note: U _T is the a.c. mains voltage prior to applications of the test level. | | | |

| Guidance and Manufacturer's Declaration - Electromagnetic Immunity | | | |
|---|--|--|---|
| SV2 is suited for use in the electromagnetic environment specified below. The customer or the user of SV2 should assure that it is used in such an environment. | | | |
| Immunity Test | IEC 60601 Test Level | Compliance Level | Electromagnetic Environment - Guidance |
| Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3 | 3 Vrms 6 Vrms in ISM bands 150 kHz to 80 MHz 3 V/m 80 MHz to 2.7 GHz | 3 Vrms 6 Vrms in ISM bands 3 V/m | Portable and mobile RF communications equipment should be used no closer to any part of SV2 , including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distance $D=(1.2) (\sqrt{P})$ $D=(1.2) (\sqrt{P})$ 80 MHz to 800 MHz $D=(2.3) (\sqrt{P})$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the |

EMC information

English
EN

(Continued)

| | | | |
|--|--|--|---|
| | | | <p>transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
|--|--|--|---|

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Note 3: The ISM (industrial, scientific, and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

Note 4: Evaluated for immunity to proximity fields from RF wireless communication equipment per IEC 60601-1-2:2014 Table 9.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which **SV2** is used exceeds the applicable RF compliance level above, **SV2** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating **SV2**.

^b Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and SV2

SV2 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of **SV2** can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and **SV2** as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter W | Separation distance according to frequency of transmitter m | | |
|--|--|---|--|
| | 150 kHz to 80 MHz $D=(1.2) (\sqrt{P})$ | 80 MHz to 800 MHz $D=(1.2) (\sqrt{P})$ | 800 MHz to 2.7 GHz $D=(2.3) (\sqrt{P})$ |
| 0.01 | 1.12 | 0.12 | 0.23 |

EMC information

English
EN

(Continued)

| Recommended separation distances between portable and mobile RF communications equipment and SV2 | | | |
|--|------|------|------|
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |
| <p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> | | | |

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Stryker European Operations B.V.