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- E. J. Kujper, B. Cegard and P. Tull: Emergence of *Clostridium difficile* associated disease in North America and Europe, Review Clinical Microbiology and Infections, 12 (suppl), p. 2-18, Oct. 2006.
- Leverly D.M., H.C. Kivan and D.J. Wilkins: *Clostridium difficile*: its disease and toxin. Clinical Microbiology Reviews, p. 2-18, Jan. 1988.
- Ramey L. et al: Fulminant *Clostridium difficile* an underappreciated and increasing cause of death and complications. Annals of Surgery 235 (3) p. 363-372, Mar. 2002.
- Wen MW, Kinross R, Sivapalan M, Shenkin M, Shetty NR: Detection of *Clostridium difficile* infection: a suggested laboratory diagnostic algorithm. British Journal of Biomedical Science, 66(4) p. 175-179, 2009.
- Wills DH, And JA Kraft: Confirmation that the latex-reactive protein of *Clostridium difficile* is a Glutamate Dehydrogenase. Journal of clinical microbiology, 30, p. 1263-1264, May 1992.

Rev. 0, 2019-01-05 OM/SDe

- Intended Use**
The NADAL® C. difficile Toxin A/B+GDH Test is a lateral flow, immunochromatographic rapid test for the qualitative detection of *Clostridium difficile* Toxin A and Toxin B (C. difficile Toxins A/B), as well as *Clostridium difficile* glutamate dehydrogenase (C. difficile GDH) in human faecal specimens. The test is intended for use as an aid in the diagnosis of *Clostridium difficile* infection and is designed for professional use only.

- Introduction and Clinical Significance**
Clostridium difficile (C. difficile) is an anaerobic gram-positive, spore-forming bacterium. The ability of C. difficile to form spores is the key feature which enables this bacterium to survive in patients and the physical environment for long periods, thereby facilitating its transmission. C. difficile is transmitted by the faecal-oral route.

C. difficile is the principal pathogen related to antibiotic-related diarrhoea and/or pseudomembranous colitis in hospitalised patients. Mature colonic bacterial flora in a healthy adult is generally resistant to C. difficile colonisation. However, if the normal enteric flora is altered, resistance to colonisation is lost. Any factor associated with the alteration of the normal enteric flora increases the risk of C. difficile colonisation after exposure to antibiotics, especially those with broad-spectrum activity such as penicillins, cephalosporins and clindamycin.

C. difficile releases two high-molecular-weight toxins: Toxin A and B, which are responsible for the clinical manifestations ranging from mild, self-limiting, watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon and death. *Clostridium difficile* glutamate dehydrogenase (C. difficile GDH) is an enzyme produced in large quantities by all toxigenic and non-toxic strains, making it an excellent biomarker.

- Test Principle**
The NADAL® C. difficile Toxin A/B+GDH Test is a lateral flow, immunochromatographic rapid test for the qualitative detection of C. difficile Toxin A and Toxin B, as well as C. difficile GDH in human faecal samples.

The test cassette contains three different test strips for the qualitative detection of C. difficile GDH, C. difficile Toxin A and C. difficile Toxin B antigens.

Anti-C. difficile GDH, anti-C. difficile Toxin A and anti-C. difficile Toxin B antibodies are pre-coated in the corresponding test line region (T) of the membrane. During testing, the specimen reacts with other anti-C. difficile GDH, anti-C. difficile Toxin A and anti-C. difficile Toxin B antibodies conjugated to coloured particles. The mixture migrates along the membrane by capillary action and, in the case of a positive sample, reacts with the anti-C. difficile GDH, anti-C. difficile Toxin A and/or anti-C. difficile Toxin B antibodies immobilised in the corresponding test line region (T), generating a coloured line. The presence of the coloured lines indicates a positive result, while its/their absence indicates a negative result.

A coloured line should always appear in each control line region (C) and serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

- Reagents and Materials Supplied**

- 10 NADAL® C. difficile Toxin A/B+GDH test cassettes (incl. disposable pipettes for liquid faecal samples)
- 30 specimen collection tubes with buffer
- 1 package insert

- Additional Materials Required**

- Specimen collection container
- Disposable gloves
- Timer

- Storage & Stability**

Test kits should be stored at 2-30°C and used before the expiry date printed on the packaging. Test cassettes are stable until the expiry date printed on foil pouches. Test cassettes should remain in sealed pouches until use. Do not freeze tests. Do not use tests beyond the expiry date.

- Warnings and Precautions**

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiry date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures, as well as standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

- Specimen Collection and Preparation**

Specimen Collection and Storage:

A sufficient quantity of faeces (1-2 g or 1-2 mL for liquid specimens) should be collected. Faecal specimens should be collected in clean and dry containers containing no preservatives or transport media.

Specimens must be tested as soon as possible after collection. Best results will be obtained if the assay is performed within 6 hours of specimen collection. If specimens are not to be tested within 6 hours of collection, they can be stored

refrigerated at 2-8°C for up to 3 days. For longer storage, specimens must be kept frozen below -20°C.

Extracted specimens in buffer can be stored at 2-8°C for up to 1 week or at -20°C for longer periods of time. Ensure that specimens are not treated with solutions containing formaldehyde or its derivatives.

Frozen samples should be completely thawed and brought to room temperature before testing. Avoid repeated freeze-thaw cycles.

Specimen preparation: Use a separate specimen collection tube with buffer for each specimen.

Hold the specimen collection tube upright and remove the applicator stick by unscrewing the green cap. Be careful not to spill or splash buffer solution from the collection tube.

For solid specimens: Collect approx. 50 mg of faeces (equivalent to the size of ¼ of a pea) by inserting the applicator stick in at least 3 different sites of the faecal specimen. Avoid scooping the faecal specimen.

For liquid specimens: Holding a pipette vertically, collect a faecal specimen and add 2 drops (approx. 80 µL) to the specimen collection tube containing buffer.

Place the applicator stick, along with the faecal specimen (or without the faecal specimen if the specimen is liquid), back into the specimen collection tube and tightly screw the green cap closed. Be careful not to break the tip of the tube's cap. Shake the specimen collection tube vigorously to thoroughly mix the specimen and buffer. Let the solution stand for 2 minutes.

9. Test Procedure

Bring the tests and samples with buffer to room temperature (15-30°C) prior to testing. Do not open the foil pouch until you are ready to perform the assay.

- Remove the NADAL® C. difficile Toxin A/B+GDH test cassette from the sealed foil pouch and use it as soon as possible. Place the test cassette on a clean and flat surface. Be sure to label the test cassette with the specimen ID number.

- Hold the specimen collection tube upright and unscrew the transparent protective cap. Dispense 3 full drops (approx. 120 µL) of the specimen solution into each specimen well (S) of the test cassette. Avoid trapping air bubbles in the specimen wells (S).

- Start the timer.

- Wait for the coloured line(s) to appear. Read the test result 10 minutes after dispensing the specimen into the specimen wells. Do not interpret the result after more than 20 minutes.

Note: If the extracted specimen does not migrate due to the presence of particles, it should be centrifuged. Then collect 120 µL of the supernatant, add it to the specimen wells (S) of a new test cassette and start again, following the instructions mentioned above.

10. Result Interpretation

Positive for GDH: One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'T' for GDH.

Positive for Toxin A: One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'T' for Toxin A (Tox A).

Positive for Toxin B: One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'T' for Toxin B (Tox B).

Note: A combination of positive results for various parameters is possible.

Note: The colour intensity in the test line region (T) may vary depending on the concentration of C. difficile antigens present in the specimen. Therefore, any shade of colour in the test line region (T) should be considered positive.

Negative: One coloured line appears in each control line region 'C'. No coloured lines appear in the test line regions 'T' for GDH, Toxin A (Tox A) and Toxin B (Tox B).

Invalid: A coloured control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

- Quality Control**
An internal procedural control is included in the test cassette: A coloured line developing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® C. difficile Toxin A/B+GDH Test is for professional *in-vitro* diagnostic use only.
- The NADAL® C. difficile Toxin A/B+GDH Test should only be used for the qualitative detection of C. difficile antigens. The quantitative value of C. difficile antigens cannot be determined using this qualitative test.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- A positive result does not exclude the presence of other pathogens.

13. Performance Characteristics

Detection limit

The lower detection limit of the NADAL® C. difficile Toxin A/B+GDH Test is 1 ng/mL for GDH, 2 ng/mL for Toxin A and 1 ng/mL for Toxin B.

Sensitivity and specificity

The performance of the NADAL® C. difficile Toxin A/B+GDH Test has been evaluated in comparison with another leading, commercially available rapid test.

NADAL® C. difficile Toxin A/B+GDH Test (c/dH)	Other rapid test			Total
	Positive	Negative	Total	
Positive	78	2	80	
Negative	1	119	120	
Total	79	121	200	

Relative sensitivity: 98.7% (95% CI: *93.1% - 100%*)

Relative specificity: 98.3% (95% CI: *94.2% - 99.8%*)

Overall agreement: 98.5% (95% CI: *95.7% - 99.7%*)

*95% confidence interval

NADAL® C. difficile Toxin A/B+GDH Test (Toxin A/B)	Other rapid test			Total
	Positive	Negative	Total	
Positive	56	2	58	
Negative	1	141	142	
Total	57	143	200	

Relative sensitivity: 98.2% (95% CI: *90.6% - 99.9%*)

Relative specificity: 98.6% (95% CI: *95.0% - 99.8%*)

Overall agreement: 98.5% (95% CI: *95.7% - 99.7%*)

*95% confidence interval

Repeatability and reproducibility

To check intra-lot accuracy (repeatability), positive and negative specimens were tested 15 times using test kits from

the same production lot in the same experimental conditions. All obtained results were confirmed as expected.

To check inter-lot accuracy (reproducibility), positive and negative specimens were tested using test kits from 3 different production lots. All obtained results were confirmed as expected.

Cross-reactivity

An evaluation was performed to determine the cross-reactivity of the NADAL® C. difficile Toxin A/B+GDH Test. No cross-reactivity with common gastrointestinal pathogens, which are occasionally present in faeces, was observed:

<i>Campylobacter jejuni</i>	<i>Salmonella Enteritidis</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	<i>Salmonella Paratyphi</i>	<i>Shigella flexneri</i>
<i>Escherichia coli O157:H7</i>	<i>Salmonella Typhi</i>	<i>Shigella sonnei</i>
<i>Helicobacter pylori</i>	<i>Salmonella Typhimurium</i>	<i>Staphylococcus aureus</i>
<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>	<i>Yersinia enterocolitica</i>

14. References

- Ramstedt Salmann, V. Saly and Saarbrunn, S. Ramonova: Estimation of faecal carriage of *Clostridium difficile* in patients with ulcerative colitis using real time polymerase chain reaction, Indian Journal of Medical Research, p. 472-477, May 2008.
- E. J. Kujper, B. Cegard and P. Tull: Emergence of *Clostridium difficile* associated disease in North America and Europe, Review Clinical Microbiology and Infections, 12 (suppl), p. 2-18, Oct. 2006.
- Leverly D.M., H.C. Kivan and D.J. Wilkins: *Clostridium difficile*: its disease and toxin. Clinical Microbiology Reviews, p. 2-18, Jan. 1988.
- Ramey L. et al: Fulminant *Clostridium difficile*: an underappreciated and increasing cause of death and complications. Annals of Surgery 235 (3) p. 363-372, Mar. 2002.
- Wen MW, Kinross R, Sivapalan M, Shenkin M, Shetty NR: Detection of *Clostridium difficile* infection: a suggested laboratory diagnostic algorithm. British Journal of Biomedical Science, 66(4) p. 175-179, 2009.
- Wills DH, And JA Kraft: Confirmation that the latex-reactive protein of *Clostridium difficile* is a Glutamate Dehydrogenase. Journal of clinical microbiology, 30, p. 1263-1264, May 1992.

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14. Bibliografia

1. Ramakrishna Balarajun, V. Bala; and Balakrishna S. Ramakrishna. Estimation of faecal carriage of *Clostridium difficile* in patients with ulcerative colitis using real time polymerase chain reaction. Indian Journal of Medical Research, p.452-457, May 2008.
2. E. J. Kujper, B. Cogliard and P. Tiel. Emergence of *Clostridium difficile*-associated disease in North America and Europe. Review Clinical Microbiology and Infectious 22 (suppl. 5): i240-2008.
3. Leventy D.M., H.C. Knapp and D.J. Wilcox. *Clostridium difficile*: its disease and treatment. Clinical Microbiology Reviews, 1: 23-38, 1988.
4. Ramsey J, et al. Fulminant *Clostridium difficile* infection: an underappreciated and increasing cause of death and complications. Annals of Surgery 205 (3): 360-372, Mar. 2002.
5. Wren MW, Knorr A, Gopalak M, Sharma M, Shetty NC. Detection of *Clostridium difficile* infection: a suggested laboratory diagnostic algorithm. British Journal of Biomedical Science, 66(4): 175-179, 2009.
6. Wells JG, Ant J, Knorr A. Confirmation that the toxin-negative strain of *Clostridium difficile* is *Clostridium difficile*. Journal of Clinical Microbiology, 30, 3: 388-389, May 1992.

Symbol	Deutsch	English	Français	Español	Italiano	Polak
CE	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
IVD	Gebrauchsanweisung beachten	Consult instructions for use	Consultez la notice d'utilisation	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso	Przeznaczony do diagnostyki in vitro
↓	Temperaturbegrenzung	Temperature limitation	Limites de température	Limites de temperatura	Limiti di temperatura	Temperatura przechowywania
LOT	Chargeberechnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
⊗	Nicht zur Weiterverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
📄	Verwendbar bis	Use by	Utiliser jusqu'à	Fecha de caducidad	Utilizzare entro	Data ważności
REF	Bestimmnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Momento di Catalogo	Numer katalogowy
🏭	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
⚠️	Ausreichend für non-Ausweis	Sufficient for non-test	Suffisant pour non-test	Suficiente para non-utilización	Sufficiente per non-Test	Wystarczająco dla non-Przeznaczony

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© 2018 nal von minden GmbH • Carl-Zeiss-Strasse 12 • 47445 Moers • Germany • info@nal-von-minden.com • www.nal-von-minden.com 2