Product label:
Puristeril 340 (F00006606 and F00006609)

04b

2015-04-28

1 / 5

Rev.no.

Date
YYYY-MM-DD

Page

(not original size)

EAN

bg опасност
cs Nebezpečí
el Κίνδυνος
en Danger
hr Opasnost
hu Veszély
no Fare
pl Niebezpieczeństwo

ro Periculos
sk Nebezpečenstvo
sl Nevarnost
sr Opasno
tr Tehlike

CE 0123

REF
LOT

5°C 25°C

size

**Puristeril
340**

bg За студена химическа дезинфекция на апарати за хемодиализа
Съдържание: пероцетна киселина, водороден пероксид
cs Na studenou dezinfekci hemodialyzačních přístrojů
Obsahuje: kyselina peroctová, peroxid vodíku
el Για Χημική Απολύμανση Μηχανημάτων Αιμοκάθαρσης
Περιεχόμενα: Παρακετικό οξύ, υ εροξείδιο υδρογόνου
en For Cold Disinfection of Haemodialysis Machines
contents: peracetic acid, hydrogen peroxide
hr Za hladnu dezinfekciju aparata za hemodijalizu
sadrži: peracetatnu kiselinu, vodikov peroksid
hu Hemodializáló készülékek hideg sterilizálásához
összetevői: peroxiecetsav és hidrogén-peroxid.
no Til desinfeksjon av hemodialysemaskiner
Innhold: Pereddiksyre, hydrogenperoksyd.
pl Do dezynfekcji na zimno aparatów do hemodializ
Zawiera: kwas nadoctowy, nadtlenek wodoru
ro Pentru dezinfectia la rece a aparatelor de hemodializa
conține: acid peracetic, peroxid de hidrogen.
sk Prostriedok na studenú dezinfekciu hemodialyzačných prístrojov
Obsahuje: kyselinu peroctovú, peroxid vodíka.
sl Za hladno dezinfekcijo hemodializnih aparatov
Sestava: perocetna kislina, vodikov peroksid.
sr Za hladnu dezinfekciju aparata za hemodijalizu
sadrži: persirćetna kiselina, vodonik peroksid.
tr Hemodializ Cihazlarının Soğuk Dezenfeksiyonu içindir
içeriği: Perasetik asit, hidrojen peroksit.

L-F00006606-00000 (05/2015)

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Page

en

For Cold Disinfection of Haemodialysis Machines

Contents: peracetic acid, acetic acid, hydrogen peroxide, water.

Action: Bactericidal, fungicidal, sporicidal, virucidal, decalcifying, cleaning

Instructions for use:

For disinfection, cleaning and decalcification of haemodialysis systems. Aspirate Puristeril 340 through the disinfectant hose. For dosage (concentration, temperature, contact time) follow according to the machine manufacturer's instruction. After disinfection the machine is rinsed with water until the residues are no longer detected with potassiumiodide starch paper or peracetic acid test strips (REF 6299161), i. e. no coloration must be observed on the strips. For disinfection of water treatment systems proceed according to the manufacturer's instruction.

Safety precautions:

Causes severe skin burns and eye damage. May be corrosive to metals. May cause respiratory irritation. Harmful if swallowed or if inhaled. May intensify fire; oxidiser. Harmful to aquatic life with long lasting effects. Keep away from dirt, rust, chemicals, alkalis and acids as well as heavy metal salts and amines. Wear protective gloves/protective clothing/eye protection/face protection.

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

Store in a well-ventilated place. Keep cool.

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