Brilliance™ GBS Agar

REF PO5320A & PB5260E*

* This IFU is intended to be read in conjunction with the IFU for Columbia CNA agar available at www.thermofisher.com

Intended Use

Brilliance GBS Agar is a qualitative, selective, chromogenic screening medium for the presumptive identification of Group B Streptococci (*Streptococcus agalactiae*; GBS). Used for vaginal and rectal samples in a diagnostic workflow for pre-natal screening and to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Group B Streptococci (GBS), are facultatively anaerobic, oxidase-negative, catalase-negative, Gram-positive organisms.^{1,2} GBS typically exhibit 1 mm colonies that are slightly larger than other Streptococci and exhibit β -haemolysis on blood agar, which can be used as an early step in identifying clinical isolates although a small proportion of strains are non-haemolytic.¹

GBS normally colonises the gastrointestinal tract and is part of the vaginal microbiome of some women. Up to 28% of women in the UK carry GBS in the vagina or rectum without any associated symptoms.^{1,2,4}

However, GBS can be responsible for invasive disease in pregnant women (urinary tract infection, amnionitis, endometritis) and in new-borns potentially causing neonatal disease.^{1,2} Approximately 50% of women who are colonised with GBS will transmit the bacteria to their new-borns. In the absence of intrapartum antibiotic prophylaxis, Preterm neonates with GBS early-onset disease are more likely to experience apnoea, require blood pressure support, and need neonatal intensive care.²

Early-onset disease presents within 7 days after birth and occurs due to infection ascending from the maternal genital tract or, very rarely, via the placenta. It is characterised primarily by sepsis, pneumonia, or less frequently meningitis, and is most likely to manifest within the first 12-48 hours after birth.^{1,2}

Late-onset disease presents between 7 days after birth and 2-3 months of age and is associated with acquisition of the organism perinatally or through nosocomial transmission or from the external environment (e.g. hospital). It is characterised by bacteraemia or meningitis.^{1,2,4}

GBS is the leading cause of new-born infection² and as the most frequent cause of severe early-onset (less than 7 days of age) infection in new-born infants.³ Over 80% of early-onset disease cases occur in the first six days after birth with most of these occurring within 12-24 hours of birth.⁴

Principle of Method

Brilliance GBS Agar is a clear, selective chromogenic screening medium incorporating two chromogens. One chromogen yields a pink colour in the presence of GBS due to its phosphatase activity. To allow the medium to differentiate GBS accurately, it contains a second chromogen. Non-GBS colonies are either inhibited or grow as blue or purple colonies. Broad-spectrum antimicrobial components in the medium suppress the growth of group A

and C Streptococci. Also included are compounds to inhibit the growth of *Enterobacteriaceae*, Staphylococci, Enterococci, and group D Streptococci. Antifungals eliminate yeasts. *Brilliance* GBS Agar also incorporates Inhibigen™ technology for targeted inhibition of Enterococci and group D Streptococci ensuring a high level of sensitivity and specificity and a significant reduction in growth of background flora, which could otherwise mask the growth of GBS.

Typical Formula

	grams per litre
Chromogenic mix	25.5
Antibiotic cocktail	0.09
Inhibigen™ mix	4.9
Agar	15.0

Physical Appearance

Colour	Oyster White
Clarity	Transparent
Fill weight	17g \pm 1
pH	6.8 \pm 0.2

Materials Provided

- The pack contains 10 x 90mm agar plates, film wrapped.
- Each plate should only be used once.
- Each pack contains enough plates for 10 individual tests.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Material Safety Data Sheet (MSDS) for safe handling and disposal of the product available at www.thermofisher.com

Specimen Collection, Handling and Storage

Requirements for specimen collection, handling and storage of samples are described in local procedures and guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 58 (Public Health England, 2018).

Procedure

- Allow product to equilibrate to room temperature.
- Inoculate and streak the specimen or a liquid suspension (e.g. Todd Hewitt broth with colistin and nalidixic acid (LIM broth)) onto the medium, using a standard loop.
- Incubate plates aerobically for 18–24 hours at 36±1 °C. There is no need to re-incubate.
- Visually inspect plates to assess colony growth and colour under good lighting.

Interpretation

The presence of pink colonies indicates the sample is presumptive positive for GBS. Non-target organisms are either inhibited or grow as blue or purple colonies.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 - 24 h @ 36° ± 1°C aerobic

Positive Controls	
Colony count is ≥ 50% of the control medium count	
<i>Streptococcus agalactiae</i> (GBS) ATCC® 13813	0.5-1 mm, Pink Colonies
<i>Streptococcus agalactiae</i> NCTC 9993	0.5-1 mm, Pink Colonies
Negative Controls	
Inoculum 10 ⁴ – 10 ⁵ cfu, qualitative, control medium COL+SB and SAB for <i>C. albicans</i> Incubation conditions: 18 – 24 h @ 36 ± 1°C, aerobic.	
<i>Candida albicans</i> ATCC® 10231	No Growth
<i>Staphylococcus saprophyticus</i> ATCC® 15305	No Growth
<i>Streptococcus dysgalactiae</i> ATCC® 12394	No Growth
<i>Escherichia coli</i> ATCC® 25922	No Growth

Inoculum ≥ 10 ⁴ cfu, quantitative control medium TSA Incubation conditions: 18 – 24 h @ 36 ± 1°C, aerobic	
<i>Enterococcus faecium</i> ATCC® 19434	Complete Inhibition (<10 CFU)

Analytical Performance

A study of 100 GBS and 150 non-GBS isolates was undertaken⁵. Suspensions of organisms with optical density equivalent to a 0.5 McFarland standard were made in

deionised water and diluted to produce inoculum ranges of 10⁶-10⁸, 10²-10⁴, and 10¹-10³ cfu/mL for all GBS and species other than GBS (including other *Streptococcus* spp., *Staphylococcus* spp., *Enterobacteriaceae*, *Lactobacilli* and *Candida* spp.) and inoculated onto all media. Plates were incubated according to manufacturers' instructions for up to 48 hours. Plates were read by staff not involved in the development of the project. Additionally, 50 GBS and 50 non-GBS were also inoculated into LIM Broth (Todd Hewitt Broth with colistin and nalidixic acid) and incubated at 36 ± 1°C, overnight, prior to being streaked onto *Brilliance* GBS Agar. Sensitivity for *Brilliance* GBS Agar was calculated based on the presence of correctly coloured pink colonies (presumptive GBS) and any other coloured or colourless colonies were reported. Specificity was calculated based on the number of true negative plates, i.e. the number of plates with correctly coloured non-GBS colonies plus plates with no growth.

Performance of *Brilliance* GBS Agar.

Performance	Inoculum level (cfu/mL)	<i>Brilliance</i> GBS Agar	
		Incubation time (hours)	
		18	24
Sensitivity (%)	10 ⁶ -10 ⁸	98	97
	10 ² -10 ⁴	90	94
	10 ¹ -10 ³	73	76
Specificity (%)	10 ⁴ -10 ⁶	100	100
	10 ² -10 ⁴	100	100
% Inhibition	10 ⁶ -10 ⁸	NC	99.3

NC- Not calculated (% inhibition was calculated at 24 hours).

Performance of *Brilliance* GBS Agar using broth enrichment.

Performance	<i>Brilliance</i> GBS Agar	
	Time (hours)	
	18	24
Sensitivity (%)	100	100
Specificity (%)	92	90

Clinical Performance

Brilliance GBS Agar has been evaluated through a series of external trials conducted at different European laboratories and hospitals, which compared and demonstrated the performance of the device in a clinical setting. The performance characteristics sensitivity, specificity, positive predictive value (PPV), negative predictive value (NPV), and percentage inhibition have been assessed for this medium.

Brilliance[™] GBS Agar proved to be a consistently highly selective media for the isolation of GBS from clinical samples, detecting presumptive GBS within 24 hours. Additionally, it showed noteworthy inhibition of non-target organisms and the ability to detect β- and non-haemolytic GBS.

In one performance study (Trial 1), a total of 920 genital swab samples were taken from women and men during routine screening at Wexham Park Hospital, UK. In a second study (Trial 2), 300 vaginal swabs were taken from pregnant women at 35-38 weeks gestation, following routine screening at University Hospital Brno, Czech Republic. The swabs were inoculated into LIM broth prior to being sub-cultured onto the medium. In a third study (Trial 3) at the Withybush General Hospital, Wales, 250 high vaginal

swabs were taken from pregnant and non-pregnant women and streaked onto the medium. In a study in Basel, Switzerland (Trial 4), 513 vaginal swabs taken from pregnant women undergoing screening were inoculated into LIM enrichment broth before streaking onto the medium. In a study in Koblenz, Germany (Trial 5), 499 vaginal and cervical swabs were collected from pregnant and non-pregnant women undergoing screening for GBS (born between 1960 and 2000).⁵

Performance of *Brilliance* GBS Agar evaluated in clinical trials.

Performance	Brilliance GBS Agar (%)				
	Trial 1*	Trial 2**	Trial 3 ^a	Trial 4**	Trial 5**
Sensitivity	94.4	96.4	94.59	100	97.4
Specificity	99.4	98.0	96.24	81.11	96.7
PPV	96.0	91.4	81.40	70	84.3
NPV	99.1	99.2	99.03	100	99.5
% Inhibition	96	61.7	90.06	49.3	68.3

* Pre-launch trial.

** Post-launch trial.

^a Trial was begun prior to launch (to compliment the pre-launch external trial conducted at Wexham Park Hospital) and continued during and after launch of the product.

Summary of results found in the studies evaluated in a literature review.

Study	Verhoeven et al, 2014 ⁶	Salem & Anderson, 2015 ⁷	Joubrel et al, 2014 ⁸
Sensitivity (%)	100	96	94.3
Specificity (%)	90.3	100	96.2
PPV (%)	68.6	100	ND
NPV (%)	100	98.97	ND

ND – No data

The Summary of Safety and Performance (SSP) for this device will be available in the European database on medical devices (Eudamed) where it is linked to the device's Basic UDI-DI (*Brilliance* GBS Agar (5032384BrillianceGBS37) and Columbia CNA Agar / *Brilliance* GBS Agar (032384CNAGBSAgarNS)). Refer to: <https://ec.europa.eu/tools/eudamed>.

Limitations

Organisms with atypical enzyme patterns may give anomalous reactions on *Brilliance* GBS Agar. Certain strains of *Streptococcus salivarius* may appear as pink colonies.

Samples containing faecal material or blood may cause some localised discolouration within the medium. This discolouration should not be confused with a true chromogenic reaction, where coloured colonies are visible.

Identifications are presumptive and should be confirmed.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Bibliography

- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>
- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from:

<https://www.acog.org/en/Clinical/ClinicalGuidance/CommitteeOpinion/Articles/2020/02/PreventionofGroupBStreptococcalEarly-OnsetDiseaseinNewborns>

- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [cited 2020 May 29]. Available from: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis - Melin Pierrette [Internet]. [cited 2020 May 29]. Available from: <https://orbi.uliege.be/handle/2268/79253>
- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. *Journal of Clinical Microbiology*. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. *Pathology*. 2015 Oct;47(6):580–2.
- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyart C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. *Diagnostic Microbiology and Infectious Disease*. 2014 Dec;80(4):282–4.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use

	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalogue marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC®

© 2019 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, England

2797
For technical assistance please contact your local distributor.

Version	Modifications introduced
1.0	2022-06-13. New document

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Literatura

- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>
- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from: <https://www.acog.org/en/Clinical/ClinicalGuidance/CommitteeOpinion/Articles/2020/02/PreventionofGroupBStreptococcalEarly-OnsetDiseaseinNewborns>
- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [cited 2020 May 29]. Available from: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis - Melin Pierrette [Internet]. [cited 2020 May 29]. Available from: <https://orbi.uliege.be/handle/2268/79253>
- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. *Journal of Clinical Microbiology*. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. *Pathology*. 2015 Oct;47(6):580–2.
- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyart C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. *Diagnostic Microbiology and Infectious Disease*. 2014 Dec;80(4):282–4.

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Bibliografi

- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>
- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from: <https://www.acog.org/en/Clinical/ClinicalGuidance/CommitteeOpinion/Articles/2020/02/PreventionofGroupBStreptococcalEarly-OnsetDiseaseinNewborns>
- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians

- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis - Melin Pierrette [Internet]. [cited 2020 May 29]. Available from: <https://orbi.uliege.be/handle/2268/79253>
- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. *Journal of Clinical Microbiology*. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. *Pathology*. 2015 Oct;47(6):580–2.
- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyart C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. *Diagnostic Microbiology and Infectious Disease*. 2014 Dec;80(4):282–4.

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Bibliographie

- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>
- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from: <https://www.acog.org/en/Clinical/ClinicalGuidance/CommitteeOpinion/Articles/2020/02/PreventionofGroupBStreptococcalEarly-OnsetDiseaseinNewborns>
- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [cited 2020 May 29]. Available from: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis - Melin Pierrette [Internet]. [cited 2020 May 29]. Available from: <https://orbi.uliege.be/handle/2268/79253>
- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et coll. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. *Journal of Clinical Microbiology*. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. *Pathology*. 2015 Oct;47(6):580–2.

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- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaguet M, Poyart C, et coll. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. Diagnostic Microbiology and Infectious Disease. 2014 Dec;80(4):282–4.

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- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline Nr. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [zitiert am 29. Mai 2020]. Abrufbar unter: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis – Melin Pierrette [Internet]. [zitiert am 29. Mai 2020]. Abrufbar unter: <https://orbi.uliege.be/handle/2268/79253>
- (5) Daten in den Akten.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. *Journal of Clinical Microbiology*. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. *Pathology*. 2015 Oct;47(6):580–2.
- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyart C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. *Diagnostic Microbiology and Infectious Disease*. 2014 Dec;80(4):282–4.

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Bibliográfia

- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>
- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from: <https://www.acog.org/en/Clinical/ClinicalGuidance/CommitteeOpinion/Articles/2020/02/PreventionofGroupBStreptococcalEarly-OnsetDiseaseinNewborns>
- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [cited 2020 May 29]. Available from: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis -

- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. *Journal of Clinical Microbiology*. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. *Pathology*. 2015 Oct;47(6):580–2.
- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyart C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. *Diagnostic Microbiology and Infectious Disease*. 2014 Dec;80(4):282–4.

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Bibliografia

- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>
- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from: <https://www.acog.org/en/Clinical/ClinicalGuidance/CommitteeOpinion/Articles/2020/02/PreventionofGroupBStreptococcalEarly-OnsetDiseaseinNewborns>
- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [cited 2020 May 29]. Available from: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis - Melin Pierrette [Internet]. [cited 2020 May 29]. Available from: <https://orbi.uliege.be/handle/2268/79253>
- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. *Journal of Clinical Microbiology*. 2014 Mar 1;52(3):991-3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. *Pathology*. 2015 Oct;47(6):580-2.
- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyart C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. *Diagnostic Microbiology and Infectious Disease*. 2014 Dec;80(4):282-4.

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Literatūra

- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>

- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from: <https://www.acog.org/en/Clinical/ClinicalGuidance/CommitteeOpinion/Articles/2020/02/PreventionofGroupBStreptococcalEarly-OnsetDiseaseinNewborns>
- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [cited 2020 May 29]. Available from: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis - Melin Pierrette [Internet]. [cited 2020 May 29]. Available from: <https://orbi.uliege.be/handle/2268/79253>
- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. Journal of Clinical Microbiology. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. Pathology. 2015 Oct;47(6):580–2.
- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyat C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. Diagnostic Microbiology and Infectious Disease. 2014 Dec;80(4):282–4.

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Bibliografia

- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>
- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from: [https://www.acog.org/en/Clinical/Clinical Guidance/Committee Opinion/Articles/2020/02/Prevention of Group B Streptococcal Early-Onset Disease in Newborns](https://www.acog.org/en/Clinical/Clinical%20Guidance/Committee%20Opinion/Articles/2020/02/Prevention%20of%20Group%20B%20Streptococcal%20Early-Onset%20Disease%20in%20Newborns)
- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [cited 2020 May 29]. Available from: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis - Melin Pierrette [Internet]. [cited 2020 May 29]. Available from: <https://orbi.uliege.be/handle/2268/79253>
- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of *Streptococcus agalactiae* Vaginal Colonization in Pregnant Women. *Journal of Clinical Microbiology*. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B *Streptococcus* from vaginal specimens in pregnant women. *Pathology*. 2015 Oct;47(6):580–2.

- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyart C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. *Diagnostic Microbiology and Infectious Disease*. 2014 Dec;80(4):282–4.

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- (1) SMI B 58: detection of carriage of group B streptococci [Internet]. GOV.UK. [cited 2020 May 29]. Available from: <https://www.gov.uk/government/publications/smi-b-58-processing-swabs-for-group-b-streptococcal-carriage>
- (2) Prevention of Group B Streptococcal Early-Onset Disease in Newborns [Internet]. [cited 2020 May 29]. Available from: <https://www.acog.org/en/Clinical/ClinicalGuidance/CommitteeOpinion/Articles/2020/02/PreventionofGroupBStreptococcalEarly-OnsetDiseaseinNewborns>
- (3) Group B Streptococcal Disease, Early-onset (Green-top Guideline No. 36) [Internet]. Royal College of Obstetricians & Gynaecologists. [cited 2020 May 29]. Available from: <https://www.rcog.org.uk/en/guidelines-research-services/guidelines/gtg36/>
- (4) Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis - Melin Pierrette [Internet]. [cited 2020 May 29]. Available from: <https://orbi.uliege.be/handle/2268/79253>
- (5) Data on file.
- (6) Verhoeven PO, Noyel P, Bonneau J, Carricajo A, Fonsale N, Ros A, et al. Evaluation of the New Brilliance GBS Chromogenic Medium for Screening of Streptococcus agalactiae Vaginal Colonization in Pregnant Women. Journal of Clinical Microbiology. 2014 Mar 1;52(3):991–3.
- (7) Salem N, Anderson JJ. Evaluation of four chromogenic media for the isolation of Group B Streptococcus from vaginal specimens in pregnant women. Pathology. 2015 Oct;47(6):580–2.
- (8) Joubrel C, Gendron N, Dmytruk N, Touak G, Verlaquet M, Poyart C, et al. Comparative evaluation of 5 different selective media for Group B Streptococcus screening in pregnant women. Diagnostic Microbiology and Infectious Disease. 2014 Dec;80(4):282–4.

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Oxoid GmbH,