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TRACHEAL TUBE INTRODUCER **TRACHEAL TUBE GUIDE** **INSTRUCTIONS FOR USE**

These instructions for use refer to the following products:

09-3617, 09-3627, 09-3610, 09-3620, 09-3717, 09-3727, 09-3817, 09-3827, 09-3810, 09-3820, 09-3117, 09-3127, 09-3110, 09-3120, 09-5617, 09-5627, 09-5610, 09-5620, 09-5717, 09-5727, 09-5817, 09-5827, 09-5810, 09-5820, 09-5117, 09-5127, 09-5110, 09-5120, 09-1610, 09-2010, 09-3657, 09-3757, 09-3857, 09-3157, 09-4657, 09-4757, 09-4857, 09-4157, 09-5657, 09-5757, 09-5857, 09-5157, 09-3650, 09-3850, 09-4757, 09-5850, 09-5150.

Pictograms

	catalogue number		manufacturer		free from phthalates (incl. DEHP)
	serial No.		sterilized by ethylene oxide		do not use if packaging is damaged
	date of manufacture		latex-free		do not re-sterilize *
	use by		consult instructions for use		keep dry
	the product complies with the requirements of the European Union/Number of the Notified Body participating in the device conformity assessment.		caution: read all warnings and precautions in instructions for use		keep away from sunlight
			do not re-use *		store at 5-40 °C

*) doesn't apply to reusable guides and introducers

Description

1,6 mm for 2,0 – 3,0 intubation tubes 2,0 mm for 2,5 – 4,0 intubation tubes
3,3 mm for 4,0 – 5,5 intubation tubes 5,0 mm for 6,0 – 10,0 intubation tubes
Lengths: 600mm, 700mm, 800 mm and 1000mm

They can be used again after washing and sterilization, maximum 4 times (relates to the multi-use guides).

In one version, the guides have a soft distal end.

Indications

A difficult intubation

Replacement of the intubation tubes

Contraindications

The guide is not recommended for intubation through nose or upon changes of trachea.

Instructions for use

Upon difficult intubations:

- To gently remove the sterile tracheal tube from the packaging and check before use.
- If introducers with soft tip are used, check the durability of joint of the soft tip and the introducer.
- To wash the guide with a sterile solution of water-soluble lubricant right before use.
- To connect the patient to oxygen and perform endoscopy upon observance of current medical techniques.
- If there is such a necessity, to tentatively shape the guide in order to get an optimal flexure.
- To keep particular caution.
- To insert the guide into the throat and next manipulate so as to insert it into the larynx.
- It is advisable to have confirmed correct arrangement of the guide by correct laryngospasms.
- To leave an endoscope in the patient's throat, and next put a tube on the guide and carefully insert the tube through the larynx into the trachea.
- To pay attention to observance of current approved medical techniques during the intubation.
- The marker of depth of intubation may be used as assistance upon correction of arrangement of the tube.
- Keep the tube and next carefully remove the guide.

Upon replacement of an intubation tube:

- To gently remove the sterile guide from the packaging and check before use.
- If introducers with soft tip are used, check the durability of joint of the soft tip and the introducer.
- To observe approved medical techniques with regard to correct size of the devices, intubation and extubation.
- To select a point of reference, for instance teeth or nostrils and mark the depth of insertion of intubation tube.

- To insert the guide and to slip it until a moment, when a depth marker of the guide comes level with a depth marker of the tube.
- Keep the guide and remove the intubation tube.
- Take advantage of a relevant marker in order to check arrangement of the tube's guide.
- Place a replaceable intubation tube on the tube's guide and insert deeply until a depth marker of the tube comes level with a depth marker of the guide. Both of them are to be correctly inserted taking into consideration relevant limiters.
- While airing, check once more, whether the new intubation tube is arranged correctly in the trachea.
- Remove the guide from the intubation tube.

Washing, disinfection and sterilization of the multi-use guides

- Insert the used guide into the disinfectant.
- Wash the guide carefully using water solution of soap, removing precisely all visible impurities.
- In order to remove impurities, wash the product in distilled water with a disinfecting agent.
- Dry the guide carefully.
- Check whether the product is not damaged, throw aside a damaged product.
- Pack up the product in new packaging of paper – foil type.
- Close the packaging tightly.
- Carry out sterilization using ethylene oxide with application of relevant sterilization procedures.
- During intervals between procedures the product should be kept in a sterile condition, in packaging of paper – foil type, in dry rooms, upon 5-40°C temperature, in a manner protecting against direct action of solar radiation, vapours of acids and organic solvents as well as other harmful agents.
- The SUMI company does not bear responsibility for the re-sterilization process.

Precautions

A physician should take proper precautions in order to avoid potential damage of an intubation tube with a guide during intubation or during replacement of a tube.

Warnings

Owing to a danger of respiratory tract infection, mainly by treatment resistant hospital bacterial strains, one should observe current approved medical procedures and techniques. Repeated use of a guide designed for single use may constitute a danger to the patient. Observance of the current approved medical techniques and performing of procedures by skilled medical personnel reduces to a minimum possibility of occurrence of complications e.g.: excoriations of epithelium and ulcerations of lips, mouth, gullet, vocal folds, injuries of lips, tongue, gullet, larynx, vocal apparatus of larynx, palate, tonsils and teeth, paresis of hypoglossal or glossal nerves, perforation of gullet, perforation of larynx, bleedings from throat, larynx, strictures of larynx, trachea.

Notes

- Do not use in case of damage of packaging.
- Do not use after expiry date.
- The device is not intended for reesterilization – doesn't apply to reusable guides and introducers

Responsibility

The guides of SUMI make undergo physical, biological and chemical investigations. It allows their safe use and top quality. There is recommended making familiar with the instructions for use before using the guides for difficult intubations. The tracheal tubes manufactured by SUMI may be used only by qualified medical personnel. SUMI shall not bear responsibility for consequences caused by inappropriate selection of the size or using the tracheal tubes contrary to their use.

In the event of a complaint or an occurrence of a medical incident, keep the medical device with its immediate packaging. Photographic documentation is also recommended.

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