**FRESENIUS
KABI**

Fresenius Kabi AG

61346 Bad Homburg

Germany

Tel. +49 (0) 61 72 / 686-0

www.fresenius-kabi.com

CE0123

STERILE EO

REF 7755645

de - Freka® FKJ-Set CH/FR 9, ENFit.

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Components

(1 x 1 piece)

- long cannula (white grip) with blunt steel mandrin for catheterisation of the jejunum
- short cannula (blue grip) for puncturing the abdominal wall
- Freka® polyurethane tube
 - FR 9, length 75 cm
 - with distal opening
 - internal diameter 1.9 mm
 - external diameter 2.9 mm
 - with radio-opaque strip and
 - length markings
- Freka® ENFit Connector for FCJ consisting of:
 - fixation screw
 - ENFit connector
- Freka® fixation plate made of silicone rubber
 - with tube clamp for fixation
 - with 2 small holes for sewing onto the abdominal wall
- Freka® tube clamp

Table of contents

1. Applications
2. Contraindications
3. Information for use
4. Aftercare instruction for the puncture site
5. Care instructions of the Feeding Tube
6. Application of medicine
7. Application of nutrition
8. Duration of use
9. Removing the Feeding Tube
10. Order information

1. Applications

Early post-operative long-term intrajejunal feeding after

- Laparotomy or
- Laparoscopy.

2. Contraindications

Absolute contraindications

- Distal GI tract obstruction
- Intestinal atony
- Peritonitis
- Blood clotting disorders
- Sepsis

Relative contraindications

- Immune suppression
- Radiation enteritis
- Pancreatitis
- Crohn's disease (risk of fistulation)
- Ascites
- Peritoneal carcinoma

3. Information for use

- After concluding the operation, select the ideal point of entry for the tube into the intestine and the transcutaneous point of entry for the tube in the upper left side of the abdomen.

1) Laparotomy

- with splitting cannula
- penetrate the abdominal wall diagonally with the blue cannula from cranially to caudally. (Fig. 1).
- The tube is guided into the gastric cavity from the outside inwards through the split cannula.
- Then withdraw the split cannula (Fig. 2).
- Puncture in the opposite direction of mesentery into the wall of the intestine beyond the ligament of treitz (the tip of the tube should be remote from the duodenojejunal flexure).
- Stretch the jejunal loop and puncture diagonally into the submucosa, use the white cannula with the withdrawn mandrin (Fig. 3).
- After extending the blunt mandrin, tunneling approx. 5-10 cm between the layers of the intestinal wall (Fig. 4).
- Remove the mandrin and use the cannula to puncture the lumen of the intestine (Fig. 5).
- Push the tube forwards through the cannula into the lumen of the intestine and slide forwards until it is approx. 30 cm into the lumen, with digital control (Fig. 6).
- Withdraw the split cannula through the tube and remove after splitting (Fig. 7).
- Using dissolvable purse-string sutures, which are first knotted around the puncture point as invaginated sutures and which serve as a temporary fastening at the parietal peritoneum with the same needle, seal off the entry point of the tube into the intestinal wall and attach the intestinal loop to the parietal peritoneum (Fig. 8).
- The intestine should also be sutured to the peritoneum a few centimetres above and below the point of entry of the tube in order to prevent dislocation due to rotation of the intestine.
- Insert the tube into the guidance of the fixation plate and fix it with the clamp.
- Secure the position of the tube by tacking the fixing plate on the abdominal wall about 1-2 cm from the entry point

(Fig. 9).

- After attaching the tube clamp, first push the fixing screw over the tube.
- Then push the metal pin of the ENFit connection as far as possible into the tube and secure it with the fixing screw (Fig. 9).
- Pull down the screw aid (external grooved ring) and remove it.

IMPORTANT NOTE

The correct position of the tube must be checked by X-ray.

2) Laparotomy without splitting cannula

- Puncture in the opposite direction of the mesentery into the wall of the intestine beyond the ligament of treitz (the tip of the tube should be remote from the duodenojejunal flexure).
- Stretch the jejunal loop and puncture diagonally into the submucosa, use the white cannula with the withdrawn mandrin (Fig. 10).
- After extending the blunt mandrin, tunneling approx. 5-10 cm between the layers of the intestinal wall (Fig. 11).
- Remove the mandrin and use the cannula to puncture the lumen of the intestine (Fig. 12).
- Push the tube forwards through the cannula into the lumen of the intestine and slide forwards until it is approx. 30cm into the lumen, with digital control (Fig. 13).
- Withdraw the cannula over the tube and remove it.
- Penetrate through the abdominal wall diagonally with the blue cannula from cranially to caudally (Fig. 14).
- Carefully guide the tube from the inside outwards through the cannula, past the sharp-edged tip of the cannula (Fig. 15).
- Withdraw the cannula and remove.
- Using dissolvable purse-string sutures, which are first knotted around the puncture point as invaginated sutures and which serve as a temporary fastening at the parietal peritoneum with the same needle, seal off the entry point of the tube into the intestinal wall and attach the intestinal loop to the parietal peritoneum (Fig. 16).
- The intestine should also be sutured to the peritoneum a few centimetres above and below the point of entry of the probe in order to prevent dislocation due to rotation of the intestine.
- Insert the tube into the guidance of the fixation plate and fix it with the clamp.
- Secure the position of the tube by tacking the fixing plate on the abdominal wall about 1-2 cm from the entry point (Fig. 17).
- After attaching the quick-release clamp, first push the fixing screw over the tube.

- Then push the metal pin of the ENFit connection as far as possible into the tube and secure with the fixing screw (Fig. 17).
- Pull down the screw aid (external grooved ring) and remove it.

IMPORTANT NOTE

The correct position of the tube must be checked by X-ray.

3) Laparoscopy

- Fix the intestinal loop to the abdominal wall at two places (Fig. 18).
- Puncture the abdominal wall with the longer (white) cannula (Fig. 19).
- Then puncture the intestine.
- Push the tube forward through the cannula approx. 30 cm.
- Withdraw the cannula (Fig. 20).
- Attach the intestinal loop to the abdominal wall with dissolvable purse-string sutures until the tube is no longer visible (Fig. 21).
- Insert the tube into the guidance of the fixation plate and fix it with the clamp.
- Secure the position of the tube by tacking the fixing plate on the abdominal wall, (1-2 cm from the entry point) (Fig. 22).
- After attaching the tube clamp, first push the fixing screw over the tube.
- Then push the metal pin of the ENFit connection as far as possible into the tube and secure with the fixing screw (Fig. 22).
- Pull down the screw aid (external grooved ring) and remove it.

IMPORTANT NOTE

The correct position of the tube must be checked by X-ray.

4. Aftercare instruction for the puncture site

- Use an applicator and a disinfectant to clean the puncture site underneath the fixation plate and dry well without removing the fixation plate.
- It is recommended to mark the length of the tube at the exit point to identify dislocation at an early stage.
- Cover the puncture site with a compress and plaster.

In the first week after insertion the puncture site should be checked by a physician daily.

IMPORTANT NOTE

The threads of the fixing plate must not be removed or the tube will not be fixed in place.

5. Care instructions for the Feeding Tube**IMPORTANT NOTE**

The tube must not be mobilized in the stoma. The risk of embedding does not exist with the FCJ.

The tube must be carefully flushed with 20 ml of water before and after each administration of nutrition, preferably using a 20 ml ENFit syringe.

IMPORTANT NOTE

If the tube blocks, the lumen of the tube should be under no circumstances cleared by force small-volume ENFit syringe, e.g. 2,5 or 5 ml) or by a guide wire. Otherwise there is a risk of perforating the tube.

Do not use fluids containing acids, particularly fruit teas or fruit juices, in conjunction with tube feed as they can cause coagulation which may clog the tube.

6. Application of medicines

Medication should preferably be administered via the ENFit connector of the feeding tube or the medication port of the giving set in dissolved form. Preference should be given to liquid medication (drops or syrup). If tablets and capsules are to be administered, they must be crushed very fine and dissolved in water, provided that the effect of the medication will not change.

If in case of any doubt, consult the physician in charge or a pharmacist.

IMPORTANT NOTE

The tube must be flushed before and after each administration of medication using 20 ml lukewarm water preferably with a 20 ml ENFit syringe. Medication should be administered separately and not mixed with nutrition.

7. Application of nutrition

It is possible in technical terms to immediately start the administration of nutrition to the small intestine after the fine-needle catheter jejunostomy. In any case, the dosage of the nutrition should be increased gradually. Commercially available enteral feeding pumps should be used for the recommended pump-controlled continuous tube feeding.

IMPORTANT NOTE

Under no circumstances should an infusion pump designed for parenteral application be used (risk of confusion).

8. Duration of use

In the past, the jejunal tube has been recorded as having been in place for several months without complications.

9. Removing the Feeding Tube

- Remove the fixing threads of the fixing plate
- Remove tube carefully
- Apply plaster dressing
- Observe a 12-hour fasting period.

The Freka® FCJ can usually be removed without a jejunocutaneous fistula being left behind.

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