# Puritan Bennett™

## D/Flex Disposable Inspiratory Filter

<b>en</b>	Instructions For Use
<b>fr</b>	Mode d'emploi
<b>de</b>	Gebrauchsanweisung
<b>nl</b>	Gebruiksaanwijzing
<b>it</b>	Istruzioni per l'uso
<b>es</b>	Instrucciones de uso
<b>sv</b>	Bruksanvisning
<b>da</b>	Brugsanvisning
<b>no</b>	Bruksanvisning
<b>fi</b>	Käyttöohjeet
<b>pt</b>	Instruções de uso
<b>ru</b>	Инструкция по применению
<b>zh</b>	使用说明
<b>pl</b>	Instrukcja użytkowania
<b>cs</b>	Návod k použití
<b>sk</b>	Návod na použitie
<b>sl</b>	Navodilo za uporabo
<b>hu</b>	Használati útmutató
<b>el</b>	Οδηγίες χρήσης
<b>tr</b>	Kullanma Talimatı
<b>bg</b>	Инструкции за употреба

<b>ro</b>	Instrucțiuni de folosire
<b>ko</b>	사용 방법
<b>hr</b>	Upute za uporabu
<b>sr</b>	Uputstva za upotrebu

## D/Flex Disposable Inspiratory Filter

↑ Indicates flow direction; filter is designed for use only in this orientation.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

### Description

The disposable inspiratory filter minimizes contaminants delivered to the patient and reduces the possibility of contaminants reentering the ventilator.

### Use

*D/Flex* filters are ready to use without sterilization. Handle with care; abusive handling can damage filter medium. *D/Flex* filters are intended for single patient use; no sterilizing or cleaning is necessary. Remove and discard filter after each patient use, or if visibly damaged.

### WARNINGS

- **Ensure all patient connections are leak-tight. Perform patient circuit leak check every time the D/Flex disposable inspiratory filter is used. Consult your ventilator's operators manual.**
- **Discard the D/Flex filter after use and follow appropriate safety practices. Do not sterilize. Ensure all patient circuit connections are leak-tight. Perform a patient circuit leak check each time filter is installed.**
- **Soaking, rinsing, or sterilization of this filter may result in retention of harmful residuals in the device and/or increasing air flow resistance.**
- **Do not reuse disposable or non-reusable breathing circuit components or devices with the D/Flex filter. The filter does not prevent cross-contamination of circuit components. Follow your institution's normal operating procedures when using disposable or non-reusable circuit components with the D/Flex filter.**

### CAUTIONS

- To ensure proper installation, observe flow direction markings on the filter.
- To avoid rapid increase in flow resistance, do not use on exhalation limb of patient circuit, or downstream of a humidifier or nebulizer.

## NOTE

- Verify filter flow resistance prior to use. For Puritan Bennett™ ventilator applications, discontinue use if flow resistance is greater than 4 cmH<sub>2</sub>O at 100 Lpm.
- Each filter is non-sterile and assembled in accordance with Class 100,000 clean room specifications.

## Disposal

Dispose of items according to your institution's protocol. Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

## Specifications

<i>Viral/bacterial filtration efficiency:</i>	Greater than 99.999%
<i>Tidal Volume:</i>	2– 2500 ml
<i>Filtration Efficiency:</i>	≥99% @ 30 l/min, ≥99% @ 15 l/min
<i>Pressure Drop:</i>	< 1.0 cm H <sub>2</sub> O @ 30 l/min, < 0.5 cm H <sub>2</sub> O @ 15 l/min
<i>Size:</i>	10.01 cm long x 7.32 cm wide x 5.97 cm deep (4.0 in. x 3.0 in. x 2.4 in.)
<i>Internal volume:</i>	Approximately 136 mL.
<i>Leakage:</i>	<10 ml/min @ 70 cm H <sub>2</sub> O
<i>Connectors:</i>	The 22 mm connectors meet requirements of the International Standards Organizations (ISO) 5356-1 (Anesthetic and respiratory equipment - conical connectors).

## Ordering Information

REF	Description
4-074601-00	D/Flex filter, 22 mm ISO connectors (carton of 12)

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© 2011 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Inc, 15 Hampshire Street,
Covidien Ireland Limited,