PL

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pro

medica

TORUŃ

BONEWAX®

BONEWAX – sterylny воск kostny

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BONEWAX® - registered trade mark by
Promedica Toruń, Poland

Healthium Medtech Limited
#472-D, 13th Cross, 4th Phase,
Peenya Industrial Area, Bangalore,
Karnataka-560058, India
Email: care@healthiummedtech.com
Mfg. Lic. No.: MFG/MD/2021/000367

MEDDEVICES LIFESCIENCES B.V
Krajenhoffstraat 137 A, 1018RG
Amsterdam, Netherlands
Email: info@meddevices.net
Phone: +31-202254558.

ENG

INSTRUCTIONS FOR USE

pro

medica

TORUŃ

BONEWAX®

BONEWAX – sterile bone wax

Description
Bonewax is sterile bone wax, a mixture of white beeswax hard paraffin and isopropyle palmitate softening agent. Bonewax is opaque and has a waxy odour.

Indications:
Bonewax may be used to control bleeding from bone surfaces.

Action:
Bonewax helps to achieve local bone hemostasis by acting as a mechanical (tamponade) barrier. It does not act biochemically.

Contraindications:
Bonewax should not be used where rapid osseous regeneration and fusion are desired.

Warnings:
Bonewax may inhibit osteogenesis and may act as a physical barrier to the reparative process. Do not use if the package is opened or damaged. Excess bone wax should be removed from the operative site and put inside "Biomedical waste" and suitably disposed.

Precautions:
Bonewax should be used sparingly. The excess of Bonewax should be removed from the operative site. The package should be opened just prior to use to minimise the possibility of contamination and excessive drying.

Adverse reactions:
Mild inflammatory reactions have been reported in tissues immediately adjacent to the site of implantation. Studies have suggested that bone wax as a foreign body may impair the ability of cancellous bone to clear bacteria. In animal models, the local accumulation of foreign body giant cells has been observed and histological examination has revealed the appearance of macrophages and occasionally polymorphonuclear leukocytes and lymphocytes.

Dosage and administration:
Bonewax should be used immediately after removal from the package. Using aseptic technique, Bonewax should be warmed to desired consistency by manipulation with the fingers or by immersion of the unopened foil packet in a warm sterile solution.

Caution:
Single-use product. The product should not be resterilised. Discard unused bone wax.

Storage conditions and shelf-life:
The product should be stored in original packaging , under controlled temperature and moisture conditions, and protected from direct sunlight. Storage temperature below 30°C is recommended. The package should be protected against dust, moistening or damaging. The expiry date is on the labelling. The product should not be used after the expiry date.

Sterility:
In order to keep the product sterile and to retain its adequate properties it should be stored in appropriate controlled conditions provided for on the labelling. The product is sterile until the package is opened or damaged.

Presentation and contents:
The product is double-packed (inner packaging: paper pouch with coating of aluminium foil, outer packaging: paper - film pouch). Bonewax is available sterile in individual packaging, each containing 2.5 grams of wax.

Symbols used on labelling:			
	Do not reuse		Sterilisation: radiation
	Do not resterilise		Manufacturer
	Do not use if the package is damaged		Manufacturing date, year/month
	Maximum temperature		Valid until
	Protect from moisture		Batch number
	Protect from direct sunlight		Reference number
	Warning		2265 Notified Body No
	See Instructions for Use		Unique device identifier
	Medical Device		GS 1 Data matrix barcode
	Double sterile barrier system		Importer

Rev.00 Date:08/2022

Promedica Toruń Sp. z o.o.,
87-100 Toruń, POLAND
Tel./Fax: +48 56 62 30 112
www.promedica.torun.pl

BONEWAX® - registered trade mark by
Promedica Toruń, Poland

Healthium Medtech Limited
#472-D, 13th Cross, 4th Phase,
Peenya Industrial Area, Bangalore,
Karnataka-560058, India
Email: care@healthiummedtech.com
Mfg. Lic. No.: MFG/MD/2021/000367

MEDDEVICES LIFESCIENCES B.V
Krajenhoffstraat 137 A, 1018RG
Amsterdam, Netherlands
Email: info@meddevices.net
Phone: +31-202254558.