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NASAL SPONGE

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NASAL SPONGE

Nasal dressing

haemostatic nasal sponge, sterile, absorbent, single use

Description:
The Nasal Sponge nasal dressing is mainly composed of PVA expanding sponge with a hemostatic gauze cover with a string or without the cover and the string, and with the string and a tube that makes breathing easier.

Product properties:
When in contact Following the contact with water, the gauze cover can form a viscous gel and quickly stop capillary bleeding. In the meantime, the expanding sponge provides a controlled pressure on the bleeding site. This double function allows to achieve hemostasis easily.

Indications:
Nasal Sponge is used for epistaxis and postoperative hemostasis of the nasal cavity and paranasal sinus.

Contraindications:
Do not use in patients allergic to any product components.

Precautions:
The package should be checked before use. After the package is opened, the product should be kept in a dry place. The product should not be moistened with water or water-bearing liquids before use. Gloves and instruments must be dry for the product application.
When the product is used in ENT surgery, the string must be fixed to the wing of the nose or the cheek to prevent it from being aspirated into the respiratory tract. The product should not be kept on the bleeding site for more than 72 hours. The product cannot be reused.
Do not use if the package is damaged. The Nasal Sponge has to be used immediately after the package is opened.
The product may be disposed of by the same methods as other medical waste.

Adverse effects:
In rare instances the physiochemical conditions associated with nasal surgery, both with and without nasal packing, may present a risk of toxic shock syndrome (TSS). The warning signs of TSS include sudden fever, vomiting, diarrhea, dizziness, fainting (or near fainting when standing up) and/or rash resembling a sunburn. Nasal dressings exhibit no antimicrobial properties, they are not bacteriostatic towards pre-existing infections, and they are not recommended without appropriate treatment being instituted.
Rarely, displacement and ingestion of a nasal dressing may occur. The nasal tampon should be inserted in the nasal cavity in the manner that it cannot move to the nasopharyngeal cavity, as this may cause the patient to choke. The tampon is non-toxic. The patient should be monitored in the normal manner for the ingestion of any non-toxic material. If aspiration should occur, prompt medical intervention is imperative with close patient monitoring and airway support while removing the material.

Dosage and administration:
1. Select a proper type of nasal dressing and unpack the product, then take it out;
2. insert the nasal dressing into the part required along the nasal meatus;
3. inject a proper amount of sterilised water, saline or antibiotic solution into the gel expansive sponge with an injector to make the sponge expand fully;
4. fix the string to the wing of the nose;
5. take the product out after 24 to 48 hours provided there is no active hemorrhage. If not, inject, sterilized water or physiological saline into the nasal cavity to saturate the sponge.
Several minutes later, gently pull the string and take the product out.

Storage:
Make sure to store in a shady, well-ventilated and clean room without corrosive gas. The preferred storage temperature should be between 5°C and 45°C.

Shelf life:
Expiry date is given on the product label. The product should not be used after the expiry date.

Sterility:
In order to keep the product sterile and to retain its adequate properties it should be stored in appropriate controlled conditions provided for on the labelling. The product is sterile until the package is opened or damaged. THE PRODUCT SHOULD NOT BE RSTERILISED!

Presentation and contents:
The product is double-packed (inner packaging: paper pouch with coating of aluminium foil, outer packaging: paper - film pouch).

Symbols used on labelling:

Do not reuse

Do not reterilise

Do not use if the package is damaged

Maximum temperature limited from 5°C to 45°C

Protect from moisture

Protect from direct sunlight

Warning

See Instructions for Use

Single sterile barrier system

Sterilisation: radiation

Manufacturer

Manufacturing date, year/month

Valid until

Batch number

Reference number

1434 Notified Body No

Medical Device

Unique Device Identifier

Double sterile barrier system

Rev: 09. 2022

Product ref No	Product name	Product description	Product size
NS100251SGS	Nasal dressing	with gauze and string	100×25×15mm
NS100201SGS	Nasal dressing	with gauze and string	100×20×15mm
NS1001510GS	Nasal dressing	with gauze and string	100×15×10mm
NS802515GS	Nasal dressing	with gauze and string	80×25×15mm
NS802015GS	Nasal dressing	with gauze and string	80×20×15mm
NS801510GS	Nasal dressing	with gauze and string	80×15×10mm
NS552515GS	Nasal dressing	with gauze and string	55×25×15mm
NS552015GS	Nasal dressing	with gauze and string	55×20×15mm
NS551510GS	Nasal dressing	with gauze and string	55×15×10mm
NS351209GS	Nasal dressing	with gauze and string	35x12x9mm
NS351212GS	Nasal dressing	with gauze and string	35x12x12mm

Product ref No	Product name	Product description	Product size
NS100251SS	Nasal dressing	with string, no gauze	100×25×15mm
NS100201SS	Nasal dressing	with string, no gauze	100×20×15mm
NS1001510S	Nasal dressing	with string, no gauze	100×15×10mm
NS80251SS	Nasal dressing	with string, no gauze	80×25×15mm
NS80201SS	Nasal dressing	with string, no gauze	80×20×15mm
NS801510S	Nasal dressing	with string, no gauze	80×15×10mm
NS55251SS	Nasal dressing	with string, no gauze	55×25×15mm
NS55201SS	Nasal dressing	with string, no gauze	55×20×15mm
NS551510S	Nasal dressing	with string, no gauze	55×15×10mm
NS351209S	Nasal dressing	with string, no gauze	35x12x9mm
NS351212S	Nasal dressing	with string, no gauze	35x12x12mm
NS452015S	Nasal dressing	with string, no gauze	45x20x15mm
NS803015S	Nasal dressing	with string, no gauze	80x30x15mm

Product ref No	Product name	Product description	Product size
NS100251S	Nasal dressing	without gauze and string	100×25×15mm
NS100201S	Nasal dressing	without gauze and string	100×20×15mm
NS1001510	Nasal dressing	without gauze and string	100×15×10mm
NS80251S	Nasal dressing	without gauze and string	80×25×15mm
NS80201S	Nasal dressing	without gauze and string	80×20×15mm
NS801510	Nasal dressing	without gauze and string	80×15×10mm
NS55251S	Nasal dressing	without gauze and string	55×25×15mm
NS55201S	Nasal dressing	without gauze and string	55×20×15mm
NS551510	Nasal dressing	without gauze and string	55×15×10mm
NS 351209	Nasal dressing	without gauze and string	35x12x9mm
NS 351212	Nasal dressing	without gauze and string	35x12x12mm

Product ref No	Product name	Product description	Product size
NS80201SSA	Nasal dressing	with string and airway, no gauze	80x20x15mm
NS45201SSA	Nasal dressing	with string and airway, no gauze	45x20x15mm

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