# SURGICEL™ ORIGINAL and SURGICEL™ NU-KNIT

## Absorbable Haemostats (*oxidised regenerated cellulose*)

SURGICEL™ ORIGINAL and SURGICEL™ NU-KNIT  
Absorbable Haemostats (*oxidised regenerated cellulose*)

en

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## Absorbable Haemostats (*oxidised regenerated cellulose*)

### DESCRIPTION

SURGICEL™ Absorbable Haemostat is a sterile absorbable knitted fabric prepared by the controlled oxidation of regenerated cellulose. SURGICEL™ Absorbable Haemostat is 100% composed of oxidised regenerated cellulose. The fabric is white with a pale-yellow cast and has a faint, caramel-like aroma. It is strong and can be sutured or cut without fraying. A slight discolouration may occur with age, but this does not affect performance.

SURGICEL™ Absorbable Haemostat is intended for use only by healthcare professionals who are trained in the surgical procedures and techniques requiring the use of an Absorbable Haemostat.

The clinical benefit to be expected is the control of capillary, venous, and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective.

A summary of safety and clinical performance can be found at the following link (upon activation): <https://ec.europa.eu/tools/eudamed>.

### INDICATIONS / INTENDED USE

SURGICEL™ Haemostat is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective.

SURGICEL™ Haemostat can be cut to size for use in endoscopic procedures (see Figures 1, 2A, 2B and 3).

*Figure 1. SURGICEL™ Haemostat should be cut to the appropriate size for endoscopic placement. Standard endoscopic procedures should be used up to the point of placement of the absorbable haemostat. Grasp the SURGICEL™ Haemostat at one corner.*

*Figures 2A and 2B. Slowly push the grasping instrument and material into the cavity.*

*Figure 3. With the use of grasping instruments in a second and/or third auxiliary site, placement can be made, and the material repositioned as needed.*

### CONTRAINDICATIONS

- Although packing or wadding sometimes is medically necessary, SURGICEL™ Haemostat should not be used in this manner, unless it is to be removed after haemostasis is achieved.
- SURGICEL™ Haemostat should not be used for implantation in bone defects, such as fractures, since there is a possibility of interference with callus formation and a theoretical chance of cyst formation.
- When SURGICEL™ Haemostat is used to help achieve haemostasis in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, or the optic nerve and chiasm, it must always be removed after haemostasis is achieved since it will swell and could exert unwanted pressure.
- SURGICEL™ Haemostat should not be used to control haemorrhage from large arteries.
- SURGICEL™ Haemostat should not be used on non-haemorrhagic serous oozing surfaces, since body fluids other than whole blood, such as serum, do not react with SURGICEL™ Haemostat to produce satisfactory haemostatic effect.
- SURGICEL™ Haemostat should not be used as an adhesion prevention product.

## WARNINGS

- SURGICEL™ Haemostat is supplied sterile and as the material is not compatible with autoclaving or ethylene oxide sterilization, SURGICEL™ Haemostat should not be resterilized.
- SURGICEL™ Haemostat is not intended as a substitute for careful surgery and the proper use of sutures and ligatures.
- Closing SURGICEL™ Haemostat in a contaminated wound may lead to complications and should be avoided.
- The haemostatic effect of SURGICEL™ Haemostat is greater when it is applied dry; therefore, it should not be moistened with water or saline.
- SURGICEL™ Haemostat should not be impregnated with anti-infective agents or with other materials such as buffering or haemostatic substances. Its haemostatic effect is not enhanced by the addition of thrombin, the activity of which is destroyed by the low pH of the product.
- Although SURGICEL™ Haemostat may be left *in situ* when necessary, it is advisable to remove it once haemostasis is achieved. It must **always** be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm and in proximity to tubular structures that could become constricted by swelling, regardless of the type of surgical procedure because SURGICEL™ Haemostat, by swelling, may exert pressure resulting in paralysis and/or nerve damage. Dislodgement of SURGICEL™ Haemostat could possibly occur by means such as repacking, further intraoperative manipulation, lavage, exaggerated respiration, etc. There have been reports that in procedures such as lobectomy, laminectomy and repair of a frontal skull fracture and lacerated lobe that SURGICEL™ Haemostat, when left in the patient after closure, migrated from the site of application into foramina in bone around the spinal cord resulting in paralysis and, in another case, the left orbit of the eye, causing blindness. While these reports cannot be confirmed, special care must be taken by physicians, **regardless of the type of surgical procedure**, to consider the advisability of removing SURGICEL™ Haemostat after haemostasis is achieved.
- Although SURGICEL™ Haemostat is bactericidal against a wide range of pathogenic microorganisms, it is not intended as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent post-operative infections.
- Discard unintentionally opened devices, partially used devices, partially used packages, used devices and packaging.

- In the event of a product malfunction before use such as product damage, the product should be discarded and a new one retrieved to complete the procedure. In the event of a product malfunction during use, continue or discontinue use at healthcare professional's discretion. In case the product is to be disposed, disposal of the device and packaging should be according to the user's facility policies and standard operating procedures concerning biohazardous materials and waste.
- DO NOT REUSE. Reuse of this device (or portion of this device) may create a risk of product degradation, which may result in device failure and/or cross contamination, which may lead to infection or transmission of blood-borne pathogens to the patients and anyone coming in contact with the device.

## PRECAUTIONS

- Use only as much SURGICEL™ Haemostat as is necessary for haemostasis, holding it firmly in place until bleeding stops. Remove any excess before surgical closure in order to facilitate absorption and minimize the possibility of foreign body reaction, such as encapsulation of the product, which may mimic artifacts on radiographic images, resulting in diagnostic errors and possible reoperation.
- In urological procedures, minimal amounts of SURGICEL™ Haemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.
- Since absorption of SURGICEL™ Haemostat could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals.
- If SURGICEL™ Haemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped.
- Precautions should be taken in otorhinolaryngologic surgery to assure that none of the material is aspirated by the patient. (Examples: controlling haemorrhage after tonsillectomy and controlling epistaxis).
- Care should be taken not to apply SURGICEL™ Haemostat too tightly when it is used as a wrap during vascular surgery.
- **DENTAL:** SURGICEL™ Haemostat should be applied loosely against the bleeding surface.

Wadding or packing should be avoided, especially within rigid cavities, where swelling may interfere with normal function or possibly cause necrosis.

## ADVERSE REACTIONS / UNDESIRABLE SIDE EFFECTS

- “Encapsulation” of fluid and foreign body reactions has been reported.
- There have been reports of stenotic effect when SURGICEL™ Haemostat has been applied as a wrap during vascular surgery. Although it has not been established that the stenosis was directly related to the use of SURGICEL™ Haemostat, it is important to be cautious and avoid applying the material tightly as a wrapping.
- Paralysis and nerve damage have been reported when SURGICEL™ Haemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. While most of these reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures. Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when SURGICEL™ Haemostat was placed in the anterior cranial fossa.
- Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy have been reported. There has been one report of a blocked ureter after kidney resection, in which post-operative catheterization was required.
- Occasional reports of “burning” and “stinging” sensations and sneezing when SURGICEL™ Haemostat has been used as packing in epistaxis, are believed to be due to the low pH of the product.
- Burning has been reported when SURGICEL™ Haemostat was applied after nasal polyp removal and after haemorrhoidectomy. Headache, burning, stinging, and sneezing in epistaxis and other rhinological procedures have been reported. Also reports of stinging when SURGICEL™ Haemostat was applied on surface wounds (varicose ulcerations, dermabrasions, and donor sites) have been made.
- Healthcare professionals should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient and advise the patient to contact a healthcare professional in case of any deviation from the normal post-operative course.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

## MAGNETIC RESONANCE IMAGING (MRI) / CARCINOGENIC, MUTAGENIC, TOXIC TO REPRODUCTION (CMR) / ENDOCRINE DISRUPTING (ED) SAFETY INFORMATION

SURGICEL™ Haemostat is Magnetic Resonance (MR) Safe. No known CMR Category 1a/1b and ED substances are present at >0.1%. Category 1a/1b are defined as known or presumed human carcinogen (H340), mutagen (H350) or reproductive toxicant (H360) based on human evidence and animal studies.

## DOSAGE AND ADMINISTRATION

Sterile technique should be observed in removing SURGICEL™ Haemostat from its sterile container.

Minimal amounts of SURGICEL™ Haemostat in appropriate size are laid on the bleeding site or held firmly against the tissues until haemostasis is obtained.

The amount required depends on the nature and intensity of the haemorrhage to be stopped. The haemostatic effect of SURGICEL™ Haemostat is particularly pronounced when used dry. Moistening the material with water or physiological saline solution is not recommended.

### Disposal

Discard unintentionally opened / partially used / used devices and packaging according to your facility's policies and procedures concerning biohazardous materials and waste.

### Patient Information

The healthcare professional or institution should provide the “Information to be Conveyed to the Patient” to the patient and if required, print the document(s).

### Performance / Actions

After SURGICEL™ Haemostat has been saturated with blood, it swells into a brownish or black gelatinous mass, which aids in the formation of a clot, thereby serving as a haemostatic adjunct in the control of local haemorrhage. When used properly in minimal amounts, SURGICEL™ Haemostat is absorbed from the sites of implantation with practically no tissue reaction. Absorption depends upon several factors including the amount used, degree of saturation with blood, and the tissue bed.

In addition to its local haemostatic properties, SURGICEL™ Haemostat is bactericidal *in vitro* against a wide range of gram-positive and gram-negative organisms including aerobes and anaerobes. SURGICEL™ Haemostat is bactericidal *in vitro* against strains of species including those of:

<i>Staphylococcus aureus</i>	<i>Serratia marcescens</i>
<i>Staphylococcus epidermidis</i>	<i>Bacillus subtilis</i>
<i>Micrococcus luteus</i>	<i>Proteus vulgaris</i>
<i>Streptococcus pyogenes</i> Group A	<i>Corynebacterium xerosis</i>
<i>Streptococcus</i> Group B	<i>Mycobacterium phlei</i>
<i>Streptococcus salivarius</i>	<i>Clostridium tetani</i>
<i>Branhamella catarrhalis</i>	<i>Clostridium perfringens</i>
<i>Escherichia coli</i>	<i>Bacteroides fragilis</i>
<i>Klebsiella pneumoniae</i>	<i>Enterobacter cloacae</i>
<i>Lactobacillus</i> sp.	<i>Pseudomonas aeruginosa</i>
<i>Salmonella enteritidis</i>	<i>Pseudomonas stutzeri</i>
<i>Shigella dysenteriae</i>	<i>Proteus mirabilis</i>

methicillin-resistant *Staphylococcus aureus* (MRSA)

penicillin-resistant *Streptococcus pneumoniae* (PRSP)

vancomycin-resistant *Enterococcus* (VRE)

methicillin-resistant *Staphylococcus epidermidis* (MRSE)

SURGICEL™ Haemostat is not, however, an alternative to systemically applied therapeutic or prophylactic antimicrobial agents.

## STERILITY

SURGICEL™ Haemostat is sterilized using radiation. Do not resterilize. Do not use if package is opened or damaged.

## STORAGE

SURGICEL™ Haemostat should be stored at or below 30°C. Do not use after expiry date.

## HOW SUPPLIED

SURGICEL™ Haemostat is available in single-use, sterile foil pouches in a variety of sizes.

## TRACEABILITY

A traceability label that identifies the type, size, expiry date, and batch code of the device is included in every package. This label should be affixed or electronically added to the patient's permanent medical record to clearly identify the device that was implanted. The following specific information can be found on the device packaging label: Catalogue number, batch code, expiry and manufacturing date, manufacturer name, address and the website, and a Unique Device Identification bar code with the Global Trade Item Number information.

## INFORMATION TO BE CONVEYED TO THE PATIENT

### Device Description

The SURGICEL™ device is sterile and single use. The device is made from plant cellulose. It can be made into many formats. Formats include a knitted fabric, loose fibers or a powder.

The SURGICEL™ device is used to help control bleeding from different types of blood vessels.

### Device Materials

The SURGICEL™ device is made completely from plant cellulose. It is made to be used safely in surgery. It is processed in a way that makes it absorbable in the body. Actual absorption time depends on where and how much is used.

### Information for Safe Use

There is always a risk of adverse reactions with any type of surgery. Some reactions could be due to your overall health condition. Some reactions could be related to just having the surgery. Some reactions could be related to the materials being placed into the body. It is difficult to know the source of some adverse reactions. Contact your doctor about any signs of infection or adverse reactions. Be sure to keep up with your doctor visits after surgery.

The device will not cause any issues with a medical test called an MRI. MRI means magnetic resonance imaging.

### Expected Device Lifetime and Follow-up

The SURGICEL™ device is absorbed by the body over time. There is almost no tissue reaction. Absorption time depends on several things, including where and how much is used.

### Patient Information Portal

ic.jnjmedicaldevices.com

You can find any updates about the device on this website. Your implant card will contain information about your specific device. This can also be found in your medical record.


Figure 1


Figure 2A


Figure 2B


Figure 3

## SYMBOLS USED ON LABELING

<b>REF</b>	Catalogue number
<b>MD</b>	Medical Device / for implant card: Device Name
	Do not use if package is damaged and consult instructions for use
	Caution
	Do not re-use
	Do not resterilize
	Single sterile barrier system with protective packaging inside
<b>STERILE R</b>	Sterilized using irradiation
<b>LOT</b>	Batch code
	Date of manufacture
	Use-by date
<b>UDI</b>	Unique Device Identifier
	Manufacturer
<b>EC REP</b>	Authorized Representative in the European Community/ European Union
	Packaging unit
	Consult instructions for use or consult electronic instructions for use. EU: Call paper on demand help desk to get paper copies free of charge within 7 days.
<b>MR</b>	MR Safe
	Upper limit of temperature
	Absorbable Hemostat
	Patient information website / for implant card: Information website for patients
	Patient Identification / for Implant card: Patient Name or patient ID
	Date / for implant card: Date of Implantation
	Health care centre or doctor / for implant card: Name and address of the implanting healthcare institution / provider

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Johnson & Johnson Medical GmbH