IOTest Conjugated Antibody CD38-PE

REF A07779 100 tests; 2 mL, 20 µL / test

TABLE OF CONTENTS

| | |
|--|-----|
| English | 2 |
| Français (France) (fr-fr) | 7 |
| Deutsche (de-de) | 12 |
| Italiano (it-it) | 17 |
| Español (es-es) | 22 |
| Português (Portugal) (pt-pt) | 27 |
| Dansk (da-dk) | 32 |
| Svenska (sv-se) | 37 |
| Norsk (nb-no) | 42 |
| Suomi (fi-fi) | 47 |
| Ελληνικά (el-gr) | 52 |
| 日本語 (ja-jp) | 57 |
| 中文 (中国) (zh-cn) | 62 |
| Lietuvių (lt-lt) | 66 |
| Magyar (hu-hu) | 71 |
| Polski (pl-pl) | 76 |
| Čeština (cs-cz) | 81 |
| Slovák (sk-sk) | 86 |
| 한국어 (ko-kr) | 91 |
| Türkçe (tr-tr) | 96 |
| Русский (ru-ru) | 101 |
| Eestlane (et-ee) | 106 |
| Hrvatski (hr-hr) | 111 |
| Български (bg-bg) | 116 |
| 中文 (台灣) (zh-tw) | 121 |
| Română (ro-ro) | 126 |
| Slovenščina (sl-si) | 131 |
| Srpski (sr-sp) | 136 |
| Latviešu (lv-lv) | 141 |
| Українська (uk-ua) | 146 |
| Português (Brasil) (pt-br) | 151 |
| Nederlands (nl-nl) | 156 |
| Tiếng Việt (vi-vn) | 161 |
| Қазақша (kk-kz) | 166 |
| APPENDIX | 171 |
| REFERENCES | 172 |

| | Specifications |
|----------------------|---|
| Specificity | CD38 |
| Clone | LS198-4-3 |
| Hybridoma | SP2/0 x balb/c |
| Immunogen | Human T cell line HUT 78 |
| Immunoglobulin | IgG1 |
| Species | Mouse |
| Purification | Affinity Chromatography |
| Fluorochrome | R Phycoerythrin (PE) |
| Molar ratio | PE / Ig: 0.5-1.5 |
| λ excitation | 488 nm |
| Emission Peak | 575 nm |
| Buffer | PBS pH 7.2 plus 2 mg / mL BSA and 0.1% NaN ₃ |

IOTest Conjugated Antibody CD38-PE

REF A07779 100 tests; 2 mL, 20 μ L / test

For *In Vitro* Diagnostic Use

INTENDED USE

This fluorochrome-conjugated antibody allows the qualitative and non automated identification of cell populations expressing the CD38 antigen present in human biological samples using flow cytometry (see section "Samples" below).

PRINCIPLE

This test is based on the ability of specific monoclonal antibodies to bind to the antigenic determinants expressed by leucocytes.

Specific staining of the leucocytes is performed by incubating the sample with the IOTest reagent. The red cells are then removed by lysis and the leucocytes, which are unaffected by this process, are analyzed by flow cytometry.

The flow cytometer measures light diffusion and the fluorescence of cells. It makes possible the delimitation of the population of interest within the electronic window defined on a histogram, which correlates the orthogonal diffusion of light (Side Scatter or SS) and the diffusion of narrow-angle light (Forward Scatter or FS). Other histograms combining two of the different parameters available on the cytometer can be used as supports in the gating stage depending on the application chosen by the user.

The fluorescence of the delimited cells is analyzed in order to distinguish the positively-stained events from the unstained ones. The results are expressed as a percentage of positive events in relation to all the events acquired by the gating.

INTENDED USER

This product is intended for laboratory professional use.

CLINICAL RELEVANCE

The CD38-PE is a CD38 antibody used to identify and characterize cells expressing the CD38 antigen by flow cytometry. This product alone cannot and is not intended to generate any diagnostic conclusion.

When used in combination with other markers, this product can be used in one or more of the following functions:

- To aid in the differential diagnosis of hematologically abnormal patients suspected of having hematopoietic neoplasm and to monitor patients with known hematopoietic neoplasm.
- To aid in the prognosis of patients having hematopoietic neoplasm.
- To aid in diagnosis of patients with suspected immunodeficiency and monitor patients with known immunodeficiency.

See following references :

1; 2; 3; 4; 5; 6; 7; 8; 9; 10

SAMPLES

Venous blood must be taken using sterile tubes containing an EDTA salt as the anticoagulant.

The samples should be kept at room temperature (18 – 25°C) and not shaken. The samples should be homogenized by gentle agitation prior to taking the test sample.

The samples must be analyzed within 24 hours of venipuncture.

CONCENTRATION

See lot specific Certificate of Analysis at www.beckman.com.

WARNING AND PRECAUTIONS

1. Do not use the reagent beyond the expiry date.

2. Do not freeze.
3. Let it come to room temperature (18 – 25°C) before use.
4. Minimize exposure to light.
5. Avoid microbial contamination of the reagents, or false results may occur.
6. Antibody solutions containing sodium azide (NaN₃) should be handled with care. Do not take internally and avoid all contact with the skin, mucosa and eyes.
 Moreover, in an acid medium, sodium azide can form the potentially dangerous hydrazoic acid. If it needs to be disposed of, it is recommended that the reagent be diluted in a large volume of water before pouring it into the drainage system so as to avoid the accumulation of sodium azide in metal pipes and to prevent the risk of explosion.
7. All blood samples must be considered as potentially infectious and must be handled with care (in particular: the wearing of protective gloves, gowns and goggles).
8. Never pipet by mouth, and avoid all contact of the samples with the skin, mucosa and eyes.
9. Blood tubes and disposable material used for handling should be disposed of in ad hoc containers intended for incineration.
10. Reagents and waste should be eliminated according to local requirements.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at beckman.com/techdocs

STORAGE AND STABILITY

This reagent must be kept at between 2 and 8°C and protected from light, before and after the vial has been opened.

Closed vial shelf life per stability study: 1095 days.

Stability of opened vial: the reagent is stable for 180 days.

EVIDENCE OF DETERIORATION

Any change in the physical appearance of the reagents may indicate deterioration and the reagent should not be used.

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-742-2345 (USA or Canada) or contact your local Beckman Coulter Representative.

CONTENTS

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76). To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT:

- Sampling tubes and material necessary for sampling.
- Automatic pipettes with disposable tips for 20, 100 and 500 µL.
- Plastic haemolysis tubes.
- Red cell lysis reagent with washing stage after lysis. For example: VersaLyse (Ref. A09777).
- Leucocyte fixation reagent. For example : IOTest 3 Fixative Solution (Ref. A07800).
- Isotype control PE : IOTest reagent (Ref. A07796).
- Buffer (PBS: 0.01 M sodium phosphate; 0.145 M sodium chloride; pH 7.2).
- Centrifuge.
- Automatic agitator (Vortex type).
- Flow cytometer.

PROCEDURE WITH VERSALYSE REAGENT

For each sample analyzed, in addition to the test tube, one control tube can be added in which the cells are mixed in the presence of the isotype control (Ref. A07796).

1. Add 20 µL of specific IOTest conjugated antibody to each test tube and if necessary, 20 µL of the isotype control to each control tube.
2. Add 100 µL of the test sample to each tube. Vortex the tubes gently.
3. Incubate for 15 to 20 minutes at room temperature (18 – 25°C), protected from light.
4. Then perform lysis of the red cells. Refer to VersaLyse (Ref. A09777) leaflet and follow preferably the procedure called “with concomitant fixation”, which consists of adding 1 mL of the “Fix-and-Lyse” mixture prepared extemporaneously. Vortex immediately for one second and incubate for 10 minutes at room temperature, protected from light.
5. Centrifuge for 5 minutes at 150 x g at room temperature.
6. Remove the supernatant by aspiration.

7. Resuspend the cell pellet using 3 mL of PBS.
8. Repeat step 5.
9. Remove the supernatant by aspiration and resuspend the cell pellet using:
 - 0.5 mL or 1 mL of PBS plus 0.1% of formaldehyde if the preparations are to be kept less than 24 hours. (A 0.1% formaldehyde PBS can be obtained by diluting 12.5 µL of the IOTest 3 Fixative Solution (Ref. A07800) at its 10X concentration in 1 mL of PBS).
 - 0.5 mL or 1 mL of PBS without formaldehyde, if the preparations are to be analyzed within 2 hours.

Note: In all cases, keep the preparations between 2 and 8°C and protected from light.

EXPECTED VALUES

In our laboratories, the whole blood samples of 25 apparently healthy donors were treated using the reagent described above. The results obtained for the count of the positive events of interest with this reagent are given in the table below :

| Positive Target | Number | Mean (%) | SD | CV (%) |
|---------------------|--------|----------|-------|--------|
| Lymphocytes | 25 | 53.38 | 14.16 | 26.52 |
| Monocytes | 25 | 97.96 | 4.53 | 4.63 |
| Granulocytes | 25 | 99.94 | 0.27 | 0.27 |

These values are intended to be representative only. Each laboratory should establish its own expected values from the local population of normal donors.

PERFORMANCE

Performance data are obtained using the procedure described above on less than 24 hour-old blood samples previously collected in sterile tubes with EDTA salt as anticoagulant. Analysis is performed within 2 hours following immunostaining.

SPECIFICITY

The CD38 antigen, also called the T10 antigen, is a single chain type II transmembranous glycoprotein the terminal N section of which is cytoplasmic (11,12). This glycoprotein with a molecular weight close to 45 kDa is expressed on a large number of hematopoietic cells: its distribution depends on the stage of differentiation and the degree of cellular activation. The functional molecule appears as a dimer or a multimer (13,14,15). One function initially attributed to CD38 is the regulation of activation and proliferation of human T lymphocytes (16,17, and 18).

CD38 ligation by agonistic mAbs induces rapid Ca²⁺ fluxes and triggers the phosphorylation of a cascade of intracellular substrates, leading to activation of the NF-κB complex (19,20). Within the B-cell compartment, CD38 is expressed at high levels only by committed progenitor bone marrow (early BM cells are CD38-) and by B lymphocytes in germinal centers, by terminally differentiated plasma cells and in activated tonsils (21,22,23). Instead, mature virgin and memory B lymphocytes express low levels of the molecule. More precisely, at the level of B cell lymphocytic differentiation, expression of the CD38 molecule appears during the early stages and disappears during B cell maturation and reappears during the final stage of differentiation (the plasmocytic stage).

The CD38 molecule is expressed by activated T and B cell lymphocytes (24,25), monocytes, the majority of NK cells, by medullary thymocytes (26,27) and on platelets (28) and red cells (29). CD38 expression on T cell circulating lymphocytes is weak. The monoclonal antibody LS198 was assigned to CD38 during the 5th HLDA Workshop on Human Leucocyte Differentiation Antigens, Boston, USA, in 1993 (WS Code: T-CD38.06, Section T) (29).

PRECISION

The percent positive values were determined using Whole Blood. Each sample was run 4 times, twice a day for 1 day on 2 instruments using 2 lots of CD38-PE Monoclonal Antibody Reagents. Measurements (% positive) were made on Navios flow cytometer. Analysis was conducted based on the CLSI method EP5-A2: Evaluation of Precision Performance of Quantitative Measurement Methods.

Our acceptance criteria is dependent on the number of positive events measured for each population:

- If positive event <1500, CV <15%
- If positive event >1500, CV <10%

| Lymphocytes | | | | | | | |
|---|----------------|-----------------|-----------|-----------|-----------|-----------|-------|
| Number of positive events (Mean) = 7398 | | | | | | | |
| | Inter-Operator | Inter-cytometer | Intra-lot | Inter-lot | Inter-run | Intra-run | Total |
| CV (%) | 1.27 | 3.61 | 2.39 | 0.57 | 1.84 | 0.83 | 4.37 |
| RESULTS | PASS | PASS | PASS | PASS | PASS | PASS | PASS |

| Monocytes | | | | | | | |
|---|----------------|-----------------|-----------|-----------|-----------|-----------|-------|
| Number of positive events (Mean) = 1050 | | | | | | | |
| | Inter-Operator | Inter-cytometer | Intra-lot | Inter-lot | Inter-run | Intra-run | Total |
| CV (%) | 0.53 | 0.60 | 1.22 | 0.35 | 0.53 | 0.97 | 1.41 |
| RESULTS | PASS | PASS | PASS | PASS | PASS | PASS | PASS |

| Granulocytes | | | | | | | |
|--|----------------|-----------------|-----------|-----------|-----------|-----------|-------|
| Number of positive events (Mean) = 13927 | | | | | | | |
| | Inter-Operator | Inter-cytometer | Intra-lot | Inter-lot | Inter-run | Intra-run | Total |
| CV (%) | 0.03 | 0.04 | 0.10 | 0.02 | 0.03 | 0.09 | 0.11 |
| RESULTS | PASS | PASS | PASS | PASS | PASS | PASS | PASS |

ACCURACY

The accuracy of CD38-PE was assessed by comparing the results with a reference reagent as the predicate on a set of whole blood samples run on a Navios flow cytometer. The bias between test and reference reagent was determined based on the difference between test results. If the bias is within the allowable error range or the p-value indicates no significant difference (> 0.05), then test results for the two reagents are considered as equivalent.

The results obtained are summarized in the following table below:

| Number of donors = 25 | | | | |
|-----------------------|---------------|--------------------------|----------|---------|
| Positive Target | Mean Δ | Δ % Cell criteria | p-value | RESULTS |
| Lymphocytes | 2.71 | <5 | 3.70E-04 | PASS |
| Monocytes | -0.74 | <5 | 0.350 | PASS |
| Granulocytes | -0.05 | <5 | 0.370 | PASS |

LIMIT OF BLANK AND LIMIT OF DETECTION

A study was conducted in accordance with CLSI EP17-A2, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures. The Limit of Detection (LOD) is the lowest concentration of analyte that can be consistently detected. Results obtained are summarized in the following table:

| Positive Target | Limit of Blank (cell/ μ L) | Limit of Detection (cell/ μ L) |
|---------------------|--------------------------------|------------------------------------|
| Lymphocytes | 4 | 7 |
| Monocytes | 1 | 3 |
| Granulocytes | 2 | 4 |

LIMITATIONS

1. Flow cytometry may produce false results if the cytometer has not been aligned perfectly, if fluorescence leaks have not been correctly compensated for and if the regions have not been carefully positioned.
2. It is preferable to use a RBC lysis technique with a washing step, as this reagent has not been optimized for "no wash" lysis techniques.
3. Accurate and reproducible results will be obtained as long as the procedures used are in accordance with the technical insert leaflet and compatible with good laboratory practices.
4. The conjugated antibody of this reagent is calibrated so as to offer the best specific signal/non-specific signal ratio. Therefore, it is important to adhere to the reagent volume/sample volume ratio in every test.
5. In the case of a hyperleucocytosis, dilute the blood in PBS so as to obtain a value of approximately 5×10^9 leucocytes/L (30).
6. In certain disease states, such as severe renal failure or haemoglobinopathies, lysis of red cells may be slow, incomplete or even impossible. In this case, it is recommended to isolate mononucleated cells using a density gradient (Ficoll, for example) prior to staining (31).
7. In patients treated with anti-human monoclonal antibody therapies, detection of the specific targeted antigens may be diminished or absent due to partial or complete blocking by the therapeutic antibody.
8. The CD38-PE results should be interpreted in light of the total clinical presentation of the patient, including: symptoms, clinical history, data from additional tests, and other appropriate information.

See the Appendix for examples and references.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The Summary of Safety and Performance is available from the EUDAMED database: ec.europa.eu/tools/eudamed.

REVISION HISTORY

| REVISION AC: | Release date : October 2019 |
|--|--|
| REVISION | Release Date |
| AW | May 2022 |
| Updates to comply with the Beckman Coulter Global labelling Policy and according to IVD-R (EU)2017/746 requirements: | |
| Added sections | BSI 2797 Number, Intended User, Clinical Relevance, Concentration, Precision, Accuracy, Limit of blank and limit of detection, Additional Information, Revision History. |
| Added information | See sections Limitations |
| Phrasing or typographic updates | See sections Procedure, Performance, Limitations, Warning and Precautions, Storage and Stability. |

| | |
|------------------|--|
| Removed sections | Example of clinical applications, Reagents, Intra-laboratory reproducibility, Linearity |
| Updated sections | Intended Use, GHS Hazard Classification, Evidence of deterioration, Procedure, Appendix. |

| | |
|------------------|-----------------------|
| REVISION | Release Date |
| AX | |
| Updated sections | Add Kazakh |
| Updated sections | Storage and stability |

| | |
|------------------------|-----------------------------|
| Revised Version | AY |
| Revised Date | May 2024 |
| Updated Translations | Clinical relevance : |
| Updated sections | Revision History |

Symbols Key

Glossary of Symbols is available at beckman.com/techdocs (document number B60062)

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
APPENDIX

- EXAMPLES

Whole blood samples stained with IOTest CD38-PE Conjugated Antibody.

Acquisition is performed with a Beckman Coulter Navios flow cytometer equipped with the Navios analysis software.

REFERENCES

1. C. E. Gooden, P. Jones, R. Bates, W. M. Shallenberger, U. Surti, S. H. Swerdlow, C. G. Roth, "CD49d shows superior performance characteristics for flow cytometric prognostic testing in chronic lymphocytic leukemia/small lymphocytic lymphoma." *Cytometry B Clin Cytom.* 94, 129–135 (2018).
2. P. M. J. Theunissen, A. van den Branden, A. Van Der Sluijs-Gelling, V. De Haas, A. Beishuizen, J. J. M. van Dongen, V. H. J. Van Der Velden, "Understanding the reconstitution of the B-cell compartment in bone marrow and blood after treatment for B-cell precursor acute lymphoblastic leukaemia." *Br. J. Haematol.* 178, 267–278 (2017).
3. A. Al-Mawali, D. Gillis, I. Lewis, "Immunoprofiling of leukemic stem cells CD34 /CD38-/CD123 delineate FLT3/ITD-positive clones." *Journal of Hematology & Oncology.* 9 (2016), doi:10.1186/s13045-016-0292-z.
4. I. S. Nijhof, R. W. J. Groen, H. M. Lokhorst, B. van Kessel, "Upregulation of CD38 expression on multiple myeloma cells by all-trans retinoic acid improves the efficacy of daratumumab." *Leukemia.* 29, 2039–2049 (2015).
5. P. P. Nesmiyanov, B. E. Tolkachev, A. V. Strygin, "ZO-1 expression shows prognostic value in chronic B cell leukemia." *Immunobiology.* 221, 6–11 (2016).
6. A. A. Basabaeen, E. A. Abdelgader, O. S. BaHashwan, E. A. Babekir, N. M. Abdelateif, S. A. Bamusa, S. O. Abdelrahim, O. A. Altayeb, E. A. Fadul, I. K. Ibrahim, "Combined analysis of ZAP-70 and CD38 expression in sudanese patients with B-cell chronic lymphocytic leukemia." *BMC Res Notes.* 12, 282 (2019).
7. E. A. Abdelgader, N. H. Eltayeb, T. A. Eltahir, O. A. Altayeb, "Evaluation of CD38 expression in Sudanese patients with chronic lymphocytic leukemia." *BMC Research Notes.* 11 (2018), doi:10.1186/s13104-018-3926-0.
8. G. d'Ettorre, G. Ceccarelli, S. Serafino, N. Giustini, "Dominant enrichment of phenotypically activated CD38 HLA-DR CD8 T cells, rather than CD38 HLA-DR CD4 T cells, in HIV/HCV coinfecting patients on antiretroviral therapy." *Journal of Medical Virology.* 88, 1347–1356 (2016).
9. M. Spaan, S. J. Hulleger, B. J. B. Beudeker, K. Kreeft, G. W. van Oord, Z. M. A. Groothuisink, M. van Tilborg, B. Rijnders, R. J. de Knecht, M. A. A. Claassen, A. Boonstra, "Frequencies of Circulating MAIT Cells Are Diminished in Chronic HCV, HIV and HCV/HIV Co-Infection and Do Not Recover during Therapy." *PLoS ONE.* 11, e0159243 (2016).
10. B. H. Davis, J.T. Holden, M. C. Bene, et al. "2006 Bethesda International Consensus Recommendations on the Flow Cytometric Immunophenotypic Analysis of Hematolymphoid Neoplasia: Medical Indications." *Cytometry Part B (Clinical Cytometry)* 72B:S5-S13 (2007).
11. Mehta, K., Shahid, U., Malavasi, F., "Human CD38, a cell-surface protein with multiple functions", 1996, *The FASEB Journal*, 10, 1408-1417.
12. Alessio M, et al. (1990) CD38 molecule: structural and biochemical analysis on human T lymphocytes, thymocytes, and plasma cells. *J. Immunol.* 145:878-84.
13. Hara-Yokoyama M, et al. (2012) Tetrameric interaction of the ectoenzyme CD38 on the cell surface enables its catalytic and raft-association activities. *Structure.* 20:1585-95.
14. Mallone R, et al. (1998) Characterization of a CD38-like 78-kilodalton soluble protein released from B cell lines derived from patients with Xlinked agammaglobulinemia. *J. Clin. Invest.* 101:2821-30.
15. Zhao YJ, Lam CM, Lee HC. (2012) The membrane-bound enzyme CD38 exists in two opposing orientations. *Sci. Signal.* 5:ra67.
16. Funaro A, et al. (1990) Involvement of the multilineage CD38 molecule in a unique pathway of cell activation and proliferation. *J. Immunol.* 145:2390-6.
17. Quarona V1, Zaccarello G, Chillemi A, et al. CD38 and CD157: a long journey from activation markers to multifunctional molecules. *Cytometry B Clin Cytom.* 2013 Jul-Aug;84(4):207-17.
18. Moodley K1, Coetzee LM, Glencross DK. Decentralised CD38 activation monitoring: aspects of practical implementation and standardisation. *J Immunol Methods.* 2012 Apr 30;378(1-2):121-7.
19. Fedele G, et al. (2004) CD38 is expressed on human mature monocyte-derived dendritic cells and is functionally involved in CD83 expression and IL-12 induction. 34:1342-50.
20. Buggins AG, et al. (2010) Interaction with vascular endothelium enhances survival in primary chronic lymphocytic leukemia cells via NFkappaB activation and de novo gene transcription. *Cancer Res.* 70:7523-33.
21. Damle RN, et al. (1999) Ig V gene mutation status and CD38 expression as novel prognostic indicators in chronic lymphocytic leukemia. *Blood.* 94:1840-7.
22. Chiorazzi N. (2012) Implications of new prognostic markers in chronic lymphocytic leukemia. *Hematology Am. Soc. Hematol. Educ. Program.* 2012:76-87.
23. Malavasi F, et al. (2011) CD38 and chronic lymphocytic leukemia: a decade later. *Blood.* 118:3470-8.
24. Terhorst C, et al. (1981) Biochemical studies of the human thymocyte cell-surface antigens T6, T9 and T10. *Cell.* 23:771-80.
25. Fedele G, et al. (2004) CD38 is expressed on human mature monocyte-derived dendritic cells and is functionally involved in CD83 expression and IL-12 induction. *Eur. J. Immunol.* 34:1342-50.
26. Funaro, A., Horenstein, A.L., Malavasi, F., "Human CD38: a versatile leukocyte molecule with emerging clinical perspectives", 1995, *Fund. Clin. Immunol.*, 3, 3, 101-113.
27. Ramaschi, G., Torti, M., Festetics, E.T., Sinigaglia, F., Malavasi, F., Balduini, C., "Expression of cyclic ADP-Ribose-synthetizing CD38 molecule on human platelet membrane", 1996, *Blood*, 6, 87, 2308- 2313.
28. Zocchi, E., Franco, L., Guida, L., Benatti, U., Bargellesi, A., Malavasi, F., Lee, H.C., DeFlora, A., " A single protein immunologically identified as CD38 display NAD+ Glycohydrolase, ADP-Ribosyl Cyclase and cyclic ADP-Ribose Hydrolase activities at the outer surface erythrocytes", 1993, *Biochem. Biophys. Res. Com.*, 3, 196, 1459-1465.
29. Boumsell, L., "T-cell antigens: section report", 1995, *Leucocyte Typing V, White Cell Differentiation Antigens*, Schlossman, S.F., et al., Eds., Oxford Univ. Press, 241-279.

30. H42-A2 Vol.27 No. 16. P30. Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; Approved Guideline - Second Edition.
31. Ibrahim FF, Ghannam MM, Ali FM., "Effect of dialysis on erythrocyte membrane of chronically hemodialyzed patients"., 2002. Ren Fail., Nov;24(6):779-90.

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Тел. +7 (495) 228 67 00, e-mail: beckman.ru@beckman.com

IMMUNOTECH SAS A Beckman Coulter Company, 130, Avenue de Lattre de Tassigny,

BP 177, 13276 Marseille Cedex 9, France, +(33) 4 91 17 27 27

www.beckman.com