IOTest Conjugated Antibody CD56-PC7

REF A21692 100 tests; 1 mL, 10 µL / test

TABLE OF CONTENTS

English	2
Français (France) (fr-fr)	6
Deutsche (de-de)	10
Italiano (it-it)	14
Español (es-es)	18
Português (Portugal) (pt-pt)	22
Dansk (da-dk)	26
Svenska (sv-se)	30
Norsk (nb-no)	34
Suomi (fi-fi)	38
Ελληνικά (el-gr)	42
日本語 (ja-jp)	46
中文 (中国) (zh-cn)	50
Lietuvių (lt-lt)	54
Magyar (hu-hu)	58
Polski (pl-pl)	62
Čeština (cs-cz)	66
Slovák (sk-sk)	70
한국어 (ko-kr)	74
Türkçe (tr-tr)	78
Русский (ru-ru)	82
Eestlane (et-ee)	86
Hrvatski (hr-hr)	90
Български (bg-bg)	94
中文 (台灣) (zh-tw)	98
Română (ro-ro)	102
Slovenščina (sl-si)	106
Srpski (sr-sp)	110
Latviešu (lv-lv)	114
Українська (uk-ua)	118
Português (Brasil) (pt-br)	122
Nederlands (nl-nl)	126
Tiếng Việt (vi-vn)	130
Қазақша (kk-kz)	134
APPENDIX	138
REFERENCES	139

	Specifications
Specificity	CD56
Clone	N901 (NKH-1)
Hybridoma	NS1 x balb/c
Immunogen	Human chronic myeloid leukemia (hCML) cells
Immunoglobulin	IgG1
Species	Mouse
Purification	Affinity Chromatography
Fluorochrome	R Phycoerythrin-Cyanine 7 (PC7)
Molar ratio	PC7 / Ig: 0.5 - 1.5
λ excitation	488 nm
Emission Peak	770 nm
Buffer	PBS pH 7.2 plus 2 mg / mL BSA and 0.1% NaN ₃

IOTest Conjugated Antibody CD56-PC7

REF A21692 100 tests; 1 mL, 10 μ L / test

For *In Vitro* Diagnostic Use

INTENDED USE

This fluorochrome-conjugated antibody allows the qualitative and non automated identification of cell populations expressing the CD56 antigen present in human biological samples using flow cytometry (see section "Samples" below).

PRINCIPLE

This test is based on the ability of specific monoclonal antibodies to bind to the antigenic determinants expressed by leucocytes.

Specific staining of the leucocytes is performed by incubating the sample with the IOTest reagent. The red cells are then removed by lysis and the leucocytes, which are unaffected by this process, are analyzed by flow cytometry.

The flow cytometer measures light diffusion and the fluorescence of cells. It makes possible the delimitation of the population of interest within the electronic window defined on a histogram, which correlates the orthogonal diffusion of light (Side Scatter or SS) and the diffusion of narrow-angle light (Forward Scatter or FS). Other histograms combining two of the different parameters available on the cytometer can be used as supports in the gating stage depending on the application chosen by the user.

The fluorescence of the delimited cells is analyzed in order to distinguish the positively-stained events from the unstained ones. The results are expressed as a percentage of positive events in relation to all the events acquired by the gating.

INTENDED USER

This product is intended for laboratory professional use.

CLINICAL RELEVANCE

The CD56-PC7 is a CD56 antibody used to identify and characterize cells expressing the CD56 antigen by flow cytometry. This product alone cannot and is not intended to generate any diagnostic conclusion.

When used in combination with other markers, this product can be used in one or more of the following functions:

- To aid in the differential diagnosis of hematologically abnormal patients suspected of having hematopoietic neoplasm and to monitor patients with known hematopoietic neoplasm.
- To aid in diagnosis of patients with suspected immunodeficiency and monitor patients with known immunodeficiency.
- To aid in the prognosis of patients having hematopoietic neoplasm.
- To monitor transplantation process or results.

See following references :

1; 2; 3; 4; 5; 6; 7; 8; 9

SAMPLES

Venous blood must be taken using sterile tubes containing an EDTA salt as the anticoagulant.

The samples should be kept at room temperature (18 – 25°C) and not shaken. The samples should be homogenized by gentle agitation prior to taking the test sample.

The samples must be analyzed within 24 hours of venipuncture.

CONCENTRATION

See lot specific Certificate of Analysis at www.beckman.com.

WARNING AND PRECAUTIONS

1. Do not use the reagent beyond the expiry date.
2. Do not freeze.
3. Let it come to room temperature (18 – 25°C) before use.

4. Minimize exposure to light.
5. Avoid microbial contamination of the reagents, or false results may occur.
6. Antibody solutions containing sodium azide (NaN_3) should be handled with care. Do not take internally and avoid all contact with the skin, mucosa and eyes.
 Moreover, in an acid medium, sodium azide can form the potentially dangerous hydrazoic acid. If it needs to be disposed of, it is recommended that the reagent be diluted in a large volume of water before pouring it into the drainage system so as to avoid the accumulation of sodium azide in metal pipes and to prevent the risk of explosion.
7. All blood samples must be considered as potentially infectious and must be handled with care (in particular: the wearing of protective gloves, gowns and goggles).
8. Never pipet by mouth, and avoid all contact of the samples with the skin, mucosa and eyes.
9. Blood tubes and disposable material used for handling should be disposed of in ad hoc containers intended for incineration.
10. Reagents and waste should be eliminated according to local requirements.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at beckman.com/techdocs

STORAGE AND STABILITY

This reagent must be kept at between 2 and 8°C and protected from light, before and after the vial has been opened.

Closed vial shelf life per stability study: 365 days.

Stability of opened vial: the reagent is stable for 180 days.

EVIDENCE OF DETERIORATION

Any change in the physical appearance of the reagents may indicate deterioration and the reagent should not be used.

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-742-2345 (USA or Canada) or contact your local Beckman Coulter Representative.

CONTENTS

Sodium azide preservative may form explosive compounds in metal drain lines. See NIOSH Bulletin: Explosive Azide Hazard (8/16/76).

To avoid the possible build-up of azide compounds, flush wastepipes with water after the disposal of undiluted reagent. Sodium azide disposal must be in accordance with appropriate local regulations.

MATERIALS NEEDED BUT NOT SUPPLIED WITH KIT:

- Sampling tubes and material necessary for sampling.
- Automatic pipettes with disposable tips for 10, 100 and 500 µL.
- Plastic haemolysis tubes.
- Red cell lysis reagent with washing stage after lysis. For example: VersaLyse (Ref. A09777).
- Leucocyte fixation reagent. For example : IOTest 3 Fixative Solution (Ref. A07800).
- Isotype control PC7 : IOTest reagent (Ref. 737662).
- Buffer (PBS: 0.01 M sodium phosphate; 0.145 M sodium chloride; pH 7.2).
- Centrifuge.
- Automatic agitator (Vortex type).
- Flow cytometer.

PROCEDURE WITH VERSALYSE REAGENT

For each sample analyzed, in addition to the test tube, one control tube can be added in which the cells are mixed in the presence of the isotype control (Ref. 737662).

1. Add 10 µL of specific IOTest conjugated antibody to each test tube and if necessary, 10 µL of the isotype control to each control tube.
2. Add 100 µL of the test sample to each tube. Vortex the tubes gently.
3. Incubate for 15 to 20 minutes at room temperature (18 – 25°C), protected from light.
4. Then perform lysis of the red cells. Refer to VersaLyse (Ref. A09777) leaflet and follow preferably the procedure called “with concomitant fixation”, which consists of adding 1 mL of the “Fix-and-Lyse” mixture prepared extemporaneously. Vortex immediately for one second and incubate for 10 minutes at room temperature, protected from light.
5. Centrifuge for 5 minutes at 150 x g at room temperature.
6. Remove the supernatant by aspiration.
7. Resuspend the cell pellet using 3 mL of PBS.
8. Repeat step 5.
9. Remove the supernatant by aspiration and resuspend the cell pellet using:
 - 0.5 mL or 1 mL of PBS plus 0.1% of formaldehyde if the preparations are to be kept less than 24 hours. (A 0.1% formaldehyde PBS can be obtained by diluting 12.5 µL of the IOTest 3 Fixative Solution (Ref. A07800) at its 10X concentration in 1 mL of PBS).

- 0.5 mL or 1 mL of PBS without formaldehyde, if the preparations are to be analyzed within 2 hours.

Note: In all cases, keep the preparations between 2 and 8°C and protected from light.

EXPECTED VALUES

In our laboratories, the whole blood samples of 25 apparently healthy donors were treated using the reagent described above. The results obtained for the count of the positive events of interest with this reagent are given in the table below :

Positive Target	Number	Mean (%)	SD	CV (%)
Lymphocytes	25	19.10	8.53	44.66

These values are intended to be representative only. Each laboratory should establish its own expected values from the local population of normal donors.

PERFORMANCE

Performance data are obtained using the procedure described above on less than 24 hour-old blood samples previously collected in sterile tubes with EDTA salt as anticoagulant. Analysis is performed within 2 hours following immunostaining.

SPECIFICITY

The molecular weight of the CD56 heavily glycosylated protein (NKH1 molecule) is 200 – 220 kDa (10,11). It is expressed on a subpopulation of peripheral blood lymphocytes that demonstrate non-major histocompatibility complex restricted cytotoxicity (10,12). The N901 (NKH-1) monoclonal antibody reacts with the majority of NK cells (10,11). It also reacts with a minor subpopulation of CD3 positive T cells that mediates reduced cytotoxic activity (12). This antibody does not react with monocytes, granulocytes, erythrocytes or B lymphocytes. The N901 (NKH-1) antibody has been assigned to the CD56 cluster of differentiation during the 4th International Workshop on Human Leucocyte Differentiation Antigens in Vienna, Austria, in 1989 (13).

PRECISION

The percent positive values were determined using Whole Blood. Each sample was run 4 times, twice a day for 1 day on 2 instruments using 2 lots of CD56-PC7 Monoclonal Antibody Reagents. Measurements (% positive) were made on Navios flow cytometer. Analysis was conducted based on the CLSI method EP5-A2: Evaluation of Precision Performance of Quantitative Measurement Methods.

Our acceptance criteria is dependent on the number of positive events measured for each population:

- If positive event <1500, CV <15%
- If positive event >1500, CV <10%

Lymphocytes							
Number of positive events (Mean) = 1090							
	Inter-Operator	Inter-cytometer	Intra-lot	Inter-lot	Inter-run	Intra-run	Total
CV (%)	2.45	0.83	4.35	0.83	2.73	2.34	4.51
RESULTS	PASS	PASS	PASS	PASS	PASS	PASS	PASS

ACCURACY

The accuracy of CD56-PC7 was assessed by comparing the results with a reference reagent as the predicate on a set of whole blood samples run on a Navios flow cytometer. The bias between test and reference reagent was determined based on the difference between test results. If the bias is within the allowable error range or the p-value indicates no significant difference (> 0.05), then test results for the two reagents are considered as equivalent.

The results obtained are summarized in the following table below:

Number of donors = 25				
Positive Target	Mean Δ	Δ % Cell criteria	p-value	RESULTS
Lymphocytes	-0.95	<5	0.000	PASS

LIMIT OF BLANK AND LIMIT OF DETECTION

A study was conducted in accordance with CLSI EP17-A2, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures. The Limit of Detection (LOD) is the lowest concentration of analyte that can be consistently detected. Results obtained are summarized in the following table:

Positive Target	Limit of Blank (cell/μL)	Limit of Detection (cell/μL)
Lymphocytes	6	10

LIMITATIONS

1. Flow cytometry may produce false results if the cytometer has not been aligned perfectly, if fluorescence leaks have not been correctly compensated for and if the regions have not been carefully positioned.
2. It is preferable to use a RBC lysis technique with a washing step, as this reagent has not been optimized for "no wash" lysis techniques.
3. Accurate and reproducible results will be obtained as long as the procedures used are in accordance with the technical insert leaflet and compatible with good laboratory practices.
4. The conjugated antibody of this reagent is calibrated so as to offer the best specific signal/non-specific signal ratio. Therefore, it is important to adhere to the reagent volume/sample volume ratio in every test.
5. In the case of a hyperleucocytosis, dilute the blood in PBS so as to obtain a value of approximately 5×10^9 leucocytes/L (14).
6. In certain disease states, such as severe renal failure or haemoglobinopathies, lysis of red cells may be slow, incomplete or even impossible. In this case, it is recommended to isolate mononucleated cells using a density gradient (Ficoll, for example) prior to staining (15).
7. Due to the tandem structure of the fluorochrome, PC7 also emits light at 575 nm. This secondary emission peak varies from lot-to-lot of PC7. Therefore, for multi-color analysis, the compensation matrix should be carefully checked when changing the lot of a -PC7 conjugate.
8. In patients treated with anti-human monoclonal antibody therapies, detection of the specific targeted antigens may be diminished or absent due to partial or complete blocking by the therapeutic antibody.

9. The CD56-PC7 results should be interpreted in light of the total clinical presentation of the patient, including: symptoms, clinical history, data from additional tests, and other appropriate information.

See the Appendix for examples and references.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

ADDITIONAL INFORMATION

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/746/EU on In vitro Diagnostic Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

The Summary of Safety and Performance is available from the EUDAMED database: ec.europa.eu/tools/eudamed.

REVISION HISTORY

REVISION AD:	Release date : September 2019
REVISION	Release Date
AW	March 2022
Updates to comply with the Beckman Coulter Global labelling Policy and according to IVD-R (EU)2017/746 requirements:	
Added sections	BSI 2797 Number, Intended User, Clinical Relevance, Concentration, Precision, Accuracy, Limit of blank and limit of detection, Additional Information, Revision History.
Added information	See sections Limitations
Phrasing or typographic updates	See sections Procedure, Performance, Limitations, Warning and Precautions, Storage and Stability.
Removed sections	Example of clinical applications, Reagents, Intra-laboratory reproducibility, Linearity
Updated sections	Intended Use, GHS Hazard Classification, Evidence of deterioration, Procedure, Appendix.
REVISION	Release Date
AX	
Updated sections	Add Kazakh
Updated sections	Storage and stability
Revised Version	AY
Revised Date	May 2024
Updated Translations	Clinical relevance :
Updated sections	Revision History

Symbols Key

Glossary of Symbols is available at beckman.com/techdocs (document number B60062)

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
APPENDIX
- EXAMPLES

Whole blood samples stained with IOTest CD56-PC7 Conjugated Antibody.
Acquisition is performed with a Beckman Coulter Navios flow cytometer equipped with the Navios analysis software.

REFERENCES

1. R. Suzuki, S. Ohtake, J. Takeuchi, M. Nagai, Y. Kodera, M. Hamaguchi, S. Miyawaki, T. Karasuno, S. Shimodaira, R. Ohno, S. Nakamura, T. Naoe, "The clinical characteristics of CD7+ CD56+ acute myeloid leukemias other than M0.", *Int. J. Hematol.* 91, 303–309 (2010).
2. S. Bhardwaj, F. Ahmad, H. Wedemeyer, M. Cornberg, J. Schulze zur Wiesch, J. van Lunzen, S. K. Sarin, R. E. Schmidt, D. Meyer-Olson, "Increased CD56bright NK cells in HIV-HCV co-infection and HCV mono-infection are associated with distinctive alterations of their phenotype.", *Virology* 13, 67 (2016).
3. P. Mantegani, G. Tambussi, L. Galli, C. T. Din, A. Lazzarin, C. Fortis, "Perturbation of the natural killer cell compartment during primary human immunodeficiency virus 1 infection primarily involving the CD56 bright subset.", *Immunology* 129, 220–233 (2010).
4. F. Angelot-Delettre, S. Biichle, C. Ferrand, E. Seilles, B. Gaugler, V. Harivel, M. A. Rosenthal-Allieri, E. Deconinck, P. Saas, F. Garnache-Ottou, "Intracytoplasmic detection of TCL1—but not ILT7-by flow cytometry is useful for blastic plasmacytoid dendritic cell leukemia diagnosis.", *Cytometry A* 81, 718–724 (2012).
5. W. Cui, D. Zhang, M. T. Cunningham, L. Tilzer, "Leukemia-associated aberrant immunophenotype in patients with acute myeloid leukemia: changes at refractory disease or first relapse and clinicopathological findings.", *Int J Lab Hematol.* 36, 636–649 (2014).
6. L. Shang, X. Chen, Y. Liu, X. Cai, "The immunophenotypic characteristics and flow cytometric scoring system of acute myeloid leukemia with t(8;21) (q22;q22); RUNX1-RUNX1T1.", *International Journal of Laboratory Hematology*, 41, 23–31 (2019).
7. Theodorakos I, Paterakis G, Papadakis V, Vicha A, Topakas G, Jencova P, Karchilaki E, Taparkou A, Tsagarakis NJ, Polychronopoulou S, "Interference of bone marrow CD56(+) mesenchymal stromal cells in minimal residual disease investigation of neuroblastoma and other CD45(-) /CD56(+) pediatric malignancies using flow cytometry.", *Pediatr Blood Cancer*. 2019 Aug;66(8):e27799. doi: 10.1002/pbc.27799. Epub 2019 May 7.
8. Hong HS, Eberhard JM, Keudel P, Bollmann BA, Ahmad F, Ballmaier M, Bhatnagar N, Zielinska-Skowronek M, Schmidt RE, Meyer-Olson D, "Phenotypically and functionally distinct subsets contribute to the expansion of CD56-/CD16+ natural killer cells in HIV infection.", *AIDS*, 2010 Jul 31;24(12):1823-34. doi: 10.1097/QAD.0b013e32833b556f.
9. B. H. Davis, J.T. Holden, M. C. Bene, et al. "2006 Bethesda International Consensus Recommendations on the Flow Cytometric Immunophenotypic Analysis of Hematolymphoid Neoplasia: Medical Indications." *Cytometry Part B (Clinical Cytometry)* 72B:S5-S13 (2007).
10. Griffin, J.D., Hercend, T., Beveridge, R., Schlossman, S.F., "Characterization of an antigen expressed by human natural killer cells", 1983, *J. Immunol.*, 130, 2947-2951.
11. Hercend, T., Griffin, J.D., Bensussan, A., Schmidt, R.E., Edson, M.A., Brennan, A., Murray, C., Daley, J.F., Schlossman, S.F., Ritz, J., "Generation of monoclonal antibodies to a human natural killer clone: characterization of two natural killer-associated antigens, NKH1A and NKH2, expressed on subsets of large granular lymphocytes", 1985, *J. Clin. Invest.*, 75, 932-943.
12. Lanier, L.L., Le, A.M., Civin, C.I., Loken, M.R., Phillips, J.H., "The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes", 1986, *J. Immunol.*, 136, 4480-4486.
13. Schubert, J., Lanier, L.L., Schmidt, R.E., "Cluster report: CD56", 1989, *Leucocyte Typing IV, White Cell Differentiation Antigens*. W. Knapp, et al., Eds., Oxford University Press, 699-702.
14. H42-A2 Vol.27 No. 16. P30. Enumeration of Immunologically Defined Cell Populations by Flow Cytometry; Approved Guideline - Second Edition.
15. Ibrahim FF, Ghannam MM, Ali FM., "Effect of dialysis on erythrocyte membrane of chronically hemodialyzed patients", 2002. *Ren Fail.*, Nov;24(6):779-90.

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Тел. +7 (495) 228 67 00, e-mail: beckman.ru@beckman.com

IMMUNOTECH SAS A Beckman Coulter Company, 130, Avenue de Lattre de Tassigny,

BP 177, 13276 Marseille Cedex 9, France, +(33) 4 91 17 27 27

www.beckman.com