**BECKMAN  
COULTER**

## Table of Contents

English — English .....	2
Brazilian Portuguese — Português do Brasil .....	3
Bulgarian — Български .....	4
Chinese Simplified — 简体中文 .....	5
Chinese Traditional — 繁體中文 .....	6
Croatian — Hrvatski .....	7
Czech — Čeština .....	8
Danish — Dansk .....	9
Dutch — Nederlands .....	10
Estonian — Eesti .....	11
Finnish — Suomi .....	12
French — Français .....	13
German — Deutsch .....	14
Greek — Ελληνικά .....	15
Hungarian — Magyar .....	16
Italian — Italiano .....	17
Japanese — 日本語 .....	18
Korean — 한국어 .....	19
Latvian — Latviski .....	20
Lithuanian — Lietuvių k. ....	21
Norwegian — Norsk .....	22
Polish — Polski .....	23
Portuguese — Português .....	24
Romanian — Română .....	25
Russian — Русский .....	26
Serbian — Srpski .....	27
Slovak — Slovensky .....	28
Slovenian — Slovenščina .....	29
Spanish — Español .....	30
Swedish — Svenska .....	31
Turkish — Türkçe .....	32
Ukrainian — Українська .....	33
Vietnamese — Tiếng Việt .....	34

## Declaration of Conformity

Immunotech SAS hereby ensures and declares that the product(s) listed below comply with the requirement of the *In-Vitro* Diagnostic Medical Devices Regulation 2017/746.

This EU Declaration of Conformity is issued under the sole responsibility of the manufacturer.

<p><b>Product(s):</b> CD117-PC7, ref. B49221</p> <p><b>EMDN Code:</b> W010308</p> <p><b>BUDI- DI:</b> 150995903B49221SG</p> <p><b>Risk Class:</b> Class C, Rule 3 h) (Article 47 in accordance with Annex VIII)</p> <p><b>Intended Purpose:</b> Intended use: The CD117-PC7 antibody allows the qualitative and non automated identification of cell populations expressing the CD117 antigen present in human biological samples using flow cytometry (see section "Samples" below).</p> <p>Intended user: This product is intended for laboratory professional use.</p> <p>Clinical relevance: The CD117-PC7 is a CD117 antibody used to identify and characterize cells expressing CD117 antigen by flow cytometry. This product alone cannot and is not intended to generate any diagnostic conclusion. When used in combination with other markers, this product can be used in one or more of the following functions:</p> <ul style="list-style-type: none"> <li>To aid in the differential diagnosis of hematologically abnormal patients suspected of having hematopoietic neoplasm and to monitor patients with known hematopoietic neoplasm.</li> <li>To aid in the prognosis of patients having hematopoietic neoplasm.</li> </ul> <p>Sample: Venous blood (as described in IFU)</p> <p><b>Common Specification(s)</b> None</p>	<p><b>Conformity Assessment Procedure</b></p> <p>Conformity Assessment is based on a Quality Management System and on Assessment of Technical Documentation (Article 48 in accordance with Annex IX)</p>								
<p>Signed for and on behalf of Immunotech SAS, the Legal Manufacturer</p>  <p>Date <u>8 Nov 2021</u></p> <p>Legal manufacturer BU RA delegate <b>Name:</b> Sophie ROQUES-VIOUIN <b>Title:</b> PRRC: Person Responsible for Regulatory Compliance Senior Manager Quality Assurance &amp; Regulatory Affairs <b>Place of Issue:</b> Marseille, France</p>	<p><b>Notified Body</b> BSI Group The Netherlands B.V. Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, Netherlands Phone: +31 (0)20 346 07 80 Email: info.nl@bsigroup.com</p> <p>Notified Body number: 2797</p>								
 <p>Immunotech SAS A Beckman Coulter Company 130 Avenue de Lattre de Tassigny B.P. 177 - 13276 Marseille, Cedex 9 France +(33) 4 91 17 27 27 <b>Manufacturer SRN:</b> FR-MF-000011121 <b>IVDR Certificate Number:</b> IVDR-738451</p>	<p><b>Document Control</b></p> <table border="1"> <tr> <td><b>Issue Date:</b></td> <td>03 Nov 2021</td> </tr> <tr> <td><b>Revision Level:</b></td> <td>1</td> </tr> <tr> <td><b>Starting Lot:</b></td> <td>200501</td> </tr> <tr> <td><b>DoC Filename:</b></td> <td>B49221-TF-810</td> </tr> </table>	<b>Issue Date:</b>	03 Nov 2021	<b>Revision Level:</b>	1	<b>Starting Lot:</b>	200501	<b>DoC Filename:</b>	B49221-TF-810
<b>Issue Date:</b>	03 Nov 2021								
<b>Revision Level:</b>	1								
<b>Starting Lot:</b>	200501								
<b>DoC Filename:</b>	B49221-TF-810								

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Bu AB Uygunluk Beyanı yalnızca üreticinin sorumluluğunda verilmektedir.