NADAL® Rota-Adenovirus Test (test cassette)

REF 481015

de	Gebrauchsanweisung	2	pt	Instruções de Utilização	22
en	Instructions for use	6	cs	Návod k použití	26
fr	Instructions d'utilisation	10	fi	Käyttöohje	34
es	Instrucciones de uso	14	no	Bruksanvisning	38
it	Istruzioni per l'uso	18		Symbols	43
pl	Sposób użycia	20		Our Teams	44

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Moers
Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Regensburg
Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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Staphylococcus aureus

14. Referenzen

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelm I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

Rev. 2, 2019-05-23 OM

1. Intended Use

The NADAL® Rota-Adenovirus Test is a rapid visual immune-assay for the qualitative detection of rotavirus and adenovirus in human faecal specimens. This test is intended to be used as an aid in the diagnosis of rotavirus and adenovirus infections. The test is designed for professional *in-vitro* diagnostic use only.

2. Introduction and Clinical Significances

Rotavirus is the most common pathogen responsible for acute gastroenteritis, mainly in young children. Its discovery in 1973 and association with infantile gastroenteritis represented a very important advancement in the study of acute non-bacterial gastroenteritis. Rotavirus is transmitted by the oro-faecal route with an incubation period of 1-3 days. Although specimens collected between the second and the fifth day of illness are ideal for antigen detection, rotavirus may still be detected while diarrhoea continues. Rotaviral gastroenteritis may result in mortality for populations at risk such as infants, the elderly and immunocompromised patients. In moderate climates, rotaviral infections occur mainly in the winter months. Endemics and epidemics affecting thousands of people have been reported. Up to 50% of the analysed specimens from hospitalised children suffering from acute enteric disease were positive for rotavirus. The viruses replicate in the cell nucleus and tend to produce a characteristic cytopathic effect (CPE). Because rotavirus is extremely difficult to culture, it is unusual to isolate the virus for the diagnosis of infection. Instead, a variety of techniques have been developed to detect rotavirus in feces.

Acute diarrheal disease in young children is a major cause of morbidity worldwide and is a leading cause of mortality in developing countries. Research has shown that enteric adenoviruses, such as Ad40 and Ad41, are also one of the primary causes of diarrhoea in children, second only to rotaviruses. These viral pathogens have been identified throughout the world and can cause diarrhoea in children all year round. Infections are most frequently seen in children less than two years of age, but have also been found in patients of all ages. Further studies indicate that adenoviruses are associated with 4-15% of all hospitalized cases of viral gastroenteritis.

Rapid and accurate diagnosis of gastroenteritis caused by adenovirus is helpful in establishing the etiology of gastroenteritis and related patient management. Other diagnostic techniques such as electron microscopy (EM) and nucleic acid hybridization are expensive and labor-intensive. Due to the self-limited nature of adenoviral infection, such expensive and labor-intensive tests may not be necessary.

3. Test Principle

The NADAL® Rota-Adenovirus Test is designed to detect rotavirus and adenovirus through visual interpretation of color development in the internal strip.

In this assay Rotavirus and Adenovirus are detected with the aid of specific antibodies. After the addition of sample (feces diluted in buffer) the color-labelled antibodies will specifically bind to the respective virus if it is present in the sample. When these antibody-virus complexes migrate along the membrane by capillary action, they are captured by another antibody

specific against Rotavirus or Adenovirus in the test line regions for the respective virus. If rotavirus is present in the sample a red test line is developed in the "R" marked region. If adenovirus is present in the sample a red test line will develop in the "A" marked region. If both viruses are present (mixed infection) two test lines will develop. If no virus is present, the color-labelled antibody will not bind to the virus specific antibodies in the test line regions. No red test lines are formed. So the presence of a red test line indicates a positive result, while its absence indicates a negative result.

To serve as a procedural control, a red line will always develop in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 individually pouched NADAL® Rota-Adenovirus test cassettes
- 10 specimen collection tubes with specimen diluent buffer for sample collection and dilution
- 10 disposable pipettes for sample collection of extremely liquid specimen
- 1 package insert

5. Additional Materials Required

Patient:

- Means to collect a stool sample.
It should be ensured that samples have no contact with the water in the toilet bowl in order to avoid dilution or contamination with e. g. detergents. Special stool collection units can be provided by nal von minden GmbH on request.

Clinic or lab:

- Tissue paper for breaking the tip of the specimen collection tube
- Timer

6. Storage & Stability

Store the test as packaged either at room temperature or refrigerated (2-30°C).

Under these conditions the test cassette and specimen diluent buffer are stable until the expiration date. The test cassette should remain in the sealed pouch containing a desiccant until use.

Do not freeze.

Do not use beyond the expiration date.

7. Warnings and Precautions

- For *in-vitro* diagnostic use by professionals only.
- For single use. Do not reuse tests.
- Do not interchange or mix reagents from different lots.
- Do not use the test if its foil pouch is damaged.
- Do not use the test after the expiration date.
- This test contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled by observing usual safety precautions (e.g. do not ingest or inhale).

- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not eat, drink or smoke in the area where the specimens or kits are handled. Protective clothing such as laboratory coats, disposable gloves and eye protection are recommended. Observe established precautions against microbiological hazards throughout testing. Follow standard procedures for proper disposal of specimens in accordance with local regulations.
- The specimen diluent buffer contains a small amount of sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of extraction solution or extracted samples, always flush with copious quantities of water to prevent azide buildup.
- Humidity and temperature can adversely affect results.
- The components of the test (e.g. antibodies/chemicals) do not cause any danger if the test is used according to the instructions.
- Follow the instructions for use carefully. Inform the patients of the procedures of the collection and dilution of the stool sample.

8. Specimen Collection and Preparation

Conditions for optimal sample collection

The NADAL® Rota-Adenovirus Test is only intended for use with human faecal specimen that has been diluted in the provided buffer.


Viral detection is facilitated by collecting specimens at the onset of diarrheic symptoms. It has been reported that the maximum excretion of rotavirus in the feces of patients with gastroenteritis occurs 2-5 days after the onset of symptoms. For adenovirus maximum excretion is about 3-13 days after onset of symptoms. If the specimens are collected long after the onset of diarrheic symptoms, the quantity of antigen may not be sufficient to obtain a positive result or the antigens detected may not be linked to the diarrheic episode.

Stool samples should be stored at 2-8°C immediately after collection and processed within 48 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided.

Specimen collection and preparation by the patient

For stool sample collection, the patient is given one of the specimen collection tubes of the test set and a disposable pipette. Stool sample should be collected in the following way:

1. Any clean and dry container or water repellent paper can be used for sample collection. Ensure that the stool sample has no direct contact with the water of the toilet bowl to avoid dilution or contamination with detergents. An amount of 1-2 mL or 1-2 g stool is sufficient.
2. Transfer a small amount of the stool into the specimen collection tube:

For solid specimens:

Unscrew the cap of the specimen collection tube, then randomly stab the specimen collection applicator in at least 3 different sites of the faecal specimen to collect

approximately 50 mg of feces (equivalent to 1/4 of a pea).

For liquid specimens:

If the stool is too liquid to stick to the applicator, the provided disposable pipette can be used. Hold the pipette vertically, aspirate some stool, and then transfer 2 drops (approximately 50 µL) into the specimen collection tube containing the diluent buffer.

3. Place the applicator back into the tube and screw the cap tightly.
4. Shake the specimen collection tube to mix the specimen and the dilution buffer. Be careful not to break the tip of the collection tube.
5. Wrap the specimen collection tube in a plastic bag. Store it in a cool place and return it to the clinic within 24 hours.


Note:

If the patient feels uncomfortable to dilute the stool specimen in the collection tube himself, he may return a container with the unprocessed stool sample to the clinic. The transfer of the specimen into the buffer of the specimen collection tube can then be carried out as described above by the personnel of the clinic or lab.

9. Test Procedure

1. Allow the test device and the diluted stool sample to reach room temperature (15-30°C) prior to testing.
2. Remove the test cassette from its pouch when ready to perform the test. The test cassette must be at room temperature to prevent condensation of moisture on the membrane. Label the test cassette on the provided space with patient or control identification.
3. Shake the collection tube thoroughly to ensure proper mixing of the faecal sample with the diluent buffer.
4. Using a piece of tissue paper, break the tip of the collection tube with a twisting motion.
5. Hold the collection tube vertically and dispense 2-3 drops of solution into the round sample well (S) of the test cassette by applying a gentle pressure to the walls of the tube. Avoid air bubbles in the sample well or splashes of liquid into the rectangular result window.
6. Start the timer. As the test begins to work, a reddish colored liquid front will be seen moving across the membrane.

Wait for the colored line(s) to appear. The result should be read after 10 minutes. Strong positive results may be observed sooner. Do not interpret the result after 20 minutes.

10. Result Interpretation

For reading of the test results, the colored lines that develop in the test result window are interpreted.

Positive result for Rotavirus, Adenovirus or both viruses

Rotavirus Positive:

A colored line develops in the control line region (C) and another colored line develops in test line region for Rotavirus (R).

Adenovirus Positive:

A colored line develops in the control line region (C) and another colored line develops in test line region for Adenovirus (A).

Rotavirus and Adenovirus Positive:

A colored line develops in the control line region (C) and two other colored lines develop in the test line regions for Rotavirus (R) and Adenovirus (A), respectively.

Note:

The intensity of the color in the test line regions (R/A) may vary depending on the concentration of the target antigens present in the specimen. Therefore, any shade of color in the test region should be considered positive. Besides, the amount of virus can not be determined by this qualitative test.

Negative result:

One colored line develops in the control line region (C). No line develops in the test line regions for Rotavirus and Adenovirus.

Invalid result:

The control line (C) fails to develop. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the test set immediately and contact your distributor.

Insufficient specimen volume, insufficient specimen migration, or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the presence of visible particles inhibited the specimen migration, they should be removed by centrifugation or sedimentation. Transfer a part of the sample into a tube, sediment the particles by a brief centrifugation and pipette approximately 80-120 µL of the supernatant into the sample well of a new test cassette. Alternatively allow the particles to settle in the upright specimen collection tube and use 80-120 µL from the topside of the liquid. If the problem persists, discontinue using the test immediately and contact your local distributor.

Note: When faecal samples are tested, the background may appear slightly yellowish due to the color of the faecal samples. This is acceptable as long as it does not interfere with

the interpretation of the test result. The test is invalid if the background fails to clear and obscures the reading of the result.

11. Quality Control

An internal procedural control is included in the test cassette. A colored line appearing in the control line region (C) is an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Limitations

- The NADAL® Rota-Adenovirus Test is for professional *in-vitro* diagnostic use, and should be used for qualitative detection of rotavirus and adenovirus only.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of rotavirus or adenovirus infection.
- In order to prevent pseudo-outbreaks, positive test results in newborns must be confirmed using an alternative testing method (PCR).
- If visible blood is present in faecal specimens, false-positive results cannot be ruled out and must therefore be verified using an alternative test method.

13. Performance Characteristics

The performance of the NADAL® Rota-Adenovirus Test (Adenovirus) has been evaluated with 210 clinical specimens collected from children and young adults in comparison with ELISA. The results are summarised in the following table:

Table: NADAL® Rota-Adenovirus Test (Adenovirus) vs. ELISA

		NADAL® Rota-Adenovirus Test (Adenovirus)		
		+	-	Total
ELISA	+	82	1	83
	-	0	127	127
	Total	82	128	210

Relative sensitivity: >98.8%

Relative specificity: >99.9%

Overall agreement: >99.5%

The performance of the NADAL® Rota-Adenovirus Test (Rotavirus) has been evaluated with 242 clinical specimens collected from children and young adults in comparison with ELISA. The results are summarised in the following table:

Table: NADAL® Rota-Adenovirus Test (Rotavirus) vs. ELISA

		NADAL® Rota-Adenovirus Test (Rotavirus)		
		+	-	Total
ELISA	+	79	3	82
	-	0	160	160
	Total	79	163	242

Relative sensitivity: >96.3%

Relative specificity: >99.9%

Overall agreement: >98.8%

Cross-reactivity

Cross-reactivity with the organisms listed below has been studied at 1.0×10^9 organisms/mL. These samples were found negative when tested with the NADAL® Rota-Adenovirus Test.

<i>Acinetobacter calcoaceticus</i>	<i>Group C Streptococcus</i>
<i>Acinetobacter spp</i>	<i>Hemophilus influenzae</i>
<i>Branhamella catarrhalis</i>	<i>Klebsiella pneumoniae</i>
<i>Candida albicans</i>	<i>Neisseria gonorrhoea</i>
<i>Chlamydia trachomatis</i>	<i>Neisseria meningitidis</i>
<i>E. coli</i>	<i>Proteus mirabilis</i>
<i>Enterococcus faecalis</i>	<i>Proteus vulgaris</i>
<i>Enterococcus faecium</i>	<i>Pseudomonas aeruginosa</i>
<i>Gardnerella vaginalis</i>	<i>Salmonella choleraesuis</i>
<i>Group B Streptococcus</i>	<i>Staphylococcus aureus</i>

14. References

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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14. Bibliografia

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.

2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.

3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.

4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.

5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.

6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954

7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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14. Reference

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003; vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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14. Lähteet

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April, 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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