
## Capillary High-Flux Dialysers

FX classix		50	60	80	100
Art.-No.		F00002385	F00002386	F00002387	F00002388
$Q_B$ [mL/min]		200 300	200 300 400	200 300 400	300 400 500
Clearances [mL/min] $Q_D = 500$ mL/min, UF = 0	Cytochrome C	53 55	69 74 76	82 89 92	100 105 108
	Inulin	68 72	88 95 99	103 113 119	122 129 133
	Vitamin B <sub>12</sub>	121 137	139 162 175	155 185 202	201 222 236
	Phosphate	167 204	178 225 252	187 244 279	253 291 316
	Creatinine	178 224	186 243 277	193 259 300	264 309 337
	Urea	191 253	195 266 312	198 279 334	280 336 370
K <sub>0</sub> A Urea	mL/min	866	1068	1394	1429
UF-coefficient (at $Q_B$ max)	mL/h/mmHg	27	38	53	68
S (sieving coefficient)	Albumin Myoglobin $\beta_2$ -Microglobulin	< 0.001 0.1 0.7			
Max. TMP	mmHg	600			
V (blood priming volume)	mL	53	74	95	116
$\Delta P$ (pressure drop blood, $Q_B = 300$ mL/min)	mmHg	121	87	67	55
Max. dialysate flow	mL/min	800	1000	1000	1000
Recommended blood flow range	mL/min	100 – 300	150 – 400	200 – 500	250 – 600
Membrane		Helixone®			
A (effective surface area)	m <sup>2</sup>	1.0	1.4	1.8	2.2
Sterilisation method		INLINE steam			

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
# FX classix Capillary Dialyser

### GENERAL NOTES

Refer to product or carton label for:

	Single use only		Blood pathway sterile. Steam sterilised
	Expiry date		Refer to instructions for use
	Batch		Order number
	Date of manufacture		Storage temperature range
	Units		Use only on machines with exact UF-control

**Indications:** FX classix dialysers are designed for single use in chronic haemodialysis.

**Contraindications:** Special contraindications are unknown. Generally, contraindications for haemodialysis are applicable.

**Side-effects:** In rare cases hypersensitivity reactions may occur during haemodialysis treatment. In severe cases dialysis must be discontinued and the appropriate medication initiated.

The dialyser is steam sterilised and thus contains no sterilisation residues. For patients with a known ethylene oxide hypersensitivity **all** components of the extracorporeal circuit should be ethylene oxide-free.

**Anticoagulation:** It is recommended to introduce an anticoagulant to the extracorporeal circuit. Nature, amount and method of application of an anticoagulant must be prescribed by the responsible physician (e.g. initial heparin bolus of 2000 IU followed by a continuous dose of 1000 IU/h up until the last hour of treatment). Coagulation should be monitored by a standard clotting time test.

**Materials:** Membrane: Helixone® (Polysulfone - PVP blend), Housing: Polypropylene, Potting material: Polyurethane, Sealing-ring: Silicone, Sterility caps: Polypropylene.

Further information may be obtained on request.

### WARNINGS

Due to the high hydraulic permeability of the membrane, the dialysers must only be used with dialysis machines which enable precise volumetric fluid control. Refer to dialysis machine manufacturer's technical and safety instructions for the use of high permeability membranes.

Use only if unit package is intact, sealing caps are in place and the dialyser is undamaged.

Dialysers must not be used after expiry date (see label).

Each dialyser is checked for integrity prior to leaving the factory. If a blood leak should arise, the dialyser must be exchanged.

The dialyser is intended for single-use only. Re-use may be hazardous to both the patient and operator. Cleansing solutions and disinfectants may damage materials employed for the housing, potting and membrane. Safety of use can no longer be guaranteed and the manufacturer assumes no liability.

### HAEMODIALYSIS

#### Priming

Affix dialyser in a vertical position. Aseptically connect the arterial and venous blood lines to the lower and upper dialyser ports respectively. Attach the dialysate tubing connectors with the inlet to the upper and the outlet to the lower dialysate port. Fill dialyser according to the instructions provided with the machine (**there is no need to turn the dialyser**). We recommend bag-free priming with ONLINE<sup>plus</sup> (priming volume: 500 mL). Alternatively the blood compartment can be filled and de-aerated in recirculation mode with isotonic saline solution (e.g. 500 mL bag). Ensure complete de-aeration of the blood compartment.

#### Patient Connection

Connect arterial line to the patient's circulation. Allow blood to flow in the blood line and dialyser (pump speed approx. 100 mL/min). Connect venous blood line to the patient.

Adjust treatment parameters to desired setting.

The weight loss should be monitored.

#### Termination

Perform reinfusion of the blood at the end of the treatment according to the instructions provided with the machine. We recommend the bag-free ONLINE<sup>plus</sup> procedure. Alternatively isotonic saline solution (e.g. 500 mL bag) can be used for reinfusion. The blood should be reinfused completely.

### WARRANTY

Products with manufacturing defects will be replaced if the defect is reported with details of the lot number.

The manufacturer will not be liable for any misuse, improper handling, non compliance with instructions for use and cautionary notes and for any damage incurred subsequent to the manufacturer's delivery of the dialyser.

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