TECHLAB®

C.DIFF QUIK CHEK COMPLETE®

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Staphylococcus aureus

14. Referenzen

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelm I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413

Rev. 2, 2019-05-23 OM

1. Intended Use

The NADAL® Rota-Adenovirus Test is a rapid visual immune-assay for the qualitative detection of rotavirus and adenovirus in human faecal specimens. This test is intended to be used as an aid in the diagnosis of rotavirus and adenovirus infections. The test is designed for professional *in-vitro* diagnostic use only.

2. Introduction and Clinical Significance

Rotavirus is the most common pathogen responsible for acute gastroenteritis, mainly in young children. Its discovery in 1973 and association with infantile gastroenteritis represented a very important advancement in the study of acute non-bacterial gastroenteritis. Rotavirus is transmitted by the oro-faecal route with an incubation period of 1-3 days. Although specimens collected between the second and the fifth day of illness are ideal for antigen detection, rotavirus may still be detected while diarrhoea continues. Rotaviral gastroenteritis may result in mortality for populations at risk such as infants, the elderly and immunocompromised patients. In moderate climates, rotaviral infections occur mainly in the winter months. Endemics and epidemics affecting thousands of people have been reported. Up to 50% of the analysed specimens from hospitalised children suffering from acute enteric disease were positive for rotavirus. The viruses replicate in the cell nucleus and tend to produce a characteristic cytopathic effect (CPE). Because rotavirus is extremely difficult to culture, it is unusual to isolate the virus for the diagnosis of infection. Instead, a variety of techniques have been developed to detect rotavirus in feces.

Acute diarrheal disease in young children is a major cause of morbidity worldwide and is a leading cause of mortality in developing countries. Research has shown that enteric adenoviruses, such as Ad40 and Ad41, are also one of the primary causes of diarrhoea in children, second only to rotaviruses. These viral pathogens have been identified throughout the world and can cause diarrhoea in children all year round. Infections are most frequently seen in children less than two years of age, but have also been found in patients of all ages. Further studies indicate that adenoviruses are associated with 4-15% of all hospitalized cases of viral gastroenteritis.

Rapid and accurate diagnosis of gastroenteritis caused by adenovirus is helpful in establishing the etiology of gastroenteritis and related patient management. Other diagnostic techniques such as electron microscopy (EM) and nucleic acid hybridization are expensive and labor-intensive. Due to the self-limited nature of adenoviral infection, such expensive and labor-intensive tests may not be necessary.

3. Test Principle

The NADAL® Rota-Adenovirus Test is designed to detect rotavirus and adenovirus through visual interpretation of color development in the internal strip.

In this assay Rotavirus and Adenovirus are detected with the aid of specific antibodies. After the addition of sample (feces diluted in buffer) the color-labelled antibodies will specifically bind to the respective virus if it is present in the sample. When these antibody-virus complexes migrate along the membrane by capillary action, they are captured by another antibody

specific against Rotavirus or Adenovirus in the test line regions for the respective virus. If rotavirus is present in the sample a red test line is developed in the "R" marked region. If adenovirus is present in the sample a red test line will develop in the "A" marked region. If both viruses are present (mixed infection) two test lines will develop. If no virus is present, the color-labelled antibody will not bind to the virus specific antibodies in the test line regions. No red test lines are formed. So the presence of a red test line indicates a positive result, while its absence indicates a negative result.

To serve as a procedural control, a red line will always develop in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 individually pouched NADAL® Rota-Adenovirus Test cassettes
- 20 specimen collection tubes with specimen diluent buffer for sample collection and dilution
- 20 disposable pipettes for sample collection of extremely liquid specimen
- 1 package insert

5. Additional Materials Required

Patient:

Means to collect a stool sample.

It should be ensured that samples have no contact with the water in the toilet bowl in order to avoid dilution or contamination with e. g. detergents. Special stool collection units can be provided by nal von minden GmbH on request.

Clinic or lab:

- Tissue paper for breaking the tip of the specimen collection tube
- Timer

6. Storage & Stability

Store the test as packaged either at room temperature or refrigerated (2-30°C).

Under these conditions the test cassette and specimen diluent buffer are stable until the expiration date. The test cassette should remain in the sealed pouch containing a desiccant until use.

Do not freeze.

Do not use beyond the expiration date.

7. Warnings and Precautions

- For *in-vitro* diagnostic use by professionals only.
- For single use. Do not reuse tests.
- Do not interchange or mix reagents from different lots.
- Do not use the test if its foil pouch is damaged.
- Do not use the test after the expiration date.
- This test contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled by observing usual safety precautions (e.g. do not ingest or inhale).

- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not eat, drink or smoke in the area where the specimens or kits are handled. Protective clothing such as laboratory coats, disposable gloves and eye protection are recommended. Observe established precautions against microbiological hazards throughout testing. Follow standard procedures for proper disposal of specimens in accordance with local regulations.
- The specimen diluent buffer contains a small amount of sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of extraction solution or extracted samples, always flush with copious quantities of water to prevent azide buildup.
- Humidity and temperature can adversely affect results.
- The components of the test (e.g. antibodies/chemicals) do not cause any danger if the test is used according to the instructions.
- Follow the instructions for use carefully. Inform the patients of the procedures of the collection and dilution of the stool sample.

8. Specimen Collection and Preparation

Conditions for optimal sample collection

The NADAL® Rota-Adenovirus Test is only intended for use with human faecal specimen that has been diluted in the provided buffer.

Viral detection is facilitated by collecting specimens at the onset of diarrheic symptoms. It has been reported that the maximum excretion of rotavirus in the feces of patients with gastroenteritis occurs 2-5 days after the onset of symptoms. For adenovirus maximum excretion is about 3-13 days after onset of symptoms. If the specimens are collected long after the onset of diarrheic symptoms, the quantity of antigen may not be sufficient to obtain a positive result or the antigens detected may not be linked to the diarrheic episode.

Stool samples should be stored at 2-8°C immediately after collection and processed within 48 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided.

Specimen collection and preparation by the patient

For stool sample collection, the patient is given one of the specimen collection tubes of the test set and a disposable pipette. Stool sample should be collected in the following way:

1. Any clean and dry container or water repellent paper can be used for sample collection. Ensure that the stool sample has no direct contact with the water of the toilet bowl to avoid dilution or contamination with detergents. An amount of 1-2 ml or 1-2 g stool is sufficient.
2. Transfer a small amount of the stool into the specimen collection tube:

For solid specimens:

Unscrew the cap of the specimen collection tube, then randomly insert the specimen collection applicator into at least 3 different sites of the faecal specimen to collect

approximately 50 mg of feces (equivalent to 1/4 of a pea).

For liquid specimens:

If the stool is too liquid to stick to the applicator, the provided disposable pipette can be used. Hold the pipette vertically, aspirate some stool, and then transfer 2 drops (approximately 50 µl) into the specimen collection tube containing the diluent buffer.

3. Place the applicator back into the tube and screw the cap tightly.
4. Shake the specimen collection tube to mix the specimen and the dilution buffer. Be careful not to break the tip of the collection tube.
5. Wrap the specimen collection tube in a plastic bag. Store it in a cool place and return it to the clinic within 24 hours.


Note:

If the patient feels uncomfortable to dilute the stool specimen in the collection tube himself, he may return a container with the unprocessed stool sample to the clinic. The transfer of the specimen into the buffer of the specimen collection tube can then be carried out as described above by the personnel of the clinic or lab.

9. Test Procedure

1. Allow the test device and the diluted stool sample to reach room temperature (15-30°C) prior to testing.
2. Remove the test cassette from its pouch when ready to perform the test. The test cassette must be at room temperature to prevent condensation of moisture on the membrane. Label the test cassette on the provided space with patient or control identification.
3. Shake the collection tube thoroughly to ensure proper mixing of the faecal sample with the diluent buffer.
4. Using a piece of tissue paper, break the tip of the collection tube with a twisting motion.
5. Hold the collection tube vertically and dispense 2-3 drops of solution into the round sample well (S) of the test cassette by applying a gentle pressure to the walls of the tube. Avoid air bubbles in the sample well or splashes of liquid into the rectangular result window.
6. Start the timer. As the test begins to work, a reddish colored liquid front will be seen moving across the membrane.

Wait for the colored line(s) to appear. The result should be read after 10 minutes. Strong positive results may be observed sooner. Do not interpret the result after 20 minutes.

10. Result Interpretation

For reading of the test results, the colored lines that develop in the test result window are interpreted.

Positive result for Rotavirus, Adenovirus or both viruses

Rotavirus Positive:

A colored line develops in the control line region (C) and another colored line develops in test line region for Rotavirus (R).

Adenovirus Positive:

A colored line develops in the control line region (C) and another colored line develops in test line region for Adenovirus (A).

Rotavirus and Adenovirus Positive:

A colored line develops in the control line region (C) and two other colored lines develop in the test line regions for Rotavirus (R) and Adenovirus (A), respectively.

Note:

The intensity of the color in the test line regions (R/A) may vary depending on the concentration of the target antigens present in the specimen. Therefore, any shade of color in the test region should be considered positive. Besides, the amount of virus can not be determined by this qualitative test.

Negative result:

One colored line develops in the control line region (C). No line develops in the test line regions for Rotavirus and Adenovirus.

Invalid result:

The control line (C) fails to develop. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the test set immediately and contact your distributor.

Insufficient specimen volume, insufficient specimen migration, or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the presence of visible particles inhibited the specimen migration, they should be removed by centrifugation or sedimentation. Transfer a part of the sample into a tube, sediment the particles by a brief centrifugation and pipette approximately 80-120 µl of the supernatant into the sample well of a new test cassette. Alternatively allow the particles to settle in the upright specimen collection tube and use 80-120 µl from the topside of the liquid. If the problem persists, discontinue using the test immediately and contact your local distributor.

Note: When faecal samples are tested, the background may appear slightly yellowish due to the color of the faecal samples. This is acceptable as long as it does not interfere with

the interpretation of the test result. The test is invalid if the background fails to clear and obscures the reading of the result.

11. Quality Control

An internal procedural control is included in the test cassette. A colored line appearing in the control line region (C) is an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Limitations

- The NADAL® Rota-Adenovirus Test is for professional *in-vitro* diagnostic use, and should be used for qualitative detection of rotavirus and adenovirus only.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of rotavirus or adenovirus infection.
- In order to prevent pseudo-outbreaks, positive test results in newborns must be confirmed using an alternative testing method (PCR).
- If visible blood is present in faecal specimens, false-positive results cannot be ruled out and must therefore be verified using an alternative test method.

13. Performance Characteristics

The performance of the NADAL® Rota-Adenovirus Test (Adenovirus) has been evaluated with 210 clinical specimens collected from children and young adults in comparison with ELISA. The results are summarised in the following table:

Table: NADAL® Rota-Adenovirus Test (Adenovirus) vs. ELISA

		NADAL® Rota-Adenovirus Test (Adenovirus)		
		+	-	Total
ELISA	+	82	1	83
	-	0	127	127
	Total	82	128	210

Relative sensitivity: >98.8%

Relative specificity: >99.9%

Overall agreement: >99.5%

The performance of the NADAL® Rota-Adenovirus Test (Rotavirus) has been evaluated with 242 clinical specimens collected from children and young adults in comparison with ELISA. The results are summarised in the following table:

Table: NADAL® Rota-Adenovirus Test (Rotavirus) vs. ELISA

		NADAL® Rota-Adenovirus Test (Rotavirus)		
		+	-	Total
ELISA	+	79	3	82
	-	0	160	160
	Total	79	163	242

Relative sensitivity: >96.3%

Relative specificity: >99.9%

Overall agreement: >98.8%

Cross-reactivity

Cross-reactivity with the organisms listed below has been studied at 1.0×10^9 organisms/ml. These samples were found negative when tested with the NADAL® Rota-Adenovirus Test.

<i>Acinetobacter calcoaceticus</i>	Group C Streptococcus
<i>Acinetobacter spp</i>	<i>Hemophilus influenzae</i>
<i>Branhamella catarrhalis</i>	<i>Klebsiella pneumoniae</i>
<i>Candida albicans</i>	<i>Neisseria gonorrhoea</i>
<i>Chlamydia trachomatis</i>	<i>Neisseria meningitidis</i>
<i>E. coli</i>	<i>Proteus mirabilis</i>
<i>Enterococcus faecalis</i>	<i>Proteus vulgaris</i>
<i>Enterococcus faecium</i>	<i>Pseudomonas aeruginosa</i>
<i>Gardnerella vaginalis</i>	<i>Salmonella choleraesuis</i>
Group B Streptococcus	<i>Staphylococcus aureus</i>

14. References

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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14. Bibliografia

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.

2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003; vol.9:247-262.

3. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.

4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.

5. Thomas, Eva E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.

6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954

7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411–1413

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14. Bibliografia

1. Wadell, G. Laboratory Diagnosis of Infectious Diseases: Principles and Practices. New York: Springer-Verlag, Volume II, 1988: 284-300.
2. Wilhelmi I, Roman E, Sanchez-Fauquier A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262.
3. Cubitt, WD (1982) Rotavirus infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38.
4. Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." Journal of Clinical Microbiology, June 1989; 27(6): 1155-1158.
5. Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." Am. J. Clin. Pathol. 1994; 101:742-746.
6. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954
7. Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. American Journal of Infection Control 2018; 46(12):1411-1413

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- Wood, D. J., K. Bijlsma, J. C. de Jong, and C. Tonkin. "Evaluation of a Commercial Monoclonal Antibody-Based Enzyme Immunoassay for Detection of Adenovirus Types 40 and 41 in Stool Specimens." *Journal of Clinical Microbiology*, June 1989; 27(6): 1155-1158.
- Thomas, Eva. E., D. Roscoe, L. Book, B. Bone, L. Browne, and V. Mah. "The Utility of Latex Agglutination Assays in the Diagnosis of Pediatric Viral Gastroenteritis." *Am. J. Clin. Pathol.* 1994; 101:742-746.
- Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. *Journal of Microbiology, Immunology and Infection* 2016; 49:947e954
- Jones N K et al. Adenovirus pseudo-outbreak in a large UK neonatal intensive care unit. *American Journal of Infection Control* 2018; 46(12):1411-1413

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15. BIBLIOGRAFIA

1. VIEIRA, A. et al., "Inflammatory bowel disease activity assessed by fecal calprotectin and lactoferrin: correlation with laboratory parameters, clinical, endoscopic and histological indexes", BMC Research Notes 2009, 2:221.
2. HANAI, H. et al. "Relationship Between Fecal Calprotectin, Intestinal Inflammation, and Peripheral Blood Neutrophils in Patients with Active Ulcerative Colitis" Digestive Diseases and Sciences, Sept. 2004, Vol 49, No 9, pp 1438-1443.
3. BONNIN TOMAS, A. et al. "Calprotectina fecal como marcador diferencia entre patología gastrointestinal orgánica y funcional". Rev. Esp. de Enf. Dig. 2007, Vol 99, No 12, pp. 689-693.

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15. Referenzen

- Centers for Disease Control (CDC). Chlamydia trachomatis infections. Policy guidelines for prevention and control. MMWR Morb Mortal Wkly Rep. 1985 Aug 23; 34 Suppl 3: 53S-74S.
- Tichonova L, Borisenko K, Ward H, Meheus A, Gromyko A, Renton A. Epidemics of syphilis in the Russian Federation: trends, origins, and priorities for control. Lancet. 1997 Jul 19; 350(9072): 210-3.
- Norgard MV, Chamberlain NR, Swancutt MA, Goldberg MS. Cloning and expression of the major 47-kilodalton surface immunogen of *Treponema pallidum* in *Escherichia coli*. Infect Immun. 1986 Nov; 54(2): 500-6.

Rev. 1, 2023-02-23 OM

1. Intended Use

The NADAL® Syphilis Test is a lateral flow chromatographic immunoassay for the qualitative detection of anti-*Treponema pallidum* (*T. pallidum*) IgM and IgG in human whole blood, serum or plasma specimens. The test is intended for use as an aid in the diagnosis of syphilis in patients suspected to have a syphilis infection (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® Syphilis Test is designed for professional use only.

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
Treponema pallidum (*T. pallidum*), a spirochete bacterium with an outer cell membrane and a cytoplasmic membrane, is the causative agent of the venereal disease known as syphilis. Although syphilis rates are declining in the United States after an epidemic between 1986 and 1990, the incidence of syphilis in Europe has been increasing since 1992, especially in the states of the Russian Federation, where peaks of 263 cases per 100,000 have been reported. In addition, the positive rate of serological test results for syphilis in HIV-infected individuals has recently been on the rise.

The serological detection of specific antibodies to *T. pallidum* has been long recognised in the diagnosis of syphilis since the natural course of the infection is characterised by periods without clinical manifestations. The antibody response to *T. pallidum* can be detected within 4 to 7 days after the syphilis chancre appears, allowing the early detection and the diagnosis of a syphilis infection.

A variety of antigens, such as Cardiolipin (RPR test), VDRL antigen and *T. pallidum* extracts derived from *in-vitro* culture or inoculated rabbit testes, have been used in syphilis serological tests. However, RPR and VDRL antigens are not treponemal specific and whole *T. pallidum* extracts are not reproducible and contain a certain amount of contaminating materials such as flagella, which may lead to a non-specific reaction in assays of test serum.

3. Test Principle

The NADAL® Syphilis Test enables the detection of anti-*T. pallidum* IgM and IgG through the visual interpretation of colour development on the test strip. Specific recombinant *T. pallidum* antigens are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with the *T. pallidum* antigens which are conjugated to coloured particles and precoat onto the conjugate pad of the test strip. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there are a sufficient number of anti-*T. pallidum* antibodies in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result. The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 40 NADAL® Syphilis test strips, incl. disposable pipettes (40 µL)
- 2 buffer bottles (4 mL each)*
- 1 package insert

*containing the following preservative: ProClin™ 300: <0.03%.

No hazard labelling for ProClin™ 300 is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below the exemption threshold of <0.03%.

5. Additional Materials Required

- Specimen collection containers (suitable for specimen material to be tested)
- Centrifuge (for serum or plasma specimens only)
- Alcohol pads
- Lancets (for fingerstick whole blood specimens only)
- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test strips are stable until the expiry date printed on the foil pouches. Test strips must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).

- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

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The NADAL® Syphilis Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.

Fingerstick whole blood should be tested immediately.

Venipuncture whole blood specimens

Containers containing anticoagulants, such as EDTA, sodium citrate or sodium heparin should be used for the preparation of venous whole blood or plasma specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time.

If the test is to be run within 2 days of specimen collection, whole blood collected by venipuncture should be stored at 2-8°C.

Do not freeze whole blood specimens.

Serum and plasma specimens

Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.


Icteric, lipemic, haemolysed, viscous, heat-treated and contaminated specimens may lead to inaccurate test results.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.


1. Remove the test strip from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test strip with the patient or control identification.
2. Place the test strip on a clean and level surface.
3. **a) For serum, plasma or venous whole blood specimens**

- Holding a pipette vertically, transfer 2 drops (approximately 80 µL) of the specimen to the sample pad of the test strip.
- Holding the buffer bottle vertically, transfer 1 drop (approximately 40 µL) of buffer to the sample pad.

b) For fingerstick whole blood specimens using hanging drops:

- Position the patient's finger so that a drop of blood is exactly above the sample pad of the test strip.
- Allow 2 hanging drops of fingerstick whole blood to fall into the centre of the sample pad of the test strip.
- Holding the buffer bottle vertically, transfer 1 drop (approximately 40 µL) of buffer to the sample pad. **Avoid squeezing the finger, as this might lead to inaccurate test results.**

4. Start the timer.

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

5. Wait for the coloured line(s) to appear. Read the test result after 10 minutes. Do not interpret the result after more than 30 minutes.

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region (C). No line develops in the test line region (T).

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

The internal procedural control is included in the test strip:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Syphilis Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of anti-*T. pallidum* IgM and IgG in human whole blood, serum or plasma specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of *T. pallidum* antibodies can be determined using this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- Some specimens containing unusually high titers of heterophile antibodies or rheumatoid factor may affect test results.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.

- The NADAL® Syphilis Test only detects the presence of *T. pallidum* antibodies in specimens and should not be used as the sole criterion for a diagnosis of *T. pallidum* infection.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a *T. pallidum* infection.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® Syphilis Test was evaluated using clinical whole blood specimens in comparison to another commercially available TPHA Syphilis rapid test.

The results are presented in the following table:

NADAL® Syphilis Test	Another TPHA Syphilis rapid test			
		Positive	Negative	Total
	Positive	411	8	419
	Negative	6	529	535
	Total	417	537	954

Diagnostic sensitivity: 98.6% (96.9% - 99.5%)*

Diagnostic specificity: 98.5% (97.1% - 99.4%)*

Overall agreement: 98.5% (97.6% - 99.2%)*

*95% confidence interval

Analytical performance

Analytical specificity

Cross-reactivity study

Negative whole blood specimens were spiked with the following, potentially cross-reacting substances and tested using the NADAL® Syphilis Test:

Substance	Concentration
HAMA	1000 ng/mL
HBsAb	Positive
HbCAb	Positive
anti- <i>H. pylori</i> IgG	Positive
anti-rubella virus IgM	Positive
rheumatoid factors (RF)	8400 IU/mL
anti-HBeAg IgM	Positive
anti-HCV IgM	Positive
anti-EBV Ab	Positive
anti- <i>Toxoplasma gondii</i> IgM	Positive
HBsAg	Positive
HBeAb	Positive
anti-HIV IgM	Positive
anti-CMV IgM	Positive

No cross-reactivity with the specimens was observed when tested using the NADAL® Syphilis Test.

Interference study

Negative whole blood and low positive serum specimens spiked with the following potentially interfering substances at the concentrations listed below showed no interference with the NADAL® Syphilis Test.

Substance	Concentration	Substance	Concentration
Acetaminophen	0.2 mg/mL	Caffeine	0.2 mg/mL
Acetylsalicylic acid	0.2 mg/mL	Gentisic acid	0.2 mg/mL
Ascorbic acid	20 mg/mL	Albumin	20 mg/mL
Creatine	2 mg/mL	Haemoglobin	0.011mg/mL
Bilirubin	10 mg/mL	Oxalic acid	6 mg/mL

Precision

Repeatability

Repeatability was established by testing 20 replicates of negative, low, medium and high positive specimens using 3 lots of the NADAL® Syphilis tests. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 95.5-100%). The NADAL® Syphilis Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of negative, low, medium and high positive specimens. Testing was performed by 4 operators using 3 independent NADAL® Syphilis test lots at 3 sites on 5 separate days. >99% of the specimens were correctly identified (100/100 correct tests per concentration, 95% confidence interval: 99.3%-100%). The NADAL® Syphilis Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® Syphilis Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

- Centers for Disease Control (CDC). Chlamydia trachomatis infections. Policy guidelines for prevention and control. MMWR Morb Mortal Wkly Rep. 1985 Aug 23; 34 Suppl 3: 535-745.
- Tichonova L, Borisenko K, Ward H, Meheus A, Gromyko A, Renton A. Epidemics of syphilis in the Russian Federation: trends, origins, and priorities for control. Lancet. 1997 Jul 19; 350(9072): 210-3.
- Norgard MV, Chamberlain NR, Swancutt MA, Goldberg MS. Cloning and expression of the major 47-kilodalton surface immunogen of Treponema pallidum in Escherichia coli. Infect Immun. 1986 Nov; 54(2): 500-6.

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15. Bibliografia

- Centers for Disease Control (CDC). Chlamydia trachomatis infections. Policy guidelines for prevention and control. MMWR Morb Mortal Wkly Rep. 1985 Aug 23; 34 Suppl 3: 53S-74S.
- Tichonova L, Borisenko K, Ward H, Meheus A, Gromyko A, Renton A. Epidemics of syphilis in the Russian Federation: trends, origins, and priorities for control. Lancet. 1997 Jul 19; 350(9072): 210-3.
- Norgard MV, Chamberlain NR, Swancutt MA, Goldberg MS. Cloning and expression of the major 47-kilodalton surface immunogen of Treponema pallidum in Escherichia coli. Infect Immun. 1986 Nov; 54(2): 500-6.

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
15. Reference

1. Centers for Disease Control (CDC). Chlamydia trachomatis infections. Policy guidelines for prevention and control. MMWR Morb Mortal Wkly Rep. 1985 Aug 23; 34 Suppl 3: 535-745.
2. Tichonova L, Borisenko K, Ward H, Meheus A, Gromyko A, Renton A. Epidemics of syphilis in the Russian Federation: trends, origins, and priorities for control. Lancet. 1997 Jul 19; 350(9072): 210-3.
3. Norgard MV, Chamberlain NR, Swancutt MA, Goldberg MS. Cloning and expression of the major 47-kilodalton surface immunogen of Treponema pallidum in Escherichia coli. Infect Immun. 1986 Nov; 54(2): 500-6.

Rev. 1, 2023-02-23 JV

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnosticc <i>in vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico- diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Limite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas européias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> - diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur- begränsning	Temperatuurlimiet	Temperatur- begrænsning	Temperatur begrænsning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merkning
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igen
	Prazo de validade	Spotřebujete do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for<n> tester

Our Teams

Germany:**Regensburg**

Tel: +49 941 290 10-0

Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0

Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29

Free Tel: 0800 291 565

Fax: +49 290 10-50

Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18

Free Tel – UK: 0808 234 1237

Free Tel – IRE: 1800 555 080

Fax: +49 290 10-50

France:

France Tel: 0800 915 240

France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720

Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82

Belgium Fax: 0800 747 07

Luxembourg

Lux, Tel: 800 211 16

Lux, Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759

Free Tel: 900 938 315

Fax: +49 941 290 10-50

Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34

Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44

Free Tel: 00 800 491 15 95

Fax: +49 941 290 10-50

Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735

Tel, Verde: 800 849 230

Fax: +49 941 290 10-50

Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600

Free Tel: 0800 0222 890

Fax: +31 70 30 30 775

Free Fax: 0800 024 9519

Nordic countries:**Denmark**

Tel: +31 703075 605

Free Tel: 808 887 53

Finland

Tel: +31 703075 606

Free Tel: 0800 918 263

Free Fax: 0800 918 262

Norway

Tel: +31 703075 605

Free Tel: 800 16 731

Sweden

Tel: +31 703075 605

Free Tel: 020 79 09 06

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Tel: +49 2841 99820-0 • Fax: +49 2841 99820-1

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© 2010 TECHLAB®, Inc. Wszystkie prawa zastrzeżone