
Revisions

BALTSO0191 Version 8.0 Template 4
Inserts

Rev from	Rev to	JOB #
02	03	7564-15

NOTES:

1. BD Catalog Number: 261185
2. Blank (Sheet) Size: Length: 11" Width: 22.5"
3. Number of Pages: 10 Number of Sheets: 1
4. Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" x 1 7/8"
5. Ink Colors: No. of Colors: 3 PMS#: 2755 Blue; 032 Red; Standard Black
6. Printed two sides: Yes ☒ No ☐
7. Style (see illustrations below): # 4

8. See specification control no. NA for material information.
9. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	<div>REVISED BY</div> <div>By Nancy Carlson at 3:04 pm, Jul 23, 2015</div>	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	<div> BD</div> <div>Becton, Dickinson and Company</div> <div>7 Loveton Circle</div> <div>Sparks, MD 21152 USA</div>	
Proofer	<div>PROOFING APPROVED BY</div> <div>By Helaine Jones at 10:37 am, Jul 24, 2015</div>			
Checked By	<div>THIRD EYE BY</div> <div>By Nichole Graham at 2:28 pm, Jul 27, 2015</div>			
Part Number:	L001142	Category and Description Package Insert, Indole Reagent Droppers	Sheet: 1 of 11	A
			Scale: N/A	

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INTENDED USE

BD BBL™ Indole Reagent Droppers (modified Kovacs' reagent) are used in determining the ability of bacteria to produce indole by the deamination of tryptophan.

SUMMARY AND EXPLANATION

The production of indole has been widely used as an aid in differentiating certain genera and species within the same genus.

Indole Reagent may be used for tube or extraction procedures or with a variety of biochemical kits requiring Kovacs' or Ehrlich's reagents.^{1,2}

Each dropper is good for one day's use after the ampule has been broken.

PRINCIPLES OF THE PROCEDURE

A positive reaction indicates the presence of the enzyme tryptophanase that reacts with tryptophan to produce indole. The indole produced reacts in the acid medium with *p*-dimethylaminobenzaldehyde of the indole test reagent to form a quinoidal red-violet compound.

REAGENTS

Indole Reagent Droppers contain 0.5 mL of 5% *p*-dimethylaminobenzaldehyde dissolved in a solution of 25% hydrochloric acid and 75% isobutyl alcohol.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper laboratory procedures in handling and disposing of infectious materials.

Causes burns. Irritating to respiratory system. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Danger

H226 Flammable liquid and vapour. **H302** Harmful if swallowed. **H314** Causes severe skin burns and eye damage. **H335+H336** May cause respiratory irritation. May cause drowsiness or dizziness.

P260 Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P403+P233** Store in a well-ventilated place. Keep container tightly closed. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store at room temperature 15 – 30 °C (59 – 86 °F). Protect from light.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until the expiration date. Change in color of the reagent from light yellow to brown indicates improper storage, which may cause weaker reactions.

PROCEDURE

Material Provided: Indole Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedures and Results:

To use Indole Reagent Droppers, hold upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break glass ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate dropper any further as the plastic may puncture and injury may occur.** Tap the bottom on the tabletop a few times. When ready for use, invert for convenient drop-by-drop dispensing of reagent for the required test procedure.

A. Rapid Tube Method (Arnold and Weaver)³

1. Heavily inoculate a pure culture of organisms to be tested into 1 mL of prewarmed broth of 0.03% tryptophan, 0.1% peptone, and 0.5% dipotassium phosphate (K_2HPO_4) in purified water. Broth should be adjusted to pH 7.4 – 7.8.
2. Incubate tubes in 35 – 37 °C water bath for 2 h.
3. Add contents of one reagent dropper.
4. The appearance of a red ring at the surface of the broth within 30 sec indicates a positive reaction for indole production. A yellow color is negative.

B. Standard Tube Method

- 1. Inoculate a pure culture of the organism to be tested into a culture tube containing 4 mL of **BD Trypticase™** Soy Broth.
- 2. Incubate for 18 – 24 h at 35 – 37 °C.
- 3. Add the contents of one reagent dropper. Shake gently.
- 4. The appearance of a red ring at the surface of the broth within 30 sec indicates a positive reaction for indole production. A yellow color is negative. *Note:* For a more sensitive reaction, 1 – 2 mL of the culture should be aseptically removed from the culture tube and incubated for an additional 24 h. Retest. A negative reaction at 48 h indicates indole has not been produced.

C. Extraction with Xylene

- 1. Inoculate a pure culture of the organism to be tested into 4 mL of appropriate test medium. For *Enterobacteriaceae* use Tryptophan (**BD Trypticase™**) 1% Solution. For others use **BD BBL Indole Nitrite Medium**.
- 2. Incubate for 18 – 24 h at 35 – 37 °C.
- 3. Add 1 mL of xylene. Mix well.
- 4. When layers have begun to separate, add the contents of one Indole Reagent Dropper.
- 5. The appearance of a red ring at the surface of the xylene within 30 sec indicates a positive reaction for indole production. A yellow color is negative. *Note:* For a more sensitive reaction, 1 – 2 mL of the culture should be aseptically removed from the culture tube and incubated for an additional 24 h. Retest. A negative reaction at 48 h indicates indole has not been produced.

D. Other Methods: Indole Reagent Droppers may be used in modifications of the indole test which employ Kovacs' or Ehrlich's indole reagent. For spot tests use **BD BBL DMACA Indole Reagent Droppers**.

User Quality Control:

- 1. Examine the reagent for signs of deterioration (See "Product Deterioration").
- 2. Positive (*Escherichia coli* ATCC® 25922) and negative (*Klebsiella pneumoniae* ATCC 33495) controls should be run simultaneously with the organism to be tested.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

LIMITATIONS OF THE PROCEDURE

Allow at least 30 sec for the color to develop before considering the indole reagent test negative. Pigmented organisms may yield ambiguous results.

PERFORMANCE CHARACTERISTICS⁴

Miller and Wright tested Kovacs indole reagent, *p*-dimethylaminobenzaldehyde (DMABA), Ehrlich indole reagent and *p*-dimethylaminocinnamaldehyde (DMACA) as spot indole reagents to test 359 strains of gram-negative rods growing on sheep blood agar (SBA), **Trypticase** Soy Agar (TSA) and MacConkey Agar (MAC). All reagents were an accurate detector of indole production by members of the *Enterobacteriaceae*, with the exception of *Providencia alcalifaciens*, when organisms were taken from SBA or TSA, but not from MAC. In addition, the overall reactivity of the spot indole procedure with 183 known indole-positive control strains was evaluated. These results are summarized in the table below.

Reagent	No. of positive detected / No. of positive controls (% positive)		
	SBA	TSA	MAC
Kovacs	156/183 (85)	168/183 (92)	41/183 (22)
DMABA	156/183 (85)	168/183 (92)	41/183 (22)
Ehrlich	155/183 (85)	168/183 (92)	41/183 (22)
DMACA	183/183 (100)	183/183 (100)	42/183 (23)

AVAILABILITY

Cat No.	Description
261185	BD BBL™ Indole Reagent Droppers, 50.

REFERENCES

- 1. Balows, A., W.J. Hausler, Jr., K.L. Hermann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
- 2. Shigel, J. 1992. Test methods used in the identification of commonly isolated aerobic gram-negative bacteria. Part 1, Indole test, p. 1.19.13-1.19.16. *In* H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- 3. Arnold, W.M., and R.H. Weaver. 1948. Quick microtechniques for the identification of cultures. J. Clin. Med. 33:1334-1337.
- 4. Miller, J.M., and J.W. Wright. 1982. Spot indole test: evaluation of four reagents. J. Clin. Microbiol. 4:589-592.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Becton Dickinson Pty Ltd.
Benex Limited