


*In-stent restenosis stops here.
RELINE with confidence.*



THE CONTINUING EVOLUTION OF A REVOLUTIONARY DEVICE



VIABAHN®
ENDOPROSTHESIS

PROPATEN
BIOACTIVE SURFACE

PERFORMANCE
through innovation

1996

Original GORE®
HEMOBAHN®
Endoprosthesis
introduced in Europe

2008

GORE® VIABAHN®
Endoprosthesis with
Heparin Bioactive Surface
introduced in Europe

5–8 mm devices decreased
in profile by one French size

2009

Laser technology
enables the new
contoured edge
at proximal end

9–13 mm devices
introduced with
0.035" guidewire
compatibility

2011

GORE® V
Endopro
Heparin
5–8 mm
decreas
by one F

TIP to HUB
deployment
introduced on
6–8 mm devices


2003

25 cm Length:
Longest stent-graft
introduced in EUROPE

2010



*GORE® VIABAHN®
Endoprosthesis with
Bioactive Surface
in devices
in profile
French size*



*Receives CE mark for the
treatment of symptomatic
venous obstruction*

2014

2016

***Radiopaque markers
introduced on 5–8 mm
devices in Europe***

We continue

to evolve the

GORE® VIABAHN®

Endoprosthesis,

demonstrating our

commitment to



providing our

customers with

innovative products.

► Now with Radiopaque Markers for Enhanced Visibility

- Addition of four gold radiopaque markers bonded to the graft at each end of endoprosthesis
- 5–8 mm diameter devices incorporate this change
- Delivery system and profile unchanged



► The Longest Stent-Graft for Endoluminal Bypass

- 25 cm longest length available
- Covers more lesion with one device
- May reduce the need for overlapping devices



➤ Total Endoluminal Bypass

Cover with Confidence

Covers and excludes the diseased irregular tissue of the arterial wall

ePTFE Lining

Provides barrier to in-stent restenosis

Nitinol Stent



Conformable yet durable

Heparin-Bonded Surface

Intended to provide sustained thromboresistance

Lowest Profile Stent-Graft

Reduced profile delivery system makes it even easier to reach and treat challenging SFA lesions




Individual Results May Vary


▶ Contoured Proximal Edge

- Precision laser trimming technology enables manufacturing change
- Excess material at the proximal edge removed
- Contoured edge may improve flow dynamics at proximal end

Canine In Vivo IVUS Examples




30 day



90 day

Non-contoured edge



30 day



90 day

Contoured edge

Excess material removed at the device margin of the contoured edge compared to a non-contoured edge

Animal Acute Examples

Non-contoured edge



Contoured edge



Excess material removed at the device margin of the contoured edge compared to a non-contoured edge

CBAS Heparin Surface

- Intended to provide a thromboresistant surface
- Sustained bioactivity*
- Proprietary end-point covalent bonding

Gore® VIABAHN® Endoprosthesis
with PROPATEN Bioactive Surface




The bioactive luminal surface of a 5 mm diameter Gore® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface appears free of thrombus after two hours in an in vitro blood loop model.

Control Endoprosthesis



The non-bioactive luminal surface of a control endoprosthesis (5 mm diameter) appears covered with thrombus after 90 minutes in the same blood loop model (data on file).

Sustained Bioactivity



Long-term Heparin Activity of Explanted Heparin-bonded
ePTFE Vascular Grafts in a Canine Model*

* Begovac PC, Thomson RC, Fisher JL, Hughson A, Gällhagen A. Improvements in GORE-TEX® Vascular Graft Performance by Carmeda® BioActive Surface Heparin Immobilization. *European Journal of Vascular and Endovascular Surgery* 2003;25(5):432-437.

Sizing Table

TIP to HUB Device Deployment – 0.014" or 0.018" Guidewire Compatibility (With radiopaque markers)

Device Sizing		Introducer Sheath (Fr)					RECOMMENDED BALLOON DIAMETER FOR DEVICE TOUCH-UP ³ (mm)
ENDOPROSTHESIS LABELED DIAMETER ¹ (mm)	RECOMMENDED VESSEL DIAMETER ² (mm)	2.5 cm DEVICE LENGTH ¹	5 cm DEVICE LENGTH ¹	10 cm DEVICE LENGTH ¹	15 cm DEVICE LENGTH ¹	25 cm DEVICE LENGTH ¹	
5	4.0 – 4.7	6	6	6	6	6	5
6	4.8 – 5.5	6	6	6	6	6	6
7	5.6 – 6.5	7	7	7	7	7	7
8	6.6 – 7.5	7	7	7	7	7 ⁵	8

TIP to HUB Device Deployment – 0.035" Guidewire Compatibility (Radiopaque markers on 5–8 mm devices)

Device Sizing		Introducer Sheath (Fr)					RECOMMENDED BALLOON DIAMETER FOR DEVICE TOUCH-UP ³ (mm)
ENDOPROSTHESIS LABELED DIAMETER ¹ (mm)	RECOMMENDED VESSEL DIAMETER ² (mm)	2.5 cm DEVICE LENGTH ¹	5 cm DEVICE LENGTH ¹	10 cm DEVICE LENGTH ¹	15 cm DEVICE LENGTH ¹	25 cm DEVICE LENGTH ¹	
5	4.0 – 4.7	7	7	7	7	7	5
6	4.8 – 5.5	7	7	7	7	7	6
7	5.6 – 6.5	8	8	8	8	8	7
8	6.6 – 7.5	8	8	8	8	8	8
9	7.6 – 8.5	–	9	9	9	–	9
10	8.6 – 9.5	–	11 ⁴	11 ⁴	11 ⁴	–	10
11	9.6 – 10.5	–	11	11	–	–	12
13	10.6 – 12.0	–	12	12	–	–	14

¹ Labeled device diameters and lengths are nominal.

² Recommended endoprosthesis compression within the vessel is approximately 5 – 20%.

³ For the 11 and 13 mm diameter devices, balloon inflation pressure should not exceed 8 atm.

⁴ The 10 mm diameter device is compatible with the following 10 Fr introducer sheaths: CORDIS® AVANTI® Sheath Introducer, Boston Scientific SUPER SHEATH Introducer Sheath, B. Braun INTRADYN Tear-Away Introducer Sheath.

⁵ The 8 mm x 25 cm device is not compatible with the 7 Fr COOK® CHECK-FLO® FLEXOR® Sheath.



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)
00800.6334.4673 (Europe)
800.437.8181 (United States)
928.779.2771 (United States)

goremedical.com

The GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface is sold in some markets as the GORE® VIABAHN® Endoprosthesis with Heparin Bioactive Surface.

Products listed may not be available in all markets.

COOK®, CHECK-FLO®, and FLEXOR® are trademarks of Cook Medical, Inc.
CORDIS® and AVANTI® are trademarks of Cordis Corporation.

GORE®, GORE-TEX®, HEMOBAHN®, PERFORMANCE THROUGH INNOVATION, PROPATEN, TIP to HUB, VIABAHN®, and designs are trademarks of W. L. Gore & Associates. CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.
© 2013–2016 W. L. Gore & Associates, Inc. AL0342-EU7 APRIL 2016

