www.thermofisher.com

Thermo
SCIENTIFIC

Tryptone Soya Broth USP

REF TV5002E

Intended Use

Tryptone Soya Broth USP (TV5002E) is a highly nutritious general-purpose medium for the growth of both aerobic and anaerobic, fastidious and non-fastidious bacteria and fungi from clinical samples (wounds, skin, respiratory secretions, genito urinal secretions etc.).

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Pathogenic aerobic and anaerobic microorganisms encompass a broad range of microorganisms, which can often be found as part of the normal skin and mucosal flora of humans¹. Infections from fastidious aerobic and anaerobic microorganisms can have a range of manifestations ranging from superficial to life-threatening¹.

Principle of Method

The device is intended to be used as a highly nutritious general-purpose medium for the growth of both aerobic and anaerobic, fastidious and non-fastidious bacteria and fungi from clinical samples. Due to the inclusion of both Tryptone and Soya Peptone, the medium will support a luxuriant growth of many fastidious organisms without the addition of serum, etc. Tryptone Soya Broth may be used for the cultivation of aerobes and facultative anaerobes, including some fungi. The addition of a small amount of agar renders the broth suitable for the cultivation of obligatory anaerobes, such as *Clostridium* spp. Tryptone Soya Broth has a greater ability to resuscitate heated spores of *Geobacillus stearothermophilus* than Dextrose Tryptone Broth. Tryptone Soya Broth is recommended in the USP (Casein soya bean digest broth) for the recovery of organisms after sterilisation processes.

Typical Formula

	<u>grams per litre</u>
Pancreatic digest of casein	17.0
Papaic digest of soyabean meal	3.0
Sodium chloride	5.0
Dibasic potassium phosphate	2.5
Glucose	2.5

Physical Appearance

Colour	Maize yellow
Clarity	Transparent
Fill weight	10.0 ± 0.5g
pH	7.3 ± 0.2

Materials Provided

TV5002E: 50 tubes of Tryptone Soya Broth USP

Each tube should only be used once.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–25°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or tubes.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 1, ID 9, B 39 and B 57.

Procedure

- Allow product to equilibrate to room temperature.
- Inoculate broth.
- Incubate aerobically:
For up to 3 days at 30–35°C for *Bacillus subtilis*.
For up to 5 days at 20–25°C for *Candida albicans* and *Aspergillus brasiliensis*.
For up to 2 days at 55–57°C for *Geobacillus stearothermophilus*.
- Inspect for turbidity.

Interpretation

Turbid growth indicates *Bacillus subtilis* or *Geobacillus stearothermophilus*. Turbid, flocculent growth indicates *Candida albicans*. White mycelia indicates *Aspergillus brasiliensis*.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Positive Controls Inoculum level: 10 – 100 cfu Incubation Conditions: up to 3 days @ 30-35°C aerobic	
<i>Bacillus subtilis</i> ATCC® 6633™	Turbid growth
Incubation Conditions: up to 5 days @ 20-25°C aerobic	
<i>Candida albicans</i> ATCC® 10231™	Turbid, flocculent growth
<i>Aspergillus brasiliensis</i> ATCC® 16404™	White mycelia
Incubation Conditions: up to 2 days @ 55-57°C aerobic	
<i>Geobacillus stearothermophilus</i> ATCC® 7593™	Turbid growth

Limitations

Organisms with atypical enzyme patterns may give anomalous reactions in Tryptone Soya Broth USP.

The device is for the maintenance of fastidious aerobic and anaerobic microorganisms. Tryptone Soya Broth USP (TV5002E) is used in a diagnostic workflow as an enrichment step to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct recovery of microorganisms from clinical samples is confirmed by the inclusion of 4 well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the device. The precision of Tryptone Soya Broth USP (TV5002E) was demonstrated by an overall pass rate of 100% obtained for the product over one month of testing (25.05.2022–23.06.2022; 10 batches). This shows that the performance is reproducible.

This device is tested in-house as part of the QC process since the product was launched in 2003. For target organisms, when using 10-100 colony forming units (cfu) and the test medium is inoculated. Following incubation for up to 3 days in aerobic condition for *Bacillus subtilis* and after incubation up to 5 days, aerobic for *Candida albicans* *Aspergillus brasiliensis* and up to 2 days, aerobic for *Geobacillus stearothermophilus*, the user can recover organisms with growth as listed in this document

Bibliography

- Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>

Symbol Legend

Symbol	Definition
	Catalogue number

	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
Made in Germany	Made in Germany

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

Version	Date of modifications introduced
1.0	2022-08-31 Original Document

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Oxoid Deutschland GmbH, Am Lippeglacis