## SK

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## EN

## SUMMARY AND EXPLANATION

The OXI is detection strip designed for detection of cytochromoxidase. Detection strip shows colour change on the base of metabolism of examined samples.

## PRINCIPLE

The activity of cytochrome oxidase is attended by change of colour of detection zone to blue.

## DIAGNOSTIC KIT CONTAINS

- 50 test strips for 50 determinations
- information leaflet

## WARNINGS AND SAFETY PRECAUTIONS

- **For in vitro diagnostics.**
- **For professional use only.**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and obvious valid regulations must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

## STORAGE

Store OXI in dark at +2 to +8°C. Expiration date is placed on the label of each package.

## REQUIRED MATERIAL

Inoculation loop

## RECOMMENDED PROCEDURE

1. Use pure 24 hour bacterial culture. Pick a 1-2 colonies from the medium by an inoculation loop and spread it on diagnostic zone of strip.
2. Evaluate the test immediately, max. till one minute. Later evaluation can cause false positive result.

## INTERPRETATION OF REACTION

**Blue colour indicates positive result of the OXI test.**

**For the need of own proof function, use following recommend bacterial strains:**

Bacterial strain	CCM / ATCC	OXI	Colour
<i>Pseudomonas aeruginosa</i>	3955 / 27853	+	blue
<i>Escherichia coli</i>	3954 / 25922	-	without change

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz)

**TIP: use negative control – follow recommended procedure, but without bacterial culture.**

Kit components doesn't contain hazardous substances.

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