remel

COAGULASE PLASMA

(Rabbit Plasma w/ EDTA)

EN

REF	R21050, Coagulase Plasma	5 ml/Vial
REF	R21060, Coagulase Plasma	6 x 5 ml/Vial
REF	R21051, Coagulase Plasma	15 ml/Vial
REF	R21052, Coagulase Plasma	25 ml/Vial

INTENDED USE

Remel Coagulase Plasma is a reagent for use in qualitative procedures for the detection of coagulase enzyme in staphylococci. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

SUMMARY AND EXPLANATION

In 1940, Fairbrother and Chapman et al. reported pathogenic staphylococci could be identified based on the ability to coagulate plasma.^{1,2} Chapman found rabbit plasma to be superior to other types of plasma in terms of clotting activity.³ Bayliss and Hall recommended replacing the citrate anticoagulant with EDTA to avoid false-positive clots caused by bacteria which utilize citrate.⁴

PRINCIPLE

The enzyme coagulase acts on a constituent of rabbit plasma (coagulase reacting factor) to produce a thrombin-like substance. This substance activates fibrinogen to form fibrin which results in the formation of a fibrin clot. Coagulase is present in two forms: bound coagulase or clumping factor remains attached to the cell wall of the organism; free coagulase is an extracellular enzyme produced when the organism is cultured in broth. Bound coagulase is detected in the slide test, while the tube test will detect bound and free coagulase.

REAGENTS (CLASSICAL FORMULA)*

Sodium Chloride (CAS 7647-14-5)	4.5 g
Rabbit Plasma w/ EDTA	500.0 ml
Demineralized Water (CAS 7732-18-5)	500.0 ml

*Adjusted as required to meet performance standards.

PRECAUTIONS

This product is for *In Vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully.

Please refer to the Safety Data Sheet (SDS) on company website and product labeling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device.

STORAGE

Store lyophilized product in its original container at 2-8°C until used. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

PRODUCT DETERIORATION

This product should not be used if (1) the plasma is clotted upon rehydration, (2) the product is contaminated, (3) the expiration date has passed, or (4) there are other signs of deterioration.

SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.⁵

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Glass slides, (7) Test tubes w/ caps, (8) Pipettes, (9) Sterile physiological saline, sterile demineralized water.

REAGENT PREPARATION

Rehydrate lyophilized Coagulase Plasma with sterile demineralized water according to the volume indicated on the product label. Invert the vial and mix well. Dispense 0.5 ml aliquots of rehydrated solution into clean tubes. Aliquots may be tightly capped and frozen at -20°C or below for up to 1 month or refrigerated at 2-8°C for 5 days.⁹ Frozen aliquots of Coagulase Plasma should not be re-frozen after thawing.

PROCEDURE

Test only 18-24-hour colonies of catalase-positive, gram-positive cocci that are morphologically characteristic of staphylococci on Gram stain and plated media. The test isolate should be growing on a nonselective medium that does not have a high concentration of salt. Positive and negative control organisms should be tested concurrently with each test run of patient isolates.

Note: When performing the slide test, verify the test isolate does not auto-agglutinate by observing for clumping of the test isolate in the drop of saline or demineralized water prior to addition of Coagulase Plasma.

Slide Test: (Detects bound coagulase, only)

1. Place a drop of demineralized water or physiological saline on a clear, clean, glass slide.
2. Emulsify a loopful of the test isolate from isolated colonies growing on a nonselective medium in the drop of water or saline. (Verify the test isolate does not auto-agglutinate.)
3. Gently mix a loopful of Coagulase Plasma into the staphylococcal suspension.
4. Observe for immediate formation of macroscopic precipitate in the form of white clumps.
5. Slide test reactions must be read quickly because false-positive results may appear with reaction times greater than 10 seconds.⁵

Note: Some strains of *S. aureus* are negative with the slide coagulase test. All negative slide tests must be re-tested by the tube test.^{5,6,8}

Tube Test: (Detects bound and free coagulase)

1. Add 0.5 ml of Coagulase Plasma to a clean test tube.
2. Add 0.5 ml from a pure broth culture or a large loopful of a pure culture from a nonselective medium.

- Mix thoroughly to emulsify the test isolate.
- Incubate at 35-37°C in a water bath or incubator.
- Observe every 30 minutes for clotting by gently slanting the tube. Do not shake.
- If no clot is visible after 4 hours, leave in water bath or place in a 35-37°C incubator overnight (24 hours). An optional method is to leave the tubes overnight (up to 24 hours) at 25°C.⁵⁻⁸

INTERPRETATION

Slide Test:

- Positive Test - Marked clumping within 10 seconds
- Negative Test - No clumping, suspension remains homogenous; confirm by tube test before reporting result as negative

Tube Test:

- Positive Test - Clot formation
- Negative Test - No clot, suspension remains homogenous

QUALITY CONTROL

All lot numbers of Coagulase Plasma have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL

<i>Staphylococcus aureus</i> ATCC® 6538	Ambient, 24 h @ 35-37°C	Positive tube test at 4 and 24 hours
<i>Staphylococcus aureus</i> ATCC® 29213	Ambient, 24 h @ 35-37°C	Positive tube test at 4 and 24 hours
<i>Staphylococcus aureus</i> ATCC® 25923	Ambient, 24 h @ 35-37°C	Positive tube test at 4 and 24 hours
<i>Staphylococcus aureus</i> ATCC® 33862	Ambient, 24 h @ 35-37°C	Positive tube test at 4 and 24 hours
<i>Staphylococcus epidermidis</i> ATCC® 12228	Ambient, 24 h @ 35-37°C	Negative tube test at 4 and 24 hours

INCUBATION

RESULTS

LIMITATIONS

- False-positive results may occur if the test isolate is removed from agar containing high concentrations of salt. Use only isolates grown on nonselective media.⁵⁻⁷
- Do not shake or agitate the tube during incubation of the tube test. This may cause a breakdown of the clot, which will not reform with additional incubation.
- When performing the tube test, it is necessary to examine the tube every 30 minutes for the first 4 hours because some strains of *S. aureus* produce fibrinolysin which lyses the clot early in the incubation period.⁶
- False-negative coagulase reactions may occur if the test isolate is older than 18-24 hours or if there is scant growth.⁶
- Some strains of *Staphylococcus* that produce free coagulase (tube test) do not form bound coagulase (slide test). Therefore, all isolates with a negative slide test must be re-tested with the tube test before reporting as negative for coagulase.^{5,6,8}

BIBLIOGRAPHY

- Fairbrother, R.W. 1940. J. Pathol. Bacteriol. 50:83-88.
- Chapman, G.H., C. Berens, and M.H. Stiles. 1940. J. Bacteriol. 43:431-439.
- Chapman, G.H. 1944. J. Bacteriol. 47:211-212.
- Bayliss, B.G. and E.R. Hall. 1965. J. Bacteriol. 89:101-105.

- Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Cowan, S.T. and K.J. Steel. 1970. Manual for the Identification of Medical Bacteria. Cambridge Univ. Press, New York, NY.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Garcia, L.S. 2010. Clinical Microbiology Procedures Handbook. 3rd ed. ASM Press, Washington, D.C.

Symbol Legend

	Catalog Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
	Use By (Expiration Date)
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Authorized European Representative
	Manufacturer
	UK Conformity Assessment
	European Conformity Assessment

Remel Inc
12076 Santa Fe Trail Drive
Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
(800) 255-6730
International: (913) 888-0939

ATCC and ATCC catalogue marks are registered trademark of American Type Culture Collection.

CAS (Chemical Abstracts Service Registry No.)

All other trademarks the property of Thermo Fisher Scientific Inc and its subsidiaries.

For technical information contact your local distributor.

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.