PL

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Instruction for Use

ANTISERUM FOR SLIDE AGGLUTINATION OF *SALMONELLA*

Intended use

The antisera for *Salmonella* slide agglutination assay are an *in vitro* diagnostic medical device. The antisera are used for qualitative identification tests of *Salmonella* serotypes isolated on bacteriological media from human or animal samples or other materials.

Description

The antisera for *Salmonella* slide agglutination assay are composed of the absorbed polyvalent or monovalent rabbit antisera, which contain antibodies against group and single factors of somatic, flagellar or capsular Vi antigens. The antisera allow for serological identification of *Salmonella* serovars occurring in Poland and worldwide. Antisera for phase inversion are designed for the phase induction in 2-phase bacteria with inhibited expression of one phase of H-antigens. It can be obtained by culturing *Salmonella* on swarm agar containing antisera for phase inversion.

The ready to use antisera are sterile and sold in glass bottles with a dropper.

Antisera are designed for serological diagnostics of bacteria belonging to *Salmonella* genus.

Product for professional laboratory application, should only be used by qualified personnel following basic aseptic and antiseptic principles.

Composition

Antisera for slide agglutination of rod shaped *Salmonella* are produced from polyclonal antibodies obtained from animals vaccinated with inactivated bacterial strains. In order to contain only specific antibodies, they undergo the process of absorption and are suspended in 0,85% NaCl. Quality control is carried out with a set of over 180 *Salmonella* strains. Antisera are sterile and preserved with 0,01% thiomersal.

Properties

Identification of *Salmonella* serovars relies on the recognition of the antigen structure of the tested strain. Serological characterization is based upon the detection of 3 antigenic components: somatic "O" antigens, flagellar "H" antigens or capsular "Vi" antigen. Capsular "Vi" antigen can be found only in certain strains e.g. S. Typhi, and can be detected by slide agglutination assay with anti-Vi antiserum. Bacterial surface antigens form visible complexes with specific antibodies in the serum (agglutination reaction).

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Method

Preparation of samples

Samples should be collected in sterile containers, transported to the laboratory and processed according to the guidelines for required type of samples. Clearly distant, typical colonies with no signs of contamination should be chosen to perform the slide agglutination assay.

Preparation of tested strains

Culture the tested strains at the required temperature, ranging from +34°C to +38°C for 24±3 hours on nutrient agar media (for identification of somatic antigens use 1,5% solid medium, for flagellar antigens use 0,5% semi-solid medium). (1,5% agar medium produced by Immunolab, catalog No. PG002, PGS02 and 0,5% Sven Gard swarm agar by Immunolab, catalog No. PG001, PGS01 are recommended).

Performance of slide agglutination assay

Antisera should be brought to room temperature (from +18 to +25°C) prior to using them. If turbidity of antiserum is observed, it should be centrifuged. In order to maintain the sterility of the product (the lack of which could make it impossible to perform the assay or make correct interpretation of the results), the test should be performed in conditions that prevent any contamination of the product.

First of all, before performing slide agglutination assay, a negative control of tested strain with 3% NaCl solution should be carried out. The occurrence of agglutination in a 3% NaCl solution indicates that the tested strain is in phase R (rough strain) and it is not suitable for serological identification in slide agglutination assay.

1. Place a drop of antiserum on the degreased glass slide.
2. With the use of a sterile wire loop/baguette, collect bacteria from the medium and place them next to the drop. While spreading the strain on the slide, mix it with the serum to form a homogeneous suspension.
3. While gently rocking the slide with circular motions for **10 - 30 seconds (max. till 1 minute!)**, observe the reaction. Be careful that the drop does not drip from the slides.

Agglutination is more visible if the reaction is observed against a dark background with a magnifying glass with 5-fold magnification.

Interpretation of the results

To define the serotype of the tested strain, its antigenic formula needs to be determined with the use of antisera containing antibodies against somatic and flagellar antigens according to the White-Kauffmann Le Minor-classification scheme.

Nature of the reaction should be taken into account: agglutination visible as little lumps is characteristic for the somatic antigens and Vi antigen, as cloudlets for flagellar antigens. Presence of agglutination is a positive reaction. Lack of agglutination is a negative reaction.

Interpretation of agglutination should be made according to the following scale:

- +++ agglutination present, distributed throughout the transparent drop or on its periphery,
- ++ agglutination present, distributed in a semitransparent drop,
- + small agglutination present on the drop's periphery or at the bottom in a milky white drop. Such results cannot be interpreted as positive in serological diagnostics. Such a situation occurs mostly with biphasic strains that have one strongly developed phase and mainly dominant phase is expressed. In order to determine antigens' composition, the dominant phase should be inhibited,
 - (-) lack of agglutination, milky white drop.
 - ! Do not let the drop dry out, as it may be interpreted as a false negative result.

An instructional video with detailed information how to perform the test and how to interpret the obtained results is available on our website: <http://immunolab.com.pl/film>

Performance of a phase inversion

In the case of diphasic strains, when the tested strain gives a positive reaction with only one out of two phases (dominant phase) and the second phase cannot be determined (latent phase), it is necessary to perform phase inversion in order to impede the dominant phase and induce the latent phase.

In order to do that, the following steps should be taken:

1. Prepare the semi-solid 0,5% Sven Gard swarm agar and allow it to cool to ~45°C.
2. Place a couple of drops (100-200µl) of inhibitory antiserum into the middle of a small Petri dish (Ø approx. 5 cm). It is necessary to pay attention and not to touch the surface of a Petri dish with a dropper.
3. Add approximately 10 ml of 0,5% Sven Gard swarm agar and mix gently with an antiserum by circular motions.
4. Allow the medium to cool and solidify at room temperature. Do not dry the medium!
5. With the use of a sterile wire loop, inoculate the strain to be inverted, at a single point in the center of a Petri dish. For inversion use bacteria previously grown on 0,5% Sven Gard swarm agar.
6. Incubate the inoculated plate overnight, at temperatures ranging from +34 to +38°C. Do not place the Petri dish upside down!
7. The swarmed strain from the peripheral parts of a plate can be used to identify the second phase in slide agglutination assay.

Performance of a phase inversion - an alternative method

1. Prepare a small Petri dish (Ø approx. 5 cm) with approx. 10 ml of sterile, semi-solid 0,5% Sven Gard swarm agar. Store the plate with the lid up. Do not dry the medium!
2. Add a couple of drops (100-200µl) of inhibitory antiserum into the middle of a plate with a nutrient medium. It is necessary to pay attention and not to touch the surface of a medium with a dropper.
3. Spread the antiserum over the surface with a sterile (glass) spatula.
4. With the use of a sterile wire loop, inoculate the strain to be inverted at a single point in the center of a dish. For inversion use bacteria previously grown on 0,5% Sven Gard swarm agar.
5. Incubate the inoculated plate overnight, at temperatures ranging from +34 to +38°C. Do not place the Petri dish upside down!
6. The swarmed strain from the peripheral parts of a plate can be used to identify the second phase in slide agglutination assay.

! Phase inversion is not always achieved at the first attempt. When necessary, the procedure should be repeated before concluding that the organism has no alternative phase 2 (monophasic strains). *The bacteria for further inoculation should be collected from the plate used during the first attempt.*

An example:

S.Typhimurium 1,4,[5],12:i:1,2

Dominant H2 phase → use antiserum H:2 for phase inversion to suppress dominant phase and, enhance expression of phase H:i

Dominant H:i phase → use antiserum H:i for phase inversion to suppress dominant phase and enhance expression of phase H:2

User quality control

When using the antiserum, it is recommended to test control the method with the reference *Salmonella* strains (both positive and negative controls) in order to verify the correctness of the method and the reactions of the antiserum, proceeding in the same way as in the case of tested strains.

Limitations of the method

In exceptional cases, there is a possibility of cross-reactions with other species of Enterobacteriaceae (e.g. strains of *Citrobacter* spp., *Hafnia alvei*, *Proteus* spp., *E. coli* or *Shigella* spp.) due to their high similarity between antigens or their phylogenetic relationship. For this reason, biochemical tests are necessary to determine whether the tested strain belongs to the *Salmonella* genus. Only pure cultures of *Salmonella* should be used for serological diagnostics.

Materials required to perform the assay, which are not included in the set

Dark background, inoculation loop, magnifying glass, glass slides, solution of 3% NaCl, container for contaminated materials, sterile (glass) spatula, Petri dish, the White-Kauffmann-Le Minor classification scheme, an instructional video (available on our website: <http://immunolab.com.pl/film>).

Validation and verification of the method

sensitivity	97% (363/376) (number of true positive results/ sum of true positive and false negative results)
specificity	99% (594/598) (number of true negative results/ sum of true negative and false negative results)
repeatability	98% (795/810) (number of repeated compatible results/total number of results)
reproducibility	98% (792/810) (number of reproduced compatible results/ total number of results)

Storage conditions and precautions

- Antisera should be stored at temperatures ranging from **2°C to 8°C**.
- Transport up to 3 days in temperatures not exceeding 25°C does not deteriorate the properties of the product.
- **DO NOT FREEZE!**
- Protect from light.
- Do not use after the expiration date stated on the package.
- The antiserum is a reusable product.
- The expiration date applies to the product in intact packaging, stored in accordance with the instructions, and the product after opening, used in accordance with the intended use, in the conditions described in the instructions and stored after opening in accordance with the instructions.
- Do not use antiserum with altered physico-chemical characteristics indicating microbiological contamination.
- Do not sterilize antiserum prior to use.
- Turbidity and/or precipitates not caused by microbiological contamination may appear in the antiserum. The antiserum is still suitable to be used, but must be centrifuged (8000rpm, 15 min) in order to give correct agglutination results.

Precautions

The product does not pose a health threat to humans and the environment, when used according to the instructions.

Aseptic working conditions must be observed and precautions against microbiological hazards should be taken during all procedures.

The instructions for use must be strictly followed.

Any serious incident involving the device should be reported to the manufacturer and the competent authority of the Member State in which the user resides or has his place of business.

Precautions when working with infectious material:

- Rabbit antisera contain material of animal origin and should be handled as potentially infectious.
- Performing a slide agglutination assay involves working with pathogenic bacteria, therefore all necessary rules applicable when working with infectious microbiological material must be observed. Wear protective clothing and gloves. Samples collected to perform the assay as well as the material from the slides should be protected against access by third parties, due to the possibility of infection. Slides used in the assay containing antiserum and infectious material, as well as other contaminated materials, should be autoclaved or handled in accordance with applicable laboratory rules.

The content of the package unit

A single antiserum bottle with a dropper contains:

- 5 milliliters of antiserum allowing to perform approx. 125 tests, or
- 3 milliliters of antiserum allowing to perform approx. 75 tests, or
- 2 milliliters of antiserum allowing to perform approx. 50 tests, or
- 1 milliliter of antiserum allowing to perform approx. 25 tests.

A single phase inversion antiserum bottle with a dropper contains:

- 5 milliliters of inhibitory antiserum allowing to perform approx. 25 phase inversion tests, or
- 3 milliliters of inhibitory antiserum allowing to perform approx. 15 phase inversion tests, or
- 2 milliliters of inhibitory antiserum allowing to perform approx. 10 phase inversion tests, or
- 1 milliliter of inhibitory antiserum allowing to perform approx. 5 phase inversion tests.

Manufacturer

Research and Development Department of Salmonella Center

"IMMUNOLAB" Ltd.

Address: ul. Kładki 24, 80-822 Gdańsk, Poland

Phone/Fax: 0048 58 781-44-91

E-mail: info@immunolab.com.pl

www.immunolab.com.pl

Explanation of the symbols used

Manufacturer

Expiration date : YYYY-XX-Z (year-month-day)

Batch code (Lot)

Consult instructions for use

Catalogue number

Storage temperature limit

In vitro diagnostic medical device

Transport temperature limit up to 3 days