www.thermofisher.com

Thermo
SCIENTIFIC

| | |
|-----------------|-------|
| Sodium Chloride | 2.5 |
| L-Cysteine | 0.5 |
| Resazurin | 0.001 |
| Agar | 0.075 |

Thioglycollate Medium USP

REF TV5001D

Intended Use

Thioglycollate Medium USP (TV5001D) is a medium for the maintenance of fastidious aerobic and anaerobic microorganisms.

Thioglycollate Medium USP is used in a diagnostic workflow as an enrichment step to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections. The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Fastidious aerobic and anaerobic microorganisms encompass a broad range of microorganisms, which can often be found as part of the normal skin and mucosal flora of humans¹. *Clostridium* spp. may manifest and cause infection in the tissues of immunocompromised individuals or those undertaking a course of antibiotics² whilst *Staphylococcus aureus* is also a nosocomial pathogen, with invasive medical interventions such as catheterization or immunosuppressant treatment, contributing to the incidence of peripheral intravenous catheter (PIVC) failure and catheter-related bloodstream infections (CRBSI)^{1,3}. Moreover, the emergence of antibiotic resistance complicates treatment and impacts patient outcomes, especially in immunocompromised or hospitalized individuals⁴.

Infections from fastidious aerobic and anaerobic microorganisms can have a range of manifestations ranging from superficial to life-threatening^{1,2,5}. According to the varieties of toxins produced by *Clostridium perfringens*, a patient may suffer from symptoms comparable to that of food poisoning, or of the fatal clostridial myonecrosis². Additionally, infection with *Staphylococcus aureus*, a common commensal of the human upper respiratory tract and skin microbiomes, can cause superficial skin infections such as cellulitis and abscesses to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis¹. Similarly, the opportunistic infection of *Pseudomonas aeruginosa* in the clinical environment may lead to pneumonia and septic shock, especially in CF and neutropenic patients⁵.

Principle of Method

Thioglycollate Medium USP (TV5001D) is a medium for the maintenance of fastidious aerobic and anaerobic microorganisms in the performance of sterility tests. Thioglycollate Medium USP contains yeast extract combined with tryptone to provide a source of nitrogen, essential vitamins and amino acids necessary for the growth of bacteria. The presence of sodium thioglycollate as a reducing agent maintains a low oxygen concentration in the media.

Typical Formula

| | grams per litre |
|-----------------------|-----------------|
| Yeast extract | 5.0 |
| Tryptone | 15.0 |
| Glucose | 5.5 |
| Sodium thioglycollate | 0.5 |

Physical Appearance

| | |
|-------------|--------------|
| Colour | Lemon Yellow |
| Clarity | Transparent |
| Fill weight | 9.0 ± 0.5mL |
| pH | 7.1 ± 0.2 |

Materials Provided

TV5001D: 50 x 9ml Thioglycollate Medium USP tubes..... 50

Each tube should only be used once.

Materials Required but Not Supplied

- 1) Inoculating loops
- 2) Swabs
- 3) Collection containers
- 4) Incubators
- 5) Quality control organisms

Storage

- Store product in its original packaging at 20–25°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- If the upper portion of the medium is red because of oxidation, anaerobic conditions can be restored by reheating for 10 minutes in boiling water or steam. Do not reheat more than once. Following reheating, if more than one-third of the medium is oxidised then it should be discarded.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

incubating the device at 32°C ± 1°C for up to 3 days, the user can recover organisms with the appropriate morphology as listed in this document.

Specimen Collection, Handling and Storage

There are no specimen collection and/or transport materials provided with the accessory device. Specimens should be collected and handled following recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q5.

Procedure

- (1) Inspect the product for oxidation before use
- (2) Allow product to equilibrate to room temperature.
- (3) Inoculate the specimen onto the medium using a standard loop.
- (4) Incubate plates at 30-35°C for up to 2 weeks.
- (5) Growth of strictly aerobic bacteria can be improved by slightly loosening the cap
- (6) Visually inspect tubes to assess colony growth and colour under good lighting.

Interpretation

The presence of turbid or flocculent colonies indicate the presence of fastidious organisms

- Turbid growth indicates *Clostridium sporogenes*, *Clostridium perfringens* or *Pseudomonas aeruginosa*
- Flocculent growth indicates *Staphylococcus aureus*.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 - 34 h @ 36° ± 1°C aerobic for *Clostridium perfringens* and 3 days @ 32 ± 1°C for *Clostridium sporogenes*, *Staphylococcus aureus* and *Pseudomonas aeruginosa*

| Control Strain | Growth |
|---|-------------------|
| <i>Clostridium sporogenes</i> ATCC® 19404 | Turbid growth |
| <i>Clostridium perfringens</i> ATCC® 13124 | Turbid growth |
| <i>Staphylococcus aureus</i> ATCC® 6538 | Flocculent growth |
| <i>Pseudomonas aeruginosa</i> ATCC® 9027 | Turbid growth. |

Limitations

Organisms with atypical enzyme patterns may give anomalous reactions on Thioglycollate Medium USP.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of fastidious microorganisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of Thioglycollate medium USP (TV5001D) was demonstrated by an overall pass rate of 100% obtained for the product over 10 batches of testing (January 2022) of testing. This shows that the performance is reproducible.

The devices are tested in-house as part of the QC process. For target organisms, when using 100 cfu inoculum of *Clostridium perfringens* and incubating the device at 36°C ± 1°C for 18-34 hours, the user can recover organisms with the appropriate morphology as listed in this document. When using 100 cfu inoculum of *Clostridium sporogenes*, *Staphylococcus aureus* and *Pseudomonas aeruginosa* and

Bibliography

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4).
<https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2016. 'Identification of Clostridium species'. ID 8 (4).
<https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>
3. Blanco-Mavillard I., Rodriguez-Calero M.A., De Pedro-Gomez J., Parra-Garcia G., Fernandez-Fernandez I. and Castro-Sanchez E. 2019. 'Incidence of peripheral intravenous catheter failure among inpatients: Variability between microbiological data and clinical signs and symptoms'. *Antimicrobial Resistance and Infection Control* 8(1):124.
<https://www.scopus.com/inward/record.uri?eid=2-s2.0-85069782323&doi=10.1186%2fs13756-019-0581-8&partnerID=40&md5=09596ab961e4fc1b1c777b8eca14664c>
4. Foster, Timothy J. 1996. 'Staphylococcus'. In Medical Microbiology. 4th Edition., 4th ed. Galveston.
5. Public Health England. 2015. 'Identification of Pseudomonas species and other Non-Glucose Fermenters'. ID 17 (3).
<https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>

Symbol Legend

| Symbol | Definition |
|--------|---|
| | Catalogue number |
| | In Vitro Diagnostic Medical Device |
| | Batch code |
| | Temperature limit |
| | Use-by date |
| | Keep away from sunlight |
| | Do not re-use |
| | Consult instructions for use or consult electronic instructions for use |
| | Contains sufficient for <n> tests |

| | |
|---|--|
|  | Do not use if packaging damaged and consult instructions for use |
|  | Manufacturer |
|  | Authorized representative in the European Community/ European Union |
|  | European Conformity Assessment |
|  | UK Conformity Assessment |
|  | Unique device identifier |

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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For technical assistance please contact your local distributor.

Revision Information

| Version | Date of issue and modifications introduced |
|---------|--|
| 1.0 | 2022-07-01. New document. (LIVE) |

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Literatura

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4).
<https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2016. 'Identification of Clostridium species'. ID 8 (4).
<https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>
3. Blanco-Mavillard I., Rodriguez-Calero M.A., De Pedro-Gomez J., Parra-Garcia G., Fernandez-Fernandez I. and Castro-Sanchez E. 2019. 'Incidence of peripheral intravenous catheter failure among inpatients: Variability between microbiological data and clinical signs and symptoms'. Antimicrobial Resistance and Infection Control 8(1):124.
<https://www.scopus.com/inward/record.uri?eid=2-s2.0-85069782323&doi=10.1186%2f13756-019-0581-8&partnerID=40&md5=09596ab961e4fc1b1c777b8eca14664c>
4. Foster, Timothy J. 1996. 'Staphylococcus'. In Medical Microbiology. 4th Edition., 4th ed. Galveston.
5. Public Health England. 2015. 'Identification of Pseudomonas species and other Non-Glucose Fermenters'. ID 17 (3).
<https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>

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Bibliographie

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4).
<https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2016. 'Identification of Clostridium species'. ID 8 (4).
<https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>
3. Blanco-Mavillard I., Rodriguez-Calero M.A., De Pedro-Gomez J., Parra-Garcia G., Fernandez-Fernandez I. and Castro-Sanchez E. 2019. 'Incidence of peripheral intravenous catheter failure among inpatients: Variability between microbiological data and clinical signs and symptoms'. Antimicrobial Resistance and Infection Control 8(1):124.
<https://www.scopus.com/inward/record.uri?eid=s2-0-85069782323&doi=10.1186%2fs13756-019-0581-8&partnerID=40&md5=09596ab961e4fc1b1c777b8eca14664c>
4. Foster, Timothy J. 1996. 'Staphylococcus'. In Medical Microbiology. 4th Edition., 4th ed. Galveston.
5. Public Health England. 2015. 'Identification of Pseudomonas species and other Non-Glucose Fermenters'. ID 17 (3).
<https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>

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Bibliographie

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4).
<https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2016. 'Identification of Clostridium species'. ID 8 (4).
<https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>
3. Blanco-Mavillard I., Rodriguez-Calero M.A., De Pedro-Gomez J., Parra-Garcia G., Fernandez-Fernandez I. and Castro-Sanchez E. 2019. 'Incidence of peripheral intravenous catheter failure among inpatients: Variability between microbiological data and clinical signs and symptoms'. Antimicrobial Resistance and Infection Control 8(1):124.
<https://www.scopus.com/inward/record.uri?eid=2-s2.0-85069782323&doi=10.1186%2fs13756-019-0581-8&partnerID=40&md5=09596ab961e4fc1b1c777b8eca14664c>
4. Foster, Timothy J. 1996. 'Staphylococcus'. In Medical Microbiology. 4th Edition., 4th ed. Galveston.
5. Public Health England. 2015. 'Identification of Pseudomonas species and other Non-Glucose Fermenters'. ID 17 (3).
<https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>

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1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2016. 'Identification of Clostridium species'. ID 8 (4). <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>
3. Blanco-Mavillard I., Rodriguez-Calero M.A., De Pedro-Gomez J., Parra-Garcia G., Fernandez-Fernandez I. and Castro-Sanchez E. 2019. 'Incidence of peripheral intravenous catheter failure among inpatients: Variability between microbiological data and clinical signs and symptoms'. *Antimicrobial Resistance and Infection Control* 8(1):124. <https://www.scopus.com/inward/record.uri?eid=2-s2.0-85069782323&doi=10.1186%2f13756-019-0581-8&partnerID=40&md5=09596ab961e4fc1b1c777b8eca14664c>
4. Foster, Timothy J. 1996. 'Staphylococcus'. In Medical Microbiology. 4th Edition., 4th ed. Galveston.
5. Public Health England. 2015. 'Identification of Pseudomonas species and other Non-Glucose Fermenters'. ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>

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
Bibliografia

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4).
<https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2016. 'Identification of Clostridium species'. ID 8 (4).
<https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>
3. Blanco-Mavillard I., Rodriguez-Calero M.A., De Pedro-Gomez J., Parra-Garcia G., Fernandez-Fernandez I. and Castro-Sanchez E. 2019. 'Incidence of peripheral intravenous catheter failure among inpatients: Variability between microbiological data and clinical signs and symptoms'. *Antimicrobial Resistance and Infection Control* 8(1):124.
<https://www.scopus.com/inward/record.uri?eid=2-s2.0-85069782323&doi=10.1186%2fs13756-019-0581-8&partnerID=40&md5=09596ab961e4fc1b1c777b8ec-a14664c>
4. Foster, Timothy J. 1996. 'Staphylococcus'. In Medical Microbiology. 4th Edition., 4th ed. Galveston.
5. Public Health England. 2015. 'Identification of Pseudomonas species and other Non-Glucose Fermenters'. ID 17 (3).
<https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>

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Bibliografia

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2016. 'Identification of Clostridium species'. ID 8 (4). <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>
3. Blanco-Mavillard I., Rodriguez-Calero M.A., De Pedro-Gomez J., Parra-Garcia G., Fernandez-Fernandez I. and Castro-Sanchez E. 2019. 'Incidence of peripheral intravenous catheter failure among inpatients: Variability between microbiological data and clinical signs and symptoms'. *Antimicrobial Resistance and Infection Control* 8(1):124. <https://www.scopus.com/inward/record.uri?eid=2-s2.0-85069782323&doi=10.1186%2fs13756-019-0581-8&partnerID=40&md5=09596ab961e4fc1b1c777b8eca14664c>
4. Foster, Timothy J. 1996. 'Staphylococcus'. In *Medical Microbiology*. 4th Edition., 4th ed. Galveston.
5. Public Health England. 2015. 'Identification of Pseudomonas species and other Non-Glucose Fermenters'. ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>

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
Bibliografie

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2016. 'Identification of Clostridium species'. ID 8 (4). <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>
3. Blanco-Mavillard I., Rodriguez-Calero M.A., De Pedro-Gomez J., Parra-Garcia G., Fernandez-Fernandez I. and Castro-Sanchez E. 2019. 'Incidence of peripheral intravenous catheter failure among inpatients: Variability between microbiological data and clinical signs and symptoms'. *Antimicrobial Resistance and Infection Control* 8(1):124. <https://www.scopus.com/inward/record.uri?eid=2-s2.0-85069782323&doi=10.1186%2fs13756-019-0581-8&partnerID=40&md5=09596ab961e4fc1b1c777b8eca14664c>
4. Foster, Timothy J. 1996. 'Staphylococcus'. In Medical Microbiology. 4th Edition., 4th ed. Galveston.
5. Public Health England. 2015. 'Identification of Pseudomonas species and other Non-Glucose Fermenters'. ID 17 (3). <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>

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©2022 Thermo Fisher Scientific Inc. Reservados todos los derechos.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus filiales.
Oxoid Deutschland GmbH, Am Lippeglacis