Keycode AT07527B

www.thermofisher.com

Europe +800 135 79 135

US 1 855 2360 190

CA 1 855 805 8539

ROW +31 20 794 7071

# remel™ NITROCEFİN DISK

EN

REF R211667..... 25 Disks/Vial

## 1. INTENDED USE

The Nitrocefın Disk is a disk impregnated with nitrocefın for use in qualitative procedures for the detection of  $\beta$ -lactamase produced by bacterial isolates. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

## 2. SUMMARY AND EXPLANATION

The  $\beta$ -lactamase enzyme, originally described by Abraham and Chain, is produced by various organisms and is a mechanism of their resistance to penicillins and cephalosporins.<sup>1,2</sup> Since then, similar enzymes with somewhat different substrate specificities have been identified in various bacterial species. Some enzymes selectively hydrolyze penicillin class antibiotics (i.e., penicillinases), other enzymes hydrolyze all  $\beta$ -lactam antibiotics except carbapenems (i.e., cephalosporinases), and still other enzymes hydrolyze both.<sup>3</sup> Test methods used to detect  $\beta$ -lactamase include iodometric, acidimetric, and chromogenic procedures.<sup>4,5</sup>

## 3. PRINCIPLE

Nitrocefın, a cephalosporin first described by Glaxo Research, is the substrate used in this test.  $\beta$ -lactamase hydrolyzes the  $\beta$ -lactam ring of nitrocefın, producing cephalosporanic acid.<sup>6</sup> A distinctive color change is associated with this reaction wherein the pale yellow nitrocefın compound is converted to a pink end product upon hydrolysis.<sup>7</sup> Aerobic and anaerobic  $\beta$ -lactamase-producing bacteria effect this color change; organisms that do not produce  $\beta$ -lactamase do not alter the pale yellow color of nitrocefın within the time limits of the test.

## 4. REAGENTS

Nitrocefın Disk (25 Disks/Vial)

Reactive Ingredient: Nitrocefın

## 5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully.

Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

## 6. STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 2-8°C until used. Allow product to equilibrate to room temperature before use.

## 7. PRODUCT DETERIORATION

This product should not be used if (1) the disks are pink, (2) the expiration date has passed, (3) the desiccant has changed from blue to pink, or (4) there are other signs of deterioration. Protect disks from moisture by removing from vial only those disks necessary for testing. Promptly replace the cap and return the vial to 2-8°C.

## 8. SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following local recommended guidelines.<sup>5</sup>

## 9. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Forceps, (7) Microscope Slide, (8) Demineralized water, (9) Wooden applicator stick.

## 10. PROCEDURE

- Using forceps, place disk on a clean glass slide or in an empty petri dish lid.
- Moisten disk with a loopful of demineralized water.
- Pick 5-6 colonies from the isolate to be tested, or a sweep of confluent growth from a pure culture, using a sterile loop or wooden applicator stick.
- Smear the organism on the disk.
- Incubate the disk for 5 minutes at room temperature.

**Note:** Positive reactions may take up to 60 minutes to develop for some staphylococci. Anaerobic bacteria may require up to 30 minutes for a positive reaction to occur.

- Observe for a color change to pink or red.

## 11. INTERPRETATION

Positive Test - A color change to pink or red

Negative Test - No color change

12. EXPECTED VALUES

This test provides reliable results with *Neisseria gonorrhoeae*, *Haemophilus* spp., *Moraxella catarrhalis*, enterococci, and staphylococci. Specifically, resistance to penicillin, ampicillin, and amoxicillin is detected for *N. gonorrhoeae*, *Haemophilus* spp., and *M. catarrhalis*. Penicillin, ampicillin, amoxicillin, azlocillin, carbenicillin, mezlocillin, piperacillin, and ticarcillin resistance is detected for enterococci and staphylococci.<sup>8</sup>

Testing for resistance to penicillin and ampicillin using the chromogenic cephalosporin method has been shown to be predictive of resistance for *Prevotella*, *Porphyromonas*, *Bacteroides*, and other anaerobes.<sup>9</sup> The great majority of *Bacteroides fragilis* group isolates are  $\beta$ -lactamase producers, and therefore do not require routine  $\beta$ -lactamase testing and should be considered resistant. It should be understood, however, that resistance to  $\beta$ -lactam drugs is not always mediated by  $\beta$ -lactamase production (e.g., some strains of *Parabacteroides distasonis*, *B. fragilis*, *Bacteroides gracilis*, and *Bilophila wadsworthia*).<sup>10</sup>

Resistance due to other mechanisms will not be detected with the use of this product.

13. QUALITY CONTROL

All lot numbers of Nitrocefim Disk have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported. Quality control testing should be run with each shipment and new lot number received.

CONTROL	INCUBATION	RESULTS
<i>Staphylococcus aureus</i> ATCC® 29213	Ambient, 5 minutes @ 25°C	Positive
<i>Neisseria gonorrhoeae</i> ATCC® 43069	Ambient, 5 minutes @ 25°C	Negative

14. PERFORMANCE CHARACTERISTICS <sup>11</sup>

In a study of 103 test isolates (4 *E. faecalis*, 13 *M. catarrhalis*, 19 *H. influenzae*, 17 *S. aureus*, 21 coagulase- negative staphylococci, 12 *N. gonorrhoeae*, 11 *B. fragilis*, 5 *B. thetaiotaomicron*, and 1 *Campylobacter* spp.), there was 100% correlation between Nitrocefim Disk and two similar, commercially available products. 61.4% of the organisms included in the study were  $\beta$ -lactamase positive and 35.9% were  $\beta$ -lactamase negative.

15. LIMITATIONS

- 1. Nitrocefim Disk is intended for use as a rapid test and should not replace conventional susceptibility testing.
- 2.  $\beta$ -lactamase negative,  $\beta$ -lactam antibiotic resistant organisms have been reported. An organism's resistance due to other mechanisms will not be detected.
- 3. Do not test members of the Enterobacteriaceae, *Pseudomonas* spp., and other aerobic, Gram-negative bacilli, as results may not be predictive of susceptibility to the  $\beta$ -lactams most often used in therapy.


- 4. Some staphylococci may require induction to increase production of  $\beta$ -lactamase to detectable levels. A negative  $\beta$ -lactamase test result from a non-induced staphylococcal isolate should be confirmed after induction, which can be easily accomplished by testing the peripheral growth of the zone surrounding a 1 $\mu$ g oxacillin disk.

16. BIBLIOGRAPHY

- 1. Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
- 2. Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
- 3. Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- 4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
- 5. Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
- 6. Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
- 7. O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
- 8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
- 9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
- 10. Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
- 11. Data on file

## 17. SYMBOL LEGEND

Symbol	Definition
	Catalogue number
	<i>In Vitro</i> Diagnostic Medical Device
	Batch code
	Temperature limitation (storage temp.)
	Use by (expiration date) YYYY-MM
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Manufacturer
	Authorized representative in the European Community
	UK Conformity Assessment

Version	Date of issue and modifications introduced
1.0	2022-03-20. New document

©2022 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS  
66215, USA

For technical assistance please contact your local distributor.

[www.thermofisher.com/microbiology](http://www.thermofisher.com/microbiology)

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- 1. Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
- 2. Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
- 3. Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman’s Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- 4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
- 5. Jorgensen J.H., M.A. Pfaller et al. 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
- 6. Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
- 7. O’Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
- 8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
- 9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
- 10. Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
- 11. Данни във файла

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## 16. BIBLIOGRAFIJA

1. Abraham, E. P. and E. Chain. 1940. Nature. 146:837.
2. Del Bene, V. E. and W. E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
3. Winn, C. W., S. D. Allen, W. M. Janda, E. W. Koneman, G. W. Procop, P. C. Schreckenberger, and G. L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
4. Montgomery, K., L. Raymundo, Jr., and W. L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
5. Jorgensen J. H., M. A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
6. Lee, D. T. and J. E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
7. O'Callaghan, C. H., A. Morris, S. M. Kirby, and A. H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
10. Finegold, S. M., E. J. Baron, and H. M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
11. Podaci su dostupni na upit

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
## 16. LITTERATUR

- Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
- Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
- Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
- Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
- Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
- O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
- Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
- Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
- Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
- Data i arkiv

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[www.thermofisher.com/microbiology](http://www.thermofisher.com/microbiology)  
Tlf.: (800) 255-6730 • International: (913) 888-0939

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3. Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.

4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.

5. Jorgensen J.H., M.A. Pfaller et al. 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.

6. Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.

7. O'Callaghan, C.H., A. Morris, S.M. Kirby ja AH Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.

8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI

9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.

10. Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.

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**16. BIBLIOGRAPHIE**

1. Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
2. Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
3. Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
5. Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
6. Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
7. O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
10. Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
11. Données archivées

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- Abraham, E.P. and E. Chain. 1940. *Nature*. 146:837.
- Del Bene, V.E. and W.E. Farrar, Jr. 1973. *Antimicrob. Agents Chemother.* 3(3):369-372.
- Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. *Koneman's Color Atlas and Textbook of Diagnostic Microbiology*. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. *J. Clin. Microbiol.* 9(2):205-207.
- Jorgensen J.H., M.A. Pfaller *et al.* 2015. *Manual of Clinical Microbiology*. 11th edition. Washington ASM Press.
- Lee, D.T. and J.E. Rosenblatt. 1983. *Diagn. Microbiol. Infect. Dis.* 1:173-175.
- O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. *Antimicrob. Agents Chemother.* 1(4):283-288.
- Clinical and Laboratory Standards Institute (CLSI). 2018. *Performance Standards for Antimicrobial Disk Susceptibility Tests*. 13th Edition, M02-A13. CLSI
- Clinical and Laboratory Standards Institute (CLSI). 2018. *Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria*. 9th Edition. M11-A9. CLSI.
- Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. *A Clinical Guide to Anaerobic Infections*. Star Publishing Co., Belmont, CA.
- Δεδομένα αρχείου (data on file)

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## 16. BIBLIOGRÁFIA

1. Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
2. Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
3. Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
5. Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
6. Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
7. O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
10. Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
11. Adatok a fájlon

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16. LITERATŪRA

- 1. Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
- 2. Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
- 3. Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman’s Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- 4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
- 5. Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
- 6. Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
- 7. O’Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
- 8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
- 9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
- 10. Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
- 11. Turimi duomenys

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## 16. BIBLIOGRAFIA

- Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
- Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
- Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
- Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
- Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
- O'Callaghan, CH, A. Morris, SM Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
- Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
- Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
- Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
- Dane w pliku

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
## 16. BIBLIOGRAFIE

1. Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
2. Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
3. Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
5. Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
6. Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
7. O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
10. Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
11. Date din dosar

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## 16. ZDROJE

1. Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
2. Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
3. Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
5. Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
6. Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
7. O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
8. Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
9. Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
10. Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
11. Údaje v záznamoch

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16. BIBLIOGRAFÍA

- Abraham, E.P. and E. Chain. 1940. Nature. 146:837.
- Del Bene, V.E. and W.E. Farrar, Jr. 1973. Antimicrob. Agents Chemother. 3(3):369-372.
- Winn, C.W., S.D. Allen, W.M. Janda, E.W. Koneman, G.W. Procop, P.C. Schreckenberger, and G.L. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. J. Clin. Microbiol. 9(2):205-207.
- Jorgensen J.H., M.A. Pfaller *et al.* 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
- Lee, D.T. and J.E. Rosenblatt. 1983. Diagn. Microbiol. Infect. Dis. 1:173-175.
- O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Antimicrob. Agents Chemother. 1(4):283-288.
- Clinical and Laboratory Standards Institute (CLSI). 2018. Performance Standards for Antimicrobial Disk Susceptibility Tests. 13th Edition, M02-A13. CLSI
- Clinical and Laboratory Standards Institute (CLSI). 2018. Methods for Antimicrobial Susceptibility Testing of Anaerobic Bacteria. 9th Edition. M11-A9. CLSI.
- Finegold, S.M., E.J. Baron, and H.M. Wexler. 1992. A Clinical Guide to Anaerobic Infections. Star Publishing Co., Belmont, CA.
- Datos de archivo.

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©2022 Thermo Fisher Scientific Inc. Reservados todos los derechos.
Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, EE. UU.