



Add: No.4, Baolong 6th Road, Baolong Industrial Town,
Longgang District 518116, Shenzhen, China
TEL: 86-755-89312160 FAX: 86-755-89312239

We, SCW Medicath Ltd. located in: No.4, Baolong 6th Road, Baolong Industrial Town, Longgang District Shenzhen 518116 Guangdong, China hereby declare for we change catalog number for guidewires because of making new catalog, please see the code list as follows:

Labels = To complete – table for Old part number to new part number

	Compatible for the use of the following injector models, which meet the provisions of Directive 93/42/EEC	SCW old part number	SCW content	SCW New part number
DSA-150-A1	Medrad Mark V/Medrad Mark V ProVis Angiographic Injection System	832603	1- syringe 150 ml 1- J fill tube	42.16.10004
DSA-150-A2	Medrad Mark 7		1- syringe 150 ml 1- J fill tube	42.16.10016
CT-200-A1	Medrad Vistron CT Injection System Medrad EnVision CT Injection System	8326601	1- syringe 200 ml 1- connecting tubing 150 cm 1- J fill tube	42.16.10000
MRI-65/115-A1	Medrad Spectris Solaris EP MR Injection System	832623	1- syringes 65ml 1-syringe 115ml 2- J fill tube 1- Y(T) connecting tubing 250 cm 2- vented spike	42.16.10010
CT-200/200-A2	Medrad Stellant D CT Injection System	832602	1- syringe 200 ml 1- connecting tubing 150 cm 1- J fill tube	42.16.10002
CT-200-A2	Medrad Stellant D Dual Syringe CT Injection System	832620A	1- syringe 200 ml 1- Y(T) connecting tubing 150 cm 1- J fill tube 1- vented spike	42.16.10013
		832620B	1- syringe 200 ml 1- J fill tube 1- vented spike	42.16.10014

CT-200-B1	Libel-Flarsheim CT9000	832609	1- syringe 200 ml 1- J fill tube	
	Libel-Flarsheim CT9000 ADV Mallinckrodt Optivantage		1- connecting tubing 150 cm	42.16.20000
CT-200-B1	Libel-Flarsheim CT9000	832609Y	1- syringe 200 ml 1- J fill tube	42.16.20040
	Libel-Flarsheim CT9000 ADV Mallinckrodt Optivantage		1- Y(T) connecting tubing 150 cm	
CT-200/200-B1	Libel-Flarsheim CT9000	832621	2- syringe 200 ml 2 - J fill tube	42.16.20001
	Libel-Flarsheim CT9000 ADV Mallinckrodt Optivantage		1- Y(T) connecting tubing 150 cm 2- vented spikes	
MRI-60/60-B1	Libel-Flarsheim Optistar	832613	2- syringes 60 ml 1- Y(T) connecting tubing 250 cm	42.16.20006
	Liebel-Flarsheim MR OptiStar LE Injector System		2- J fill tubes 2- vented spikes	
MRI-60/60-NE	Nemoto Sonic Shot	832618	2- syringes 60 ml 1- Y(T) connecting tubing 250 cm 2- J fill tubes 2- vented spikes	42.16.30007
CT-200-NE	Nemoto Dual-Shot Injector A-25	832615	1- syringe 200 ml 1- connecting tubing 150 cm	
	Nemoto Dual Shot CT-Injector A-60		1- J fill tube	42.16.30002
CT-100/200-NE	Nemoto Dual-Shot Injector	832632	1- syringes 100 ml 1- syringes 200 ml 1- Y(T) connecting tubing 150 cm 2- J fill tubes 2- vented spikes	42.16.30008



CT-200-NE	Nemoto Dual-Shot Injector	832632A	1- syringe 200 ml 1- Y(T) connecting tubing 150 cm 1- J fill tube 1- vented spike	42.16.30009
CT-100-NE	Nemoto Dual Shot CT-Injector	832632B	1- syringe 100 ml 1- J fill tube 1- vented spike	42.16.30010
CT/DSA-150-A1	Angiomat Illumena	632611	1- syringe 150 ml 1- J fill tube 1- connecting tubing 150 cm	42.16.20003
Tubing - 0.12inch	Medrad Vistron CT Nemoto Dual Shot A-25 Libel-Flarsheim CT9000, CT9000 ADV, Optivantage Medrad Mark V ProVis Angiomat Illumena	832653	connecting tube 0.12", 150cm ,375psi	42.02.90001
Tubing - 0.15inch	Medrad Stellant D Dual Syringe CT Nemoto Dual Shot Libel-Flarsheim CT9000 DH, Optivantage	832651	connecting tube 0.15",150cm ,375psi	42.02.90000
Tubing - 0.15inch	Medrad Spectris Solaris EP MR Nemoto Sonic Shot Libel-Flarsheim Optistar	832654	connecting tube 0.15",250cm ,375psi	42.02.90002