

Ulrich GmbH & Co. KG | Buchbrunnenweg 12 | 89081 Ulm | Germany

Ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm, Germany

To whom it may concern

Renner, Nadine
Regulatory Affairs Manager

T. +49 (0)731 9654-2975
n.renner@ulrichmedical.com

Ulm, 2024-12-16

Declaration of compatibility of XD 10711, XD 10710, XD 10701, XD 10702, XD 10741 with the ulricheasy-INJECT Max 3

I hereby declare the compatibility between

- the XD 10701: Patient tubing, 250 cm for ulricheasyINJECT CT/MRT injectors
- the XD 10702: Patient tubing, 320 cm for ulricheasyINJECT CT/MRT injectors
- the XD 10710: Easy-Click-Cassette Max 3 - holder s for ulricheasyINJECT MRI injectors
- the XD 10711: Easy-Click-Cassette Max 3 - holder l for ulricheasyINJECT MRI injectors
- the XD 10741: Adapter ring for small vials

and the XD 10160 ulricheasy INJECT Max 3.

Ulrich medical



Corina Matzdorf
Head of Regulatory

Ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm
Germany
Tel. +49 (0) 731 96 54-0
Fax +49 (0) 731 96 54-199

ulrich GmbH & Co. KG | Buchbrunnenweg 12 | 89081 Ulm | Germany

ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm, Germany

To whom it may concern

Renner, Nadine
Regulatory Affairs Manager

T. +49 (0)731 9654-2975
n.renner@ulrichmedical.com

Ulm, 2024-12-16

**Declaration of compatibility of XD 10713, XD 10701, XD 10702, XD 10741 with the
ulricheasy-INJECT Max 2M**

I hereby declare the compatibility between

- the XD 10701: Patient tubing, 250 cm for ulricheasyINJECT CT/MRT injectors
- the XD 10702: Patient tubing, 320 cm for ulricheasyINJECT CT/MRT injectors
- the XD 10713: Easy-Click-Cassette Max 2M - holder I for ulricheasyINJECT MRI injectors
- the XD 10741: Adapter ring for small vials

and the XD 10130 ulricheasy INJECT Max 2M.

ulrich medical



Corina Matzdorf
Head of Regulatory



ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm
Germany
Tel. +49 (0) 731 96 54-0
Fax +49 (0) 731 96 54-199

ulrich GmbH & Co. KG | Buchbrunnenweg 12 | 89081 Ulm | Germany

ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm, Germany

To whom it may concern

Renner, Nadine
Regulatory Affairs Manager

T. +49 (0)731 9654-2975
n.renner@ulrichmedical.com

Ulm, 2024-12-16

Declaration of compatibility of XD 8151, XD 8003, XD 2035, XD 2040 and XD 2045 with the ulrichINJECT CT motion injector

I hereby declare the compatibility between

- the XD 8151 Pump tubing
- the XD 8003 Pump tubing
- the XD 2035 Patient tubing 150 cm, 2 check valves
- the XD 2040 Patient tubing 250 cm, 2 check valves
- the XD 2045 Patient tubing 320 cm, 2 check valves

and the XD 8000 ulrichINJECT CT motion. The two different pump tubings are compatible with different software versions of the injector:

	≤ SW 3.4	SW 4.0	SW 4.1	SW 4.3	SW 4.5	SW 4.6	SW 6.1
XD 8151	X	✓	✓	✓	✓	✓	✓
XD 8003	✓	✓	✓	✓	✓	✓	X

ulrich medical



Corina Matzdorf
Head of Regulatory



ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm
Germany
Tel. +49 (0) 731 96 54-0
Fax +49 (0) 731 96 54-199

ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm | Germany
T. +49 (0)731 9654-0
www.ulrichmedical.com
UID DE147049042

Amtsgericht Ulm HRA 222
Komplementär: ulrich Verwaltungs GmbH
Amtsgericht Ulm HRB 3209
Geschäftsführender Gesellschafter,
Managing Partner: Christoph Ulrich
Geschäftsführer, CEO: Klaus Kiesel
Geschäftsführer, Co-CEO: Friedrich von Rechteren

Bankverbindungen, Bank accounts:
Baden-Württembergische Bank AG Ulm
IBAN DE47 6005 0101 7439 5033 53
BIC (SWIFT-Code) SOLADEST600

Sparkasse Ulm
IBAN DE68 6305 0000 0000 1324 11
BIC (SWIFT-Code) SOLADES1ULM

Volksbank Ulm-Biberach eG
IBAN DE09 6309 0100 0001 4040 08
BIC (SWIFT-Code) ULMVDE66