



A/V Linie do Dialog+

steryl. E-BEAM

bez DEHP

DiaStream[®]
A/V SYSTEM LINII KRWI

DiaStream do Dialog+

NUMER KAT.	NAZWA	ILOSC W OP.	WYPELNIENIE OBJETOSC (ML)	TETNICZA OBJETOSC (ML)	ZYLNA OBJETOSC (ML)	Ø ROZM. POMPY (MM)	SREDNICA LINII KRWI A i V (MM)	SREDNICA LINII HEPAR. (MM)
DiaStream HD								
7211508	DiaStream Basic HD	25	140	71	69	8	4.3 / 4.8	1.5
<u>7211509</u>	DiaStream Classic HD	25	140	71	69	8	4.3 / 4.8	1.5
7211510	DiaStream Advanced HD	25	140	71	69	8	4.3 / 4.8	1.5
7211511	DiaStream Premium HD	25	140	71	69	8	4.3 / 4.8	1.5
DiaStream HDF								
7211512	DiaStream Classic HDF	25	139	70	69	8	4.3 / 4.8	1.5
DiaStream Universal HD								
7211513	DiaStream Classic Universal HD	25	151	76	75	8	4.3 / 4.8	1.5

NR KATALOGOWY	Ø KOMORA TETNICZA (MM)	Ø KOMORA ZYLNA (MM)	PRE/POST DYLUCJA	WOREK (ML)	SPIKE	RECYRKULACJA LACZNIK	INFUSION LINIA
DiaStream HD							
7211508	22	22	Post	–	–	–	–
<u>7211509</u>	22	22	Post	2000	■	–	–
7211510	22	22	Post	2000	■	■	–
7211511	22	22	Post	2000	■	■	■
DiaStream HDF							
7211512	22	22	Pre / Post	–	–	–	–
DiaStream Universal HD							
7211513	–	21.5 / 30	–	2000	■	–	–